INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004935

A. GENERAL Setup and Checks:

1. Dever On

Ž 5 min. countdown during initial warm-up

- X 1 min. countdown to begin test after standby mode
- 2. X Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security: Level 1 (Operator) X Level 2 (Field Inspector) X Level 3 (Lab Personnel)
- 4. X Printer Setup (Level 3, S, P)
- 5. X Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

- 1. X Internal Printer Test and Paper Feed (Level 1, P). Attach record.
- 2. X External Printer Test (Level 1,P). Attach record.
- 3. Barcode Reader (Level 1,B): Operator ID 🕅 ND License 🛱 R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

- 1. X Continuous Air Blank (Level 2,A)
- 2. XDiagnostic (Level 2,D). Attach record.
- 3. X Prelim Data Entry (Level 2,E)
- 4. Ser Number/Version Info (Level 2,V): Inst Ser#: <u>80-004935</u> DSP Ver: <u>1362-14</u> Core Ver: <u>8164-14</u> CPLD Ver: <u>7</u>

D. LEVEL 3 Setup and Checks:

- 1. DVM Monitor (Level 3,D,D): 3: $\underline{12737}$ s: $\underline{1}$ C: $\underline{47.00}$ 9: $\underline{12957}$ s: $\underline{2}$ B: $\underline{45}$ O (RFI indicator with keyed radio): $\underline{1}$ 2. Function (F) Keys: : KF5 – Tone KF9 - Fan F6 – Green LED F10 – Pump F7 – Red LED F11 – Solenoid
- Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: <u>945</u> CF: <u>1.07</u>
 Flow Monitor (Level 3,D,F): Flow Sensor, R: <u>93</u> F: <u>0.000</u> L/s

BrW-003

- 5. Temperature Monitor (Level 3,D,T): C: $\frac{47.0}{B}$ A: $\frac{46.5}{V}$ I: 0.0B: $\frac{45.0}{V}$ V: 0.0
- 6. Tank Monitor (Level 3,D,G): P (display): 236 psi P (gauge): 225 psi
- 7. General Setup (Level 3,S,G)
- 8. X Configure Start Test (Level 3,S,T)
- 9. X Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

- 1. Zero AC (Level 1,C): Sim ser #: <u>MP53·20</u> Lot #: <u>MILLI Q H₂0</u> Exp date: <u>NA</u> Known value: <u>0.000</u> AC Results: <u>0.000</u> AC <u>0.000</u> AC <u>0.000</u> AC
- 2. Low AC (Level 1,C): Sim ser #: <u>MP5289</u> Lot #: <u>201911E</u> Exp date: <u>11.19.21</u> Known value: <u>0.030</u> AC Results: <u>0-028</u> AC <u>0.030</u> AC <u>0.030</u> AC
- 3. High AC (Level 1,C): Sim ser #: <u>MP306</u>∠ Lot #: <u>202/03</u>E Exp date: <u>3.24.23</u> Known value: <u>0.400</u> AC Results: <u>0.403</u> AC <u>0.400</u> AC <u>0.401</u> AC

Configure dry gas standard before running tests (Level 1,S):

- 4. Ethanol + 0.05% Acetone (Level 1,B): Sim ser #: <u>DR3847</u> Lot #: <u>ICS6</u> Exp date: <u>NA</u> Known value: <u>0.10AC + 0.057. ACETONE</u>
 ☑ Display reads "Interferent Detect"
- 5. Acetone (Level 1,B): Sim ser #: <u>DR5188</u> Lot #: <u>ACE1</u> Exp date: <u>NA</u> Known value: <u>0.057</u> <u>ACETONE</u>
- Display reads "Interferent Detect"
 Methanol (Level 1,B): Sim ser #: <u>DR5186</u> Lot #: <u>MET3</u> Exp date: <u>NA</u> Known value: <u>0.10% MET7-IANOL</u>
 Display reads "Interferent Detect"
- 7. Isopropanol (Level 1,B): Sim ser #: <u>DR5189</u> Lot #: <u>IPAZ</u> Exp date: <u>NA</u> Known value: <u>0.10%</u> <u>Isopropavoc</u>
 X Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. X Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. X ABA (Level 1,B): Test record reads correct.
- 12. X Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. XRFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16 X Improper Sample (CMS mode): Display and test record read correct.
- 17. X Difference Too Great (CMS mode): Display and test record read correct.
- 18. E-CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3, T, S):

Ave: 0. 0782 Std Dev: 0. 0006 Rel Std Dev: 0. 76%

Maintenance/Repairs/Notes: INITIAL INSPECTION DUE TO REPEACHE CEE REAL TIME CLOCK BATTERY AND HAVING TO PPEN UP INSTICUTE FORT PERFORM THE REPAIR. LNSTRUMENT

X Instrument is acceptable to be used in the field.

Field Inspector's Signature

Reviewed by

<u>4.28.21</u> 5.3-21

Date

Page 3 of 3

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Norti Dakota Model 8000 Location = TOXL 04/25/2021	Alcohol Analyzer SN 80-004935 8164.14.00 09/16
04/20/2021	12:03

********* Printer Test ********

12 char/line abcdefghijkl ABCDEFGHIJKL

16 char per line abcdefghijklmnop ABCDEFGHIJKLMNOP

21 chars per line abcdefghijklmnopqrstu ABCDEFGHIJKLMNOPQRSTU

24 characters per line abcdefghijklmnopqrstuvwx ABCDEFGHIJKLMNOPQRSTUVWX

32 characters per line abcdefghijklmnopqrstuvwxyz123456 ABCDEFGHIJKLMNOPQRSTUVWXYZ123456

42 characters per line abcdefghijklmnopqrstuuwxyz123456789Dabcdef ABCDEFGHIJKLMNOPDRSTUUWXYZ123456789DABCDEF

Current Instrument 1 Data Entry Mode: Start Test Sequence Display Prelim RsIt Display Third Digit Inhib Printer(Y/N)? Display Volume? Disable On Memfull? H of Print Copies? Select Std (D/W/I)? Standard Value? Standard Lot H? Standard Cyl H? Standard Expiration Oper No?	. ? ?	Enabled DABACABA Yes Yes No Yes I Dry 0.080 34418060 5 02/05/20 666666	IA2
Flow Cal. Date: Slope Intercept		07/12/20 655 -510289	112
IR Calibration Date	e: 3u⊓	05/20/2	015 9cm
Oth Coef(*100):	-7871		-19684
lst Coef(*100):	26450	5	133108
2nd Coef(*100):	2865		1320
H2O adj(mg/l*lOk):	200		298

********** Test End **********

Operator Signature CHARLES EDER

Remarks: INTERNAL

PRINT TEST

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL8164.14.00 09/16 04/26/2021 12:05 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? NO Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 5 Standard Expiration? 02/05/2021 Oper No? 666666 Flow Cal. Date: 07/12/2012 Slope 655 Intercept -510289 IR Calibration Date: 05/20/2015 3um 9um _____ 0th Coef(*100): -7871 -19684 264505 2nd Coef(*100): 2865 H2O add (* 17 133108 1320 H2O adj(mq/l*10k): 200298

Operator Signature CHARLES EDER

ENAL

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 12:08 Test AC Time 01 Room Air 0.000 12:09 02 *Subject Test REF* 12:09 03 Room Air 12:10 0.000 *Subject Test Refused Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = R KEY REFUSAL Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 5

Expiration Date = 02/05/2021 County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

"R"KEY REFUSAL

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004935Location = TOXL8164.14.00 09/1604/26/202112:11

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666

tel 1 e

Operator Signature CHARLES EDER

GAUSTIC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004935
Location = TOXL	8164.14.00 09/16
04/26/2021	12:21

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.000 0.000 0.000 0.000 0.000 0.000	12:21 12:22 12:23 12:23 12:24 12:24 12:25

08 Sim Temp = 34.0°C

Simul Ser No = MP5320 Std Sol No = MILLIQ H2O County = 08Oper No. = 666666

Operator Signature CHARLES EDER

ZERO AC

0.000AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004935
Location = TOXL	8164.14.00 09/16
04/26/2021	14:59

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.028 0.000 0.030 0.000 0.030 0.000	15:00 15:01 15:02 15:03 15:03 15:04

08 Sim Temp = 34.0°C

Simul Ser No = MP5289 Std Sol No = 201911E County = 08

Oper No. = 666666

1

Operator Signature CHARLES EDER

LOW A 0.030AC

CEE Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004935
Location = TOXL	8164.14.00 09/16
04/26/2021	15:05

	WET CAL CHECK]
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air	0.000 0.403 0.000 0.400 0.000	15:05 15:06 15:07 15:07 15:08
06 Std. Sol. 07 Room Air 08 Sim Temp =	0.401 0.000 34.0°C	15:09 15:09

Simul Ser No = MP3062 Std Sol No = 202103ECounty = 08

Oper No. = 666666

UU e. Ce d 1

Operator Signature CHARLES EDER

GAAAC 0.400 M

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL8164.14.00 09/16 04/26/2021 16:34 Test AC Time 01 Room Air 0.000 16:35 02 *Subject Test INT* 16:35 03 Room Air 0.000 16:36 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER + 0.05% ACETONE pley le and

Remarks:

Form 106-18000

0.10AC

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 16:37 Test AC Time 01 Room Air 0.000 16:38 02 *Subject Test INT* 16:38 03 Room Air 0.000 16:39 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER Remarks: Inferferen Check 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 16:39 Test AC Time 01 Room Air 0.000 16:40 02 *Subject Test INT* 16:40 03 Room Air 0.000 16:41 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER

Check 0.10% METHANOL REE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 16:42 Test AC Time 01 Room Air 0.000 16:43 02 *Subject Test INT* 16:43 03 Room Air 0.000 16:44 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Weight = 150Sub Sex = Female Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER

-format Check ISOPROPANOL

0.10%

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 16:45 Test AC Time 01 Room Air AMB* 16:46 02 Room Air 0.000 16:46 *Invalid Test Check Ambient Conditions Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = AMBIENT FAIL Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER

MBIEN

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL8164.14.00 09/16 04/26/2021 16:47 Test AC Time 01 Room Air 0.000 16:48 02 Invalid Sample X.XXX 16:48 03 Room Air 0.000 16:49 *Invalid Test - Mouth Alcohol Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INVALID SAMPLE Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER

WYALIO SAMPLE

MOUTH ALCOHOL

Form 106-18000

Intoxilyzer Test NDOAG Crime Lab. D			1
CMI, Inc. Intoxily North Dakota Model Location = TOXL 04/26/2021	8000	Alcohol Analyzen SN 80-004935 8164.14.00 09/16 16:49	5
Test	AC	Time	9
01 Room Air 02 *Subject Test 03 Room Air	0.000 ABT* 0.000	16:50)
*Sequence Aborted			

Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992 Sub Sex = Female Weight = 150 Test = DUI Cit = SEQUENCE ABORT Dr. Lic. = ND/DIS921456 Lot No = 02621080A1 Cyl No = 25 Expiration Date = 03/05/2023 County = 08 Oper No. = 666666

Operator Signature CHARLES EDER SEQUENCE ABORTED

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/26/2021 16:51 Test AC Time 01 Room Air 0.000 16:51 02 Subject Test 1 0.148 16:52 03 Room Air 0.000 16:52 04 Reported AC 0.148 16:52

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992 Sub Sex = Female Weight = 150 Test = DUI Cit = ABA TEST Dr. Lic. = ND/DIS921456 Lot No = 02621080A1 Cyl No = 25 Expiration Date = 03/05/2023 County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

1

ES1

Remarks:

Form 106-I8000

Intoxilyzer Test R		
NDOAG Crime Lab. Div	., Bismarck, N	D 58501
CMI, Inc. Intoxilyze North Dakota Model 8 Location = TOXL 04/27/2021		-004935
Test	AC	Time
03 Subject Test 1	0.000	16:23 16:23 16:24 16:26 16:27 16:28
*Cal Check Out of To	lerance	
Sub Name = DISCOVER, Sub DOB = 02/01/1992 Sub Sex = Female Test = DUI Dr. Lic. = ND/DIS9214 Lot No = 02621080A1	Weight Cit = CAL CHK (z = 150 DUT TOL

Lot No = 02621080A1 Cyl No = 25 Expiration Date = 03/05/2023 County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

lack e

Operator Signature CHARLES EDER

CALIBRATION CIVECK OUT

Form 106-18000

Intoxilyzer Test Re NDOAG Crime Lab. Div.		
CMI, Inc. Intoxilyzer North Dakota Model 80 Location = TOXL 04/27/2021	000 S	ol Analyzer N 80-004935 14.00 09/16 16:29
Test	AC	Time
03 *Subject Test	OK 0.000 RFI* 0.000	16:29 16:30 16:30 16:31
Sub Name = DISCOVER, Sub DOB = 02/01/1992 Sub Sex = Female Test = DUI Dr. Lic. = ND/DIS9214 Lot No = 02621080A1 Cyl No = 25 Expiration Date = 03/ County = 08	₩ 56 05/2023	eight = 150 Cit = NA o. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Ua

Operator Signature CHARLES EDER

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL8164.14.00 09/16 04/27/2021 16:32 Test AC Time 01 Diagnostic OK 16:33 02 Room Air 0.000 16:34 03 *Subject Test 0.000* 16:37 04 Room Air 0.000 16:39 05 Std. Gas 0.078 16:40 06 Room Air 0.000 16:42 07 *Subject Test 0.000* 16:45 08 Room Air 0.000 16:45 No RFI Detected *Deficient Sample - Value Printed was Highest Obtained Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = DEFICIENT SAMP Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 66666620 minute waiting period ascertained? Y Clean Mouthpiece used and disposed? Y Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

FICIENT SAMPLES

Form 106-I8000

Intoxilyzer Test NDOAG Crime Lab. Di		
CMI, Inc. Intoxilyz North Dakota Model Location = TOXL 04/27/2021	8000 SN	l Analyzer 80-004935 4.00 09/16 16:49
Test	AC	Time
01 Diagnostic 02 Room Air 03 Reference 04 Room Air	OK 0.000 IPS* 0.000	16:50 16:51 16:51 16:51

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992 Sub Sex = Female Weight = 150 Test = DUI Cit = IMPROPER SAMPLE Dr. Lic. = ND/DIS921456 Lot No = 02621080A1 Cyl No = 25 Expiration Date = 03/05/2023 County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

lack

Operator Signature CHARLES EDER

Remarks:

NPROPER SAMPLE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL 8164.14.00 09/16 04/27/2021 16:52 Test AC Time 01 Diagnostic OK 16:53 02 Room Air 0.000 16:53 03 Subject Test 1 0.199* 16:53 04 Room Air 0.000 16:56 05 Std. Gas 0.078 16:57 06 Room Air 0.000 16:58 07 Subject Test 2 0.148* 16:59 08 Room Air 0.000 17:00 *Difference Too Great

Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150 Test = DUI Cit = DIFF TOO GREAT Dr. Lic. = ND/DIS921456 Lot No = 02621080A1Cyl No = 25 Expiration Date = 03/05/2023County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

lart

Operator Signature CHARLES EDER

100 CENCE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004935 Location = TOXL8164.14.00 09/16 04/28/2021 08:57 Test AC Time 01 Diagnostic OK 08:58 02 Room Air 0.000 08:59 03 Subject Test 1 0.147 08:59 04 Room Air 0.000 09:01 05 Std. Gas 0.077 09:02 06 Room Air 0.000 09:04 07 Subject Test 2 0.147 09:04 08 Room Air 0.000 09:05 09 Reported AC 0.147 09:04 Difference OK No RFI Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = CMS TEST Dr. Lic. = ND/DIS921456Lot No = 02621080A1Cyl No = 25Expiration Date = 03/05/2023County = 08Oper No. = 66666620 minute waiting period ascertained? Y Clean Mouthpiece used and disposed? Y Clean Mouthpiece used and disposed? Υ

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

all

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004935
Location = TOXL	8164.14.00 09/16
04/28/2021	09:06

Т	est		AC	Time
01	Room		0.000	09:07
02	Std.	Gas	0.078	09:07
03	Room	Air	0.000	09:08
04	Std.	Gas	0.077	09:08
05	Room	Air	0.000	09:08
06	Std.	Gas	0.078	09:09
07	Room	Air	0.000	09:09
08 09	Std. Room	Gas Air	0.079 0.000	09:10
10	Std.	Gas	0.079	09:10
11	Room	Air	0.009	09:10 09:11
12	Std.	Gas	0.078	09:11
13	Room	Air	0.000	09:11
14	Std.	Gas	0.078	09:12
15	Room	Air	0.000	09:12
16	Std.	Gas	0.079	09:13
17	Room	Air	0.000	09:14
18	Std.	Gas	0.078	09:14
19	Room	Air	0.000	09:14
20	Std.	Gas	0.077	09:15
21	Room	Air	0.000	09:15
22	Std.	Gas	0.078	09:16
23	Room	Air	0.000	09:16
24	Std.	Gas	0.078	09:16
25	Room	Air	0.000	09:17
26	Std.	Gas	0.078	09:17
27	Room	Air	0.000	09:18
28	Std.	Gas	0.078	09:18
29	Room	Air	0.000	09:19
30	Std.	Gas	0.079	09:19
31 32	Room Std.		0.000	09:20
33	Room	Gas Air	0.078 0.000	09:20 09:20
34	Std.	Gas	0.078	09:20
35	Room	Air	0.000	09:21
	Std.		0.078	09:21
37			0.000	09:22
38			0.079	09:22
39	Room		0.000	09:23
40	Std.		0.078	09:23
41	Room	Air	0.000	09:24
42	Cal C	heck Stats	3:	
43	Aver	-	0.0782	
44	Std		0.0006	
45	Rel	Std Dev	0.76%	

Charlet El

Operator Signature CHARLES EDER

STABILITY TEST 0,080 AC

Form 106-18000

LOT: 02621080A1 EXP: 3.5-23 CYL: 25