

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-005957 Reason for Install/Repair	MAND	
☐ Install After Receiving From Crime Laboratory ☐ Ins	tall After Location Change	
Other (Specify) After Angual Tasp.		
HTTE Aller List.		
Check When Done:		
2 1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and	Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function	on S).	
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mo	ode or Level 1, Function B or C; Key	Radio During Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place	it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form	120-G) at the Intoxilyzer® Location	at the Agency.
10. Send the Following to the Crime Laboratory: Complete 104-G), Print Test, ACA Test, and RFI Test.	ed Intoxilyzer® 8000 Installation and	d Repair Checkout (SFN59281, Form
Field Inspector Signature	anterior de la companya de la compa	Date 5/7/2021
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved the date the Field Inspector performed the installation. This record county of Burleigh, North Dakota, is certified to be a true and corre	d on file at the Office of Attorney Ge	neral, Crime Laboratory Division, in the
Reviewed/Certified By	Lo	Certified Date
(Maria Ce	Ed .	5.7.2021

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005957
Location = MAND 8164.14.00 09/16
05/07/2021 09:39

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Start Test Sequence: Enabled DABACABA Display Prelim Rslt? Yes Display Third Digit?
Inhib Printer(Y/N)? Yes No No Display Volume? Yes Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 056200 Standard Lot #? 05620080A1 Standard Cyl #? 39 Standard Expiration? 04/05/2022 05620080A1 020611 Oper No?

Flow Cal. Date: 06/21/2013 Slope 647 Intercept -613317

IR Calibration Date: 05/11/2015
3um 9um

 0th Coef(*100):
 -16666
 -18045

 1st Coef(*100):
 278478
 132949

 2nd Coef(*100):
 3303
 1397

 H2O adj(mg/l*10k):
 379
 351

******* Printer Test End *******

Operator Signature
TOM SHARP

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MAND
05/07/2021

SN 80-005957 8164.14.00 09/16 09:39

DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.080	09:40 09:40 09:41 09:41 09:42 09:42
0 / 1000111 1177		

Lot No = 05620080A1

Cyl No = 39

Exp Date = 04/05/2022

County = 30

Oper No. = 020611

Operator Signature
TOM SHARP

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005957 Location = MAND 05/07/2021

8164.14.00 09/16 09:43

Test	AC	Time
01 Room Air	RFI*	09:45
02 Room Air	RFI*	09:45

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 01/01/2000

Weight = 150 Sub Sex = Male Cit = 00000Test = DUI

Dr. Lic. = ND/TES000000

Lot No = 05620080A1

Cyl No = 39

Expiration Date = 04/05/2022

Oper No. = 020611 County = 30

> perator Signature TOM SHARP

Remarks:

Form 106-I8000