

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location (ass County 80)	Tail
Reason for Install/Repair	
Install After Receiving From Crime Laboratory 🔲 Install After Location Change	
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
☐ 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
X 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
🔀 C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	uring Test).
√7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
2 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form
Field Inspector Signature	5-13-2621
Crime Laboratory Use Only	•
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to deter the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crin County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	mine alcohol concentration from ne Laboratory Division, in the
Reviewed/Certified By	Certified Date 5 21 21

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CASS
05/13/2021

Alcohol Analyzer
SN 80-005944
8164.14.00 09/16
07:34

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Enabled Data Entry Mode: DABACABA Start Test Sequence: Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No No Display Volume? Yes Disable On Memfull? 1 # of Print Copies? Dry Select Std (D/W/I)? 0.080 Standard Value? 05620080A1 Standard Lot #? Standard Cyl #? 26 Standard Expiration? 04/05/2022 090801 Oper No?

Flow Cal. Date: 06/20/2013
Slope 677
Intercept -702620

******* Printer Test End ********

Operator Signature KATHRYN ALLEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = CASS 05/13/2021 Alcohol Analyzer SN 80-005944 8164.14.00 09/16 07:29

DRY CAL CHECK

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std. 07 Room	Gas Air Gas Air Gas	0.000 0.083 0.000 0.082 0.000 0.083 0.000	07:30 07:30 07:30 07:31 07:31 07:32

Lot No = 05620080A1

Cyl No = 26

Exp Date = 04/05/2022

County = 09

Oper No. = 090801

pperator Signature KATHRYN ALLEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

Alcohol Analyzer CMI, Inc. Intoxilyzer SN 80-005944 North Dakota Model 8000 8164.14.00 09/16 Location = CASS 07:36 05/13/2021

Test	AC	Time
01 Room Air	0.000	07:38
02 *Subject Test	RFI*	07:38
03 Room Air	0.000	07:38

^{*}Invalid Test Inhibited - RFI

Sub Name = TEST SUBJECT, TEST SUBJECT NA Sub DOB = 01/01/1901Weight = 200Sub Sex = Female Cit = 100000000000000 Test = DUI Dr. Lic. = ND/NALot No = 05620080A1Cyl No = 26Expiration Date = 04/05/2022Oper No. = 090801

> Operator Signature KATHRYN ALLEN

Remarks:

County = 09

Form 106-I8000