

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-005362

**A. GENERAL Setup and Checks:**

1.  Power On  
 5 min. countdown during initial warm-up  
 1 min. countdown to begin test after standby mode
2.  Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:  
 Level 1 (Operator)     Level 2 (Field Inspector)     Level 3 (Lab Personnel)
4.  Printer Setup (Level 3,S,P)
5.  Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

1.  Internal Printer Test and Paper Feed (Level 1,P). Attach record.
2.  External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):  
 Operator ID     ND License     R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

1.  Continuous Air Blank (Level 2,A)
2.  Diagnostic (Level 2,D). Attach record.
3.  Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):  
Inst Ser#: 80-005362    DSP Ver: 1362.14  
Core Ver: 8164.14    CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

1. DVM Monitor (Level 3,D,D): 3: 12248 s: 1 c: 47.00  
9: 12872 s: 2 B: 45 0  
(RFI indicator with keyed radio): 1
2. Function (F) Keys: :  F5 – Tone                                     F9 - Fan  
 F6 – Green LED                                     F10 – Pump  
 F7 – Red LED                                         F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 948 CF: 1.06
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 132 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 49.0 I: 0.0  
B: 45.1 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 296 psi P (gauge): 300 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C):  
Sim ser #: MP5320 Lot #: MILLI @ H<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3070 Lot #: 201906A Exp date: 6-11-21  
Known value: 0.010 AC  
Results: 0.008 AC 0.008 AC 0.008 AC
3. High AC (Level 1,C):  
Sim ser #: MP3061 Lot #: 20101 Exp date: 2-13-22  
Known value: 0.500 AC  
Results: 0.495 AC 0.493 AC 0.492 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR3847 Lot #: ICS6 Exp date: NA  
Known value: 0.10% AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5188 Lot #: ACE1 Exp date: NA  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):  
Sim ser #: DR5184 Lot #: MET3 Exp date: NA  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA  
Known value: 0.10% ISOPROPNOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):


Ave: 0.0800 Std Dev: 0.0006 Rel Std Dev: 0.82%

Maintenance/Repairs/Notes: INITIAL INSPECTION DUE TO OPENING  
UP INSTRUMENT TO REPLACE <sup>CELL</sup> REAL TIME  
CLOCK BATTERY.

Instrument is acceptable to be used in the field.

  
Field Inspector's Signature

5/20/21  
Date

  
Reviewed by

5/21/21  
Date

Intoxilyzer Test Record and Checklist  
NDODAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-005362  
Location = TOXL 8164.14.00 09/16  
05/19/2021 16:39

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMNPO

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNOPQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNOPQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNOPQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup


Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 24119080A1  
Standard Cyl #? 31  
Standard Expiration? 11/05/2021  
Oper No? 666666

Flow Cal. Date: 05/22/2020  
Slope 659  
Intercept -620990

IR Calibration Date: 05/26/2020  
3um 9um

|                    |        |        |
|--------------------|--------|--------|
| 0th Coef(*100):    | -20883 | -25370 |
| 1st Coef(*100):    | 265508 | 144960 |
| 2nd Coef(*100):    | 936    | 714    |
| H2O adj(ng/l*10K): | 391    | 379    |

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: **INTERNAL  
PRINT  
TEST**  
Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005362  
Location = TOXL                 8164.14.00 09/16  
05/19/2021                        16:42

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?


Current Instrument Setup

Data Entry Mode:                Enabled  
Start Test Sequence:            DABACABA  
Display Prelim Rslt?            Yes  
Display Third Digit?            Yes  
Inhib Printer(Y/N)?            No  
Display Volume?                No  
Disable On Memfull?            Yes  
# of Print Copies?              1  
Select Std (D/W/I)?            Dry  
Standard Value?                0.080  
Standard Lot #?                 24119080A1  
Standard Cyl #?                 31  
Standard Expiration?            11/05/2021  
Oper No?                        666666  
  
Flow Cal. Date:                 05/22/2020  
Slope                            659  
Intercept                        -620990

IR Calibration Date:            05/26/2020  
                                    3um                    9um

-----  
0th Coef(\*100):                -20883                -25370  
1st Coef(\*100):                265508                144960  
2nd Coef(\*100):                936                    714  
H2O adj(mg/l\*10k):            391                    379

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

EXTERNAL PRINT  
TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005362  
Location = TOXL                    8164.14.00 09/16  
05/19/2021                            16:44

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 16:44 |
| 02 *Subject Test | REF*  | 16:45 |
| 03 Room Air      | 0.000 | 16:45 |

\*Subject Test Refused

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                            Weight = NA  
Test = DUI                            Cit = R KEY REFUSAL  
Dr. Lic. = ND/TES851012  
Lot No = 24119080A1  
Cyl No = 31  
Expiration Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: "R" KEY REFUSAL

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/19/2021      16:46

DIAGNOSTIC

|                           |      |
|---------------------------|------|
| Voltage/Current Test      | Pass |
| RAM Test                  | Pass |
| EEPROM Checksum Test      | Pass |
| Real Time Clock Test      | Pass |
| DSP Test                  | Pass |
| Analytical Stability Test | Pass |
| Modem Test                | Pass |
| Temp Regulation Test      | Pass |

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *DIAGNOSTIC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/19/2021      16:55

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 16:56 |
| 02 Std. Sol. | 0.000 | 16:57 |
| 03 Room Air  | 0.000 | 16:57 |
| 04 Std. Sol. | 0.000 | 16:58 |
| 05 Room Air  | 0.000 | 16:58 |
| 06 Std. Sol. | 0.000 | 16:59 |
| 07 Room Air  | 0.000 | 17:00 |

08 Sim Temp = 34.0°C

Simul Ser No = MP5320  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      ZERO AC  
                    0.000 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      09:32

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 09:32 |
| 02 Std. Sol. | 0.008 | 09:33 |
| 03 Room Air  | 0.000 | 09:33 |
| 04 Std. Sol. | 0.008 | 09:34 |
| 05 Room Air  | 0.000 | 09:35 |
| 06 Std. Sol. | 0.008 | 09:35 |
| 07 Room Air  | 0.000 | 09:36 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3070  
Std Sol No = 201906A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      0.010 AC  
                    Low AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      09:43

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 09:44 |
| 02 Std. Sol. | 0.495 | 09:45 |
| 03 Room Air  | 0.000 | 09:45 |
| 04 Std. Sol. | 0.493 | 09:46 |
| 05 Room Air  | 0.000 | 09:46 |
| 06 Std. Sol. | 0.492 | 09:47 |
| 07 Room Air  | 0.000 | 09:48 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 20101  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.500 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005362  
Location = TOXL                 8164.14.00 09/16  
05/20/2021                        09:59

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 09:59 |
| 02 *Subject Test | INT*  | 10:00 |
| 03 Room Air      | 0.000 | 10:00 |

\*Invalid Test  
Interferent Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                     Weight = NA  
Test = DUI                         Cit = INTERFERENT  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005362  
Location = TOXL                 8164.14.00 09/16  
05/20/2021                        10:01

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 10:02 |
| 02 *Subject Test | INT*  | 10:03 |
| 03 Room Air      | 0.000 | 10:03 |

\*Invalid Test  
Interferent Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                     Weight = NA  
Test = DUI                         Cit = INTERFERENT  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      10:04

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 10:05 |
| 02 *Subject Test | INT*  | 10:05 |
| 03 Room Air      | 0.000 | 10:06 |

\*Invalid Test  
Interferent Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*

*0.10% METHANOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      10:07

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 10:08 |
| 02 *Subject Test | INT*  | 10:09 |
| 03 Room Air      | 0.000 | 10:09 |

\*Invalid Test  
Interferent Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10% ISOPROPAVOL*  
Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      10:10

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | AMB*  | 10:11 |
| 02 Room Air | 0.000 | 10:12 |

\*Invalid Test  
Check Ambient Conditions

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = AMBIENT FAIL  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *AMBIENT FAIL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005362  
Location = TOXL                    8164.14.00 09/16  
05/20/2021                            10:13

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Room Air       | 0.000 | 10:13 |
| 02 Invalid Sample | X.XXX | 10:14 |
| 03 Room Air       | 0.000 | 10:14 |

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                            Weight = NA  
Test = DUI                                Cit = INVALID SAMPLE  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                                Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INVALID SAMPLE*

Form 106-I8000






Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005362  
Location = TOXL                    8164.14.00 09/16  
05/20/2021                            10:17

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Room Air       | 0.000 | 10:18 |
| 02 Subject Test 1 | 0.146 | 10:18 |
| 03 Room Air       | 0.000 | 10:19 |
| 04 Reported AC    | 0.146 | 10:18 |

No RFI Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                            Weight = NA  
Test = DUI                                    Cit = ABA MODE  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                                    Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:    *ABA MODE TEST*  
*0.150 AC*

Form 106-I8000

*FOR BREATH SAMPLE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005362  
Location = TOXL                 8164.14.00 09/16  
05/20/2021                        10:20

| Test              | AC     | Time  |
|-------------------|--------|-------|
| 01 Diagnostic     | OK     | 10:21 |
| 02 Room Air       | 0.000  | 10:21 |
| 03 Subject Test 1 | 0.195  | 10:22 |
| 04 Room Air       | 0.000  | 10:24 |
| 05 Std. Gas       | 0.096* | 10:25 |
| 06 Room Air       | 0.000  | 10:26 |

\*Cal Check Out of Tolerance

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                     Weight = NA  
Test = DUI                         Cit = CAL CHK OUT TOL  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                        Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
OUT OF TOLERANCE

Form 106-I8000

0.200 AC FOR  
BREATH SAMPLE

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      10:40

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Diagnostic     | OK    | 10:41 |
| 02 Room Air       | 0.000 | 10:41 |
| 03 Subject Test 1 | 0.192 | 10:42 |
| 04 Room Air       | 0.000 | 10:44 |
| 05 Std. Gas       | 0.083 | 10:45 |
| 06 Room Air       | 0.000 | 10:46 |
| 07 *Subject Test  | RFI*  | 10:47 |
| 08 Room Air       | 0.000 | 10:47 |

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = RFI DETECTED  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-I8000

*EXPIRED 0.200 AC SIM. SOL.  
USED FOR BREATH SAMPLE ONE\*  
SIM. SOL. SITTING IN SIMULATOR FOR  
OVER 6 MONTHS  
\*SUBJECT TEST 1*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005362  
Location = TOXL                    8164.14.00 09/16  
05/20/2021                            10:49

| Test             | AC     | Time  |
|------------------|--------|-------|
| 01 Diagnostic    | OK     | 10:50 |
| 02 Room Air      | 0.000  | 10:51 |
| 03 *Subject Test | 0.000* | 10:54 |
| 04 Room Air      | 0.000  | 10:56 |
| 05 Std. Gas      | 0.081  | 10:57 |
| 06 Room Air      | 0.000  | 10:58 |
| 07 *Subject Test | 0.000* | 11:01 |
| 08 Room Air      | 0.000  | 11:02 |

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male                            Weight = NA  
Test = DUI                            Cit = DEFICIENT SAMPL  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: DEFICIENT SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      11:03

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Diagnostic     | OK    | 11:04 |
| 02 Room Air       | 0.000 | 11:04 |
| 03 Subject Test 1 | 0.141 | 11:05 |
| 04 Room Air       | 0.000 | 11:08 |
| 05 Std. Gas       | 0.082 | 11:09 |
| 06 Room Air       | 0.000 | 11:10 |
| 07 Reference      | IPS*  | 11:10 |
| 08 Room Air       | 0.000 | 11:11 |

\*Improper Sample

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = IMPROPER SAMPLE  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *IMPROPER SAMPLE*

Form 106-I8000

*EXPIRED 0.150 AC SIM. SOL. USED  
FOR SUBJECT TEST 1e SIM. SOL. SITTING  
IN SIMULATOR FOR OVER 6 MONTHS.*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      11:13

| Test              | AC     | Time  |
|-------------------|--------|-------|
| 01 Diagnostic     | OK     | 11:14 |
| 02 Room Air       | 0.000  | 11:14 |
| 03 Subject Test 1 | 0.194* | 11:15 |
| 04 Room Air       | 0.000  | 11:17 |
| 05 Std. Gas       | 0.082  | 11:18 |
| 06 Room Air       | 0.000  | 11:19 |
| 07 Subject Test 2 | 0.142* | 11:20 |
| 08 Room Air       | 0.000  | 11:21 |

\*Difference Too Great

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = DIFF TOO GREAT  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *DIFFERENCE TOO GREAT*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      11:40


| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Diagnostic     | OK    | 11:41 |
| 02 Room Air       | 0.000 | 11:41 |
| 03 Subject Test 1 | 0.000 | 11:42 |
| 04 Room Air       | 0.000 | 11:44 |
| 05 Std. Gas       | 0.080 | 11:45 |
| 06 Room Air       | 0.000 | 11:46 |
| 07 Subject Test 2 | 0.000 | 11:47 |
| 08 Room Air       | 0.000 | 11:47 |
| 09 Reported AC    | 0.000 | 11:42 |

Difference OK  
No RFI Detected

Sub Name = TEST, BOBBY NONE  
Sub DOB = 01/01/1985  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = CMS MODE TEST  
Dr. Lic. = ND/TES851012  
Lot No = 02621080A1  
Cyl No = 25  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*CMS MODE  
TEST*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = TOXL      8164.14.00 09/16  
05/20/2021      11:52

| Test                | AC     | Time  |
|---------------------|--------|-------|
| 01 Room Air         | 0.000  | 11:52 |
| 02 Std. Gas         | 0.080  | 11:53 |
| 03 Room Air         | 0.000  | 11:53 |
| 04 Std. Gas         | 0.080  | 11:53 |
| 05 Room Air         | 0.000  | 11:54 |
| 06 Std. Gas         | 0.081  | 11:54 |
| 07 Room Air         | 0.000  | 11:55 |
| 08 Std. Gas         | 0.081  | 11:55 |
| 09 Room Air         | 0.000  | 11:56 |
| 10 Std. Gas         | 0.080  | 11:56 |
| 11 Room Air         | 0.000  | 11:57 |
| 12 Std. Gas         | 0.079  | 11:57 |
| 13 Room Air         | 0.000  | 11:57 |
| 14 Std. Gas         | 0.081  | 11:58 |
| 15 Room Air         | 0.000  | 11:58 |
| 16 Std. Gas         | 0.080  | 11:59 |
| 17 Room Air         | 0.000  | 11:59 |
| 18 Std. Gas         | 0.080  | 12:00 |
| 19 Room Air         | 0.000  | 12:00 |
| 20 Std. Gas         | 0.081  | 12:00 |
| 21 Room Air         | 0.000  | 12:01 |
| 22 Std. Gas         | 0.080  | 12:01 |
| 23 Room Air         | 0.000  | 12:02 |
| 24 Std. Gas         | 0.079  | 12:02 |
| 25 Room Air         | 0.000  | 12:03 |
| 26 Std. Gas         | 0.080  | 12:03 |
| 27 Room Air         | 0.000  | 12:04 |
| 28 Std. Gas         | 0.080  | 12:04 |
| 29 Room Air         | 0.000  | 12:04 |
| 30 Std. Gas         | 0.079  | 12:05 |
| 31 Room Air         | 0.000  | 12:05 |
| 32 Std. Gas         | 0.080  | 12:06 |
| 33 Room Air         | 0.000  | 12:06 |
| 34 Std. Gas         | 0.080  | 12:07 |
| 35 Room Air         | 0.000  | 12:07 |
| 36 Std. Gas         | 0.080  | 12:07 |
| 37 Room Air         | 0.000  | 12:08 |
| 38 Std. Gas         | 0.080  | 12:08 |
| 39 Room Air         | 0.000  | 12:09 |
| 40 Std. Gas         | 0.079  | 12:09 |
| 41 Room Air         | 0.000  | 12:10 |
| 42 Cal Check Stats: |        |       |
| 43 Average          | 0.0800 |       |
| 44 Std Dev          | 0.0006 |       |
| 45 Rel Std Dev      | 0.82%  |       |



Operator Signature  
CHARLES EDER

Remarks: 0.080 AC

STABILITY TEST

Form 106-I8000

LOT: 02621080A1

EXP: 3-5-23

QYL: 25