INTOXILYZER® 8000 REPAIR AND MAINTENANCE

	_
Intoxil	yzer® 8000 Serial Number: <u>80-006496</u> Location: <u>70メレ</u>
1.	Reported Issue: CORE VERSION WAS 8164.13. SHOULD BE UPDATED TO 8164.14.
2.	Found/Determined: INSTRUMENT CORE VERSION WAS AT 8164.13, SHOULD BE UPDATED TO 8164.14. THIS WILL ADD DRIVERS LICENSE FORMATS CONSISTENT WITH UPDATES IN 2014,
3.	Repair/Maintenance Performed: UP DATED CORE VERSION TO 8164.14. USED FLASH UP DATE PROCESS IN CMI'S SOFTWARE. DEMONSTRATED PROCESS TO KALI HIEB ON 5/18/21.
	THETE WHENT NOW READS NO DRIVER'S LICENSES ISSUED BEFORE, DURING AND AFTER 2014.
4.	Check and reset if needed; time, date, and location. (Level 2, E)
5.	Tank monitor display and regulator within 50 psi. (Level 3, D, G) Display: 210 psi Regulator: 225 psi If needed, completed tare of tank sensor. (Level 3, M, C, G)
6.	Does the flow need to be calibrated? If Yes, attach paperwork. Yes of No
7.	Does the optical bench need a calibration adjustment? If Yes, attach test records. Does the instrument require further testing? Yes or No
8.	Does the instrument require further testing? If No, sign and date at end of form. If Yes, state testing required or: Deform all tests for repair. WILL USE DRIVEL'S LICENSES BEFORE AND DURING 2014. THIS WILL DEMONSTRATE READING BARCODES ON LICENSES CORECTLY. ABA INTERFERENT - PRE 2014 LICENSE. (ISSUED 2006) CMS RFI TEST - 2014 ISSUED LICENSE. (ISSUED 2014)
	CHOINT ILST USIN CHENSE. (15) MEDICOM

03/20 Issuing Authority - ARQ Page 1 of 2

OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

BrW-005

Configure simulator(s) for the following tests: (Level 1, S) Attach test records			
9. Wet Calibration Check: (ACA Test)			
X Low AC (Level 1, C): Known value ≤ 0.03 AC. 0.010 AC Sim SN: MP3070 Lot #: 201906 A Exp. Date:			
X High AC (Level 1, C): Known value ≥ 0.25 AC.			
Wet Calibration Check results are within \pm 0.005 or \pm 5% (whichever is greater) of stated AC value.			
Configure dry gas standard for the following tests: (Level 1, S) Attach test records.			
10. Print test. (Level 1, P)			
11. Interferent Check. O. O. AC ethanol plus 0.05% Acetone. (Level 1, B) Sim SN: DR3847 Lot #: <u>ICS 6</u> Display reads "Interferent Detect"			
12.RFI Check. (CMS mode) Display reads "RFI Detect"			
13. Dry Calibration Check/Stability Test (Level 3, T, S, D, 15) Known: 0.080 AC Ethanol Gas Standard. Lot No. 02621080 A1 Cylinder No. 25 Expiration Date 3.5.23			
Average AC: 0.0806 Std. Dev. 0.005 Rel. Std. Dev. 0.63%			
Remarks/Notes: APDED CMS TEST - LICENSE ISSUED IN 2019. - 2019 LICENSE IS A REAL FD.			
Instrument is acceptable to be used in the field.			
Malla E El 5/20/21			
Breath Analyst Signature Date			
Reviewed by $\frac{2/24/21}{\text{Date}}$			

03/20 Issuing Authority - ARQ Page 2 of 2

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/20/2021 Alcohol Analyzer SN 80-006496 8164.14.00 09/16 14:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:41
02 Std. Sol.	0.009	14:42
03 Room Air	0.000	14:42
04 Std. Sol.	0.009	14:43
05 Room Air	0.000	14:43
06 Std. Sol.	0.009	14:44
07 Room Air	0.000	14:45

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3070Std Sol No = 201906A

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low AC

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/20/2021 SN 80-006496 8164.14.00 09/16 14:50

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Sol.	0.500	14:52
03 Room Air	0.000	14:52
04 Std. Sol.	0.500	14:53
05 Room Air	0.000	14:54
06 Std. Sol.	0.499	14:55
07 Room Air	0.000	14:55

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3061 Std Sol No = 20101

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.500 AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006496
Location = TOXL 8164.14.00 09/16
05/20/2021 14:59

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 25 Standard Expiration? 03/05/2023 Oper No? 666666

Flow Cal. Date: 01/29/2020

Slope 663 Intercept -706121

IR Calibration Date: 08/18/2015

Oth Coef(*100): -22081 -15963
1st Coef(*100): 278490 136818
2nd Coef(*100): 2086 1282
H2O adj(mg/l*10k): 598 386

3um

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/20/2021

Alcohol Analyzer SN 80-006496 8164.14.00 09/16 15:19

Т	est	AC	Time
	Room Air	0.000	15:20
	Subject Test Room Air	INT 0.000	15:21 15:22

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Weight = 150 Sub Sex = Female Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 02621080A1

Cyl No = 25

Expiration Date = 03/05/2023

County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER

Front Check + 0.05% ACE HOWE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006496
Location = TOXL 8164.14.00 09/16
05/20/2021 15:23

Test	AC	Time
01 Diagnostic	OK	15:24
02 Room Air	0.000	15:24
03 *Subject Test	RFI*	15:24
04 Room Air	0.000	15:25

*Invalid Test Inhibited - RFI

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/TES989643

Lot No = 02621080A1

Cyl No = 25

Expiration Date = 03/05/2023

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/20/2021

Alcohol Analyzer SN 80-006496 8164.14.00 09/16 15:36

Test	AC	Time
01 Room Air	0.000	15:37
02 Std. Gas	0.081	15:37
03 Room Air	0.000	15:38
04 Std. Gas	0.080	15:38
05 Room Air	0.000	15:39
06 Std. Gas	0.081	15:39
07 Room Air	0.000	15:40
08 Std. Gas	0.081	15:40
09 Room Air	0.000	15:41
10 Std. Gas	0.080	15:41
11 Room Air	0.000	15:42
12 Std. Gas	0.081	15:42
13 Room Air	0.000	15:43
14 Std. Gas	0.081	15:43
15 Room Air	0.000	15:44
16 Std. Gas	0.081	15:44
17 Room Air	0.000	15:45
18 Std. Gas	0.080	15:45
19 Room Air	0.000	15:45
20 Std. Gas	0.080	15:46
21 Room Air	0.000	15:46
22 Std. Gas	0.081	15:47
23 Room Air	0.000	15:47
24 Std. Gas	0.080	15:48
25 Room Air	0.000	15:48
26 Std. Gas	0.080	15:49
27 Room Air	0.000	15:49
28 Std. Gas	0.081	15:50
29 Room Air	0.000	15:50
30 Std. Gas	0.081	15:51
31 Room Air	0.000	15:51
32 Cal Check Stats:		
33 Average	0.0806	
34 Std Dev	0.0005	
35 Rel Std Dev	0.63%	

Operator Signature CHARLES EDER

PALIBRATION CHECK, STABILITY TEST 0.080 AC

15:26

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006496 Location = TOXL 8164.14.00 09/16 05/20/2021

Test	AC	Time
01 Diagnostic	OK	15:28
02 Room Air	0.000	15:29
03 Subject Test 1	0.000	15:29
04 Room Air	0.000	15:31
05 Std. Gas	0.080	15:32
06 Room Air	0.000	15:33
07 Subject Test 2	0.000	15:34
08 Room Air	0.000	15:35
09 Reported AC	0.000	15:29

Difference OK No RFI Detected

Sub Name = REAL ID, REAL ID REAL ID

Sub DOB = 10/18/1920

Sub Sex = MaleWeight = NA

Test = DUI Cit = REAL ID LICENSE

Dr. Lic. = ND/XXX201234

Lot No = 02621080A1

Cyl No = 25

Expiration Date = 03/05/2023

County = 08Oper No. = 666666

20 minute waiting period ascertained? Y Clean Mouthpiece used and disposed? Y Clean Mouthpiece used and disposed?

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

> Operator Signature CHARLES EDER

REAL ID LIVENSE

CMS TEST

Remarks: