

**INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-006493

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 220 psi Regulator: 225 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3064 Lot #: 20070 Exp. Date: 2-13-22
  - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: MP5317 Lot #: 202012A Exp. Date: 12-8-22

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3847 Lot #: IC56 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 02621080A1 Cyl #: 26 Exp. Date: 3-5-23

Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.082</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

*Chandra E. Eder*  
Breath Analyst Signature

5/25/21  
Date

*Kali G. Dier*  
Reviewed by

5/26/21  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer              Alcohol Analyzer  
North Dakota Model 8000              SN 80-006493  
Location = TOXL                      8164.14.00 09/16  
05/25/2021                              14:49

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                      Enabled  
Start Test Sequence:                  DABACABA  
Display Prelim Rslt?                  Yes  
Display Third Digit?                  Yes  
Inhib Printer(Y/N)?                  No  
Display Volume?                      No  
Disable On Memfull?                  Yes  
# of Print Copies?                    1  
Select Std (D/W/I)?                  Dry  
Standard Value?                      0.080  
Standard Lot #?                      05620080A1  
Standard Cyl #?                      43  
Standard Expiration?                  04/05/2022  
Oper No?                              666666

Flow Cal. Date:                      05/23/2018  
Slope                                  640  
Intercept                              -624033

IR Calibration Date:                  05/23/2018  
    3um                                  9um

-----  
0th Coef(\*100):                      -19857                                  -14695  
1st Coef(\*100):                      263752                                  135087  
2nd Coef(\*100):                      3253                                      1458  
H2O adj(mg/l\*10k): 643                                  483

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006493  
Location = TOXL      8164.14.00 09/16  
05/25/2021      15:09

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Sol.	0.019	15:11
03 Room Air	0.000	15:11
04 Std. Sol.	0.019	15:12
05 Room Air	0.000	15:12
06 Std. Sol.	0.019	15:13
07 Room Air	0.000	15:14

08 Sim Temp = 34.0°C

Simul Ser No = MP3064

Std Sol No = 20070

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006493  
Location = TOXL      8164.14.00 09/16  
05/25/2021      15:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:16
02 Std. Sol.	0.295	15:17
03 Room Air	0.000	15:18
04 Std. Sol.	0.296	15:18
05 Room Air	0.000	15:19
06 Std. Sol.	0.296	15:20
07 Room Air	0.000	15:20

08 Sim Temp = 34.0°C

Simul Ser No = MP5317  
Std Sol No = 202012A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGH AC*  
*0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006493  
Location = TOXL      8164.14.00 09/16  
05/25/2021      15:22

Test	AC	Time
01 Room Air	0.000	15:23
02 *Subject Test	INT*	15:23
03 Room Air	0.000	15:24

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 02621080A1  
Cyl No = 26  
Expiration Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006493  
Location = TOXL                8164.14.00 09/16  
05/25/2021                      15:25

Test	AC	Time
01 Diagnostic	OK	15:26
02 Room Air	0.000	15:27
03 *Subject Test	RFI*	15:27
04 Room Air	0.000	15:27

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = DUI                        Cit = RFI CHECK  
Dr. Lic. = ND/TES989643  
Lot No = 02621080A1  
Cyl No = 26  
Expiration Date = 03/05/2023  
County = 08                      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006493  
Location = TOXL      8164.14.00 09/16  
05/25/2021      15:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:28
02 Std. Gas	0.081	15:28
03 Room Air	0.000	15:29
04 Std. Gas	0.081	15:29
05 Room Air	0.000	15:30
06 Std. Gas	0.081	15:30
07 Room Air	0.000	15:31

Lot No = 02621080A1  
Cyl No = 26  
Exp Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006493  
Location = TOXL                    8164.14.00 09/16  
05/25/2021                            15:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:32
02 Std. Gas	0.081	15:33
03 Room Air	0.000	15:33
04 Std. Gas	0.081	15:34
05 Room Air	0.000	15:34
06 Std. Gas	0.082	15:35
07 Room Air	0.000	15:35

Lot No = 02621080A1  
Cyl No = 26  
Exp Date = 03/05/2023  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006493  
Location = TOXL      8164.14.00 09/16  
05/25/2021      15:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:36
02 Std. Gas	0.081	15:36
03 Room Air	0.000	15:37
04 Std. Gas	0.081	15:37
05 Room Air	0.000	15:38
06 Std. Gas	0.081	15:38
07 Room Air	0.000	15:39

Lot No = 02621080A1  
Cyl No = 26  
Exp Date = 03/05/2023  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000