

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006490	Instrument Location			
Reason for Install/Repair				
Install After Receiving From Crime Laboratory Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
∑ 2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
🐧 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form			
Field Jospector Signature	05/28/2021			
Crime Laboratory Use Only	•			
This installation has been reviewed and the instrument is approved to I	be used for the analysis of breath to determine alcohol concentration from the at the Office of Attorney General, Crime Laboratory Division, in the copy of the documents received.			

FORM 104-G

Reviewed/Certified By

Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006490
Location = ASHL 8164.14.00 09/16
05/28/2021 08:52

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies.
Select Std (D/W/I)? Dry
0.080 # of Print Copies? Standard Value?
Standard Lot #?
Standard Cyl #? 24119080A1 28 Standard Expiration? 11/05/2021 Oper No? 132584

Flow Cal. Date: 09/04/2020 Slope 649 Intercept -598817

IR Calibration Date: 09/04/2020

3um 9um

Oth Coef(*100): -25190 -16392
1st Coef(*100): 271645 136906
2nd Coef(*100): 2533 1222
H2O adj(mg/l*10k): 460 309

******* Printer Test End *******

perator Signature DYLAN THOMES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ASHL
05/28/2021

Alcohol Analyzer SN 80-006490 8164.14.00 09/16 08:58

DRY CAL CHECK

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std. 07 Room	Gas Air Gas Air Gas	0.000 0.079 0.000 0.079 0.000 0.079	08:58 08:59 08:59 08:59 09:00 09:00

Lot No = 24119080A1

Cyl No = 28

Exp Date = 11/05/2021

County = 26

Oper No. = 132584

Operator Signature DYLAN THOMES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ASHL
05/28/2021

Alcohol Analyzer SN 80-006490 8164.14.00 09/16 09:02

DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Room Air	0.000 0.080 0.000 0.079 RFI* 0.000	09:03 09:03 09:03 09:04 09:04

*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 28

Exp Date = 11/05/2021

County = 26

Oper No. = 132584

Operator Signature DYLAN THOMES

Remarks:

Form 106-I8000