

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005937	Instrument Location CAVA - Cavalier, ND	
Reason for Install/Repair		
☑ Install After Receiving From Crime Laboratory ☐ In:	nstall After Location Change	
Other (Specify)		
Check When Done:		
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
☑ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).		
✓ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
☑ 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
☑ C. Radio Frequency Interference (RFI) Test (CMS III)	Mode or Level 1, Function B or C; Key Radio During Test).	
7. Repair and/or Maintenance Performed (if any): Not	ot Applicable	
8. Complete the Top Portion of the Intoxilyzer® Record	(SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.	
9. File Previous Intoxilyzer® Record (SFN504096, Form	m 120-G) at the Intoxilyzer® Location at the Agency.	
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	leted Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form	
Field Inspector Signature	7502 Date 05/29/2021	
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approv	ved to be used for the analysis of breath to determine alcohol concentration fr	
the date the Field Inspector performed the installation. This reco County of Burleigh, North Dakota, is certified to be a true and coi	ord on file at the Office of Attorney General, Crime Laboratory Division, in the prect copy of the documents received.	

FORM 104-G

Reviewed/Certified By

Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005937
Location = CAVA	8164.14.00 09/16
05/29/2021	16:18

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	05620080A1
Standard Cyl #?	2
Standard Expiration?	04/05/2022
Oper No?	060849

Flow Cal. Date:	11/17/2017
Slope	606
Intercept	-534563

IR	Calibration	Date: 3um	11/17/2017 9um
	n Coef(*100): : Coef(*100):		-23772 142357

1st Coef(\*100): 294336 1423 2nd Coef(\*100): -906 850 H2O adj(mg/l\*10k): 205 250

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature DANIEL C. GARVENS

Remarks: Printer Tess

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CAVA
05/29/2021

Alcohol Analyzer SN 80-005937 8164.14.00 09/16 16:20

## DRY CAL CHECK

_	est		AC	Time
	Room		0.000	16:21
02	Std.	Gas	0.079	16:21
03	Room	Air	0.000	16:22
04	std.	Gas	0.080	16:22
05	Room	Air	0.000	16:23
06	Std.	Gas	0.079	16:23
07	Room	Air	0.000	16:24

Lot No = 05620080A1

Cyl No = 2

Exp Date = 04/05/2022

County = 34

Oper No. = 060849

Operator Signature DANIEL C. GARVENS

Remarks: ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CAVA
05/29/2021

Alcohol Analyzer SN 80-005937 8164.14.00 09/16 16:24

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:25
02 Std. Gas	0.080	16:25
03 Room Air	RFI*	16:25
04 Room Air	0.000	16:26

\*Invalid Test Inhibited - RFI

Lot No = 05620080A1

Cyl No = 2

Exp Date = 04/05/2022

County = 34

Oper No. = 060849

Operator Signature DANIEL C. GARVENS

Remarks: RFI Test

Form 106-I8000