



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

|  |  |
|--|--|
| Serial Number<br><b>80-006495</b>  | Instrument Location<br><b>Fargo Police PD 3 1230 25<sup>th</sup> St S. Fargo</b> |
| Reason for Install/Repair  |  |
| <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input type="checkbox"/> Other (Specify) _____ |  |

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

|                               |                           |
|-------------------------------|---------------------------|
| Field Inspector Signature<br> | Date<br><b>05/30/2021</b> |
|-------------------------------|---------------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

|                           |                                    |
|---------------------------|------------------------------------|
| Reviewed/Certified By<br> | Certified Date<br><b>7.22.2021</b> |
|---------------------------|------------------------------------|

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006495  
Location = FAPD                8164.14.00 09/16  
05/30/2021                      22:43

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            02621080A1  
Standard Cyl #?            13  
Standard Expiration?       03/05/2023  
Oper No?                    130056

Flow Cal. Date:            08/18/2015  
Slope                      678  
Intercept                  -657244

IR Calibration Date:        08/18/2015  
                                    3um                      9um

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -22405 | -18664 |
| 1st Coef(*100):    | 270469 | 136654 |
| 2nd Coef(*100):    | 3040   | 1360   |
| H2O adj(mg/l*10k): | 699    | 443    |

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
PATRICK THOMAS

Remarks: *Print test for installation*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006495  
Location = FAPD      8164.14.00 09/16  
05/30/2021      22:45

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 22:46 |
| 02 Std. Gas | 0.081 | 22:46 |
| 03 Room Air | 0.000 | 22:47 |
| 04 Std. Gas | 0.082 | 22:47 |
| 05 Room Air | 0.000 | 22:48 |
| 06 Std. Gas | 0.082 | 22:48 |
| 07 Room Air | 0.000 | 22:49 |

Lot No = 02621080A1  
Cyl No = 13  
Exp Date = 03/05/2023  
County = 09      Oper No. = 130056

  
Operator Signature  
PATRICK THOMAS

Remarks: *ACA Calibration Test for Installation.*

Form 106-18000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006495  
Location = FAPD      8164.14.00 09/16  
05/30/2021      22:50

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 22:52 |
| 02 *Subject Test | RFI*  | 22:52 |
| 03 Room Air      | 0.000 | 22:53 |

\*Invalid Test  
Inhibited - RFI

Sub Name = RFI TEST, NA NA  
Sub DOB = 01/01/2021  
Sub Sex = Unknown      Weight = NA  
Test = DUI      Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 02621080A1  
Cyl No = 13  
Expiration Date = 03/05/2023  
County = 09      Oper No. = 130056



Operator Signature  
PATRICK THOMAS

Remarks: RFI test (detected) for Installation.

Form 106-I8000