

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-004191

**A. GENERAL Setup and Checks:**

- 1.  Power On
  - 5 min. countdown during initial warm-up
  - 1 min. countdown to begin test after standby mode
- 2.  Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
  - Level 1 (Operator)     Level 2 (Field Inspector)     Level 3 (Lab Personnel)
- 4.  Printer Setup (Level 3,S,P)
- 5.  Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

- 1.  Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- 2.  External Printer Test (Level 1,P). Attach record.
- 3. Barcode Reader (Level 1,B):
  - Operator ID     ND License     R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

- 1.  Continuous Air Blank (Level 2,A)
- 2.  Diagnostic (Level 2,D). Attach record.
- 3.  Prelim Data Entry (Level 2,E)
- 4. Ser Number/Version Info (Level 2,V):
  - Inst Ser#: 80-004191    DSP Ver: 1362.14
  - Core Ver: 8164.14    CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

- 1. DVM Monitor (Level 3,D,D): 3: 13221 s: 0 c: 47.06  
9: 13013 s: 2 B: 44 0  
(RFI indicator with keyed radio): 1
- 2. Function (F) Keys: :  F5 – Tone                                     F9 - Fan  
 F6 – Green LED                                     F10 – Pump  
 F7 – Red LED                                         F11 – Solenoid
- 3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 949 CF: 1.06
- 4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 193 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 49.0 I: 0.0  
B: 45.1 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 198 psi P (gauge): 225 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C): ~~CEE 5320~~  
Sim ser #: MP3070 Lot #: MILLI Q H<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3070 Lot #: 201906A Exp date: 6-11-21  
Known value: 0.010 AC  
Results: 0.010 AC 0.010 AC 0.009 AC
3. High AC (Level 1,C):  
Sim ser #: MP3061 Lot #: 20101 Exp date: 2.13.22  
Known value: 0.500 AC  
Results: 0.498 AC 0.497 AC 0.498 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR5114 Lot #: ICS 6 Exp date: NA  
Known value: 0.10AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5188 Lot #: ACE 1 Exp date: NA  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):  
Sim ser #: DR5186 Lot #: MET 3 Exp date: NA  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR5189 Lot #: IPA #2 Exp date: NA  
Known value: 0.10% ISOPROPANOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0800 Std Dev: 0.0008 Rel Std Dev: 1.00%

Maintenance/Repairs/Notes: INITIAL INSPECTION DUE TO OPENING UP INSTRUMENT TO REPLACE REAL TIME CLOCK BATTERY.

Instrument is acceptable to be used in the field.

Charles E. Ea  
Field Inspector's Signature

6/2/21  
Date

Kali L. Hib  
Reviewed by

6-7-21  
Date

Intoxilyzer Test Record and Checklist  
ND006G Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004191  
Location = TOXL 6154.14.00 09/16  
06/02/2021 13:13

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNOPQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNOPQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNOPQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Reslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/!)? Dry  
Standard Value? 0.000  
Standard Lot #? 1541035191  
Standard Cyl #? 4  
Standard Expiration? 07/01/2012  
Oper No? 666666

Flow Cal. Date: 01/26/2018  
Slope: 736  
Intercept: -777674

IR Calibration Date: 04/20/2017  
Sun Sun  
-----  
8th Coef(\*100): -6756 -1.9235  
1st Coef(\*100): 262162 125903  
2nd Coef(\*100): 3907 1462  
420 Adj(ng)\*100): 423 566

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

Internal  
PRINT  
TEST  
Form 106-18000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      13:18

Test	AC	Time
01 Room Air	0.000	13:20
02 *Subject Test	REF*	13:20
03 Room Air	0.000	13:20

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = R KEY REFUSAL  
Dr. Lic. = ND/DIS921456  
Lot No = 15410080A1  
Cyl No = 4  
Expiration Date = 07/01/2012  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: "R" KEY REFUSAL

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      13:22

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

DIAGNOSTIC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      13:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:38
02 Std. Sol.	0.000	13:39
03 Room Air	0.000	13:40
04 Std. Sol.	0.000	13:40
05 Room Air	0.000	13:41
06 Std. Sol.	0.000	13:42
07 Room Air	0.000	13:42

08 Sim Temp = 34.0°C

Simul Ser No = MP5320  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

ZERO AC  
0.000 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      13:44

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:44
02 Std. Sol.	0.010	13:45
03 Room Air	0.000	13:45
04 Std. Sol.	0.010	13:46
05 Room Air	0.000	13:47
06 Std. Sol.	0.009	13:47
07 Room Air	0.000	13:48

08 Sim Temp = 34.0°C

Simul Ser No = MP3070  
Std Sol No = 201906A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Low AC*  
*0.010 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      13:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:55
02 Std. Sol.	0.498	13:56
03 Room Air	0.000	13:57
04 Std. Sol.	0.497	13:57
05 Room Air	0.000	13:58
06 Std. Sol.	0.498	13:59
07 Room Air	0.000	13:59

08 Sim Temp = 34.0°C

Simul Ser No = MP3061  
Std Sol No = 20101  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.500 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      14:04

Test	AC	Time
01 Room Air	0.000	14:04
02 *Subject Test	INT*	14:05
03 Room Air	0.000	14:05

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:08

Test	AC	Time
01 Room Air	0.000	14:09
02 *Subject Test	INT*	14:09
03 Room Air	0.000	14:10

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = DUI                                Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                                Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.05% ACETONE*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:11

Test	AC	Time
01 Room Air	0.000	14:12
02 *Subject Test	INT*	14:12
03 Room Air	0.000	14:13

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                          Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10% METHANOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      14:14

Test	AC	Time
01 Room Air	0.000	14:15
02 *Subject Test	INT*	14:15
03 Room Air	0.000	14:16

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
~~0.10% METHANOL CEE~~  
*car*      Form 106-I8000      ~~CEE~~

0.10% ISOPROPANOL


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:18

Test	AC	Time
01 Room Air	AMB*	14:19
02 Room Air	0.000	14:20

\*Invalid Test  
Check Ambient Conditions

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = DUI                          Cit = AMBIENT FAIL  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                          Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *AMBIENT FAIL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      14:21

Test	AC	Time
01 Room Air	0.000	14:21
02 Invalid Sample	X.XXX	14:22
03 Room Air	0.000	14:22

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = INVALID SAMPLE  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INVALID SAMPLE*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:23

Test	AC	Time
01 Room Air	0.000	14:24
02 *Subject Test	ABT*	14:24
03 Room Air	0.000	14:25

\*Sequence Aborted

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = DUI                            Cit = SEQUENCE ABORT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: SEQUENCE ABORTED

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      14:25

Test	AC	Time
01 Room Air	0.000	14:25
02 Subject Test 1	0.194	14:27
03 Room Air	0.000	14:28
04 Reported AC	0.194	14:27

No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = ABA MODE  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: ABA MODE  
TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:29

Test	AC	Time
01 Diagnostic	OK	14:30
02 Room Air	0.000	14:30
03 Subject Test 1	0.192	14:30
04 Room Air	0.000	14:33
05 Std. Gas	0.087*	14:34
06 Room Air	0.000	14:35

\*Cal Check Out of Tolerance

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = DUI                                Cit = CAL CHK OUT TOL  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                                Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
OUT OF TOLERANCE

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            14:36

Test	AC	Time
01 Diagnostic	OK	14:37
02 Room Air	0.000	14:38
03 Subject Test 1	0.191	14:38
04 Room Air	0.000	14:41
05 Std. Gas	0.079	14:42
06 Room Air	0.000	14:43
07 *Subject Test	RFI*	14:43
08 Room Air	0.000	14:44

\*Invalid Test  
  Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                          Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:



Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004191  
Location = TOXL      8164.14.00 09/16  
06/02/2021      14:44

Test	AC	Time
01 Diagnostic	OK	14:45
02 Room Air	0.000	14:46
03 *Subject Test	0.000*	14:49
04 Room Air	0.000	14:51
05 Std. Gas	0.079	14:52
06 Room Air	0.000	14:54
07 *Subject Test	0.000*	14:57
08 Room Air	0.000	14:58


No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = DUI      Cit = DEFICIENT SAMP  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

DEFICIENT  
SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                 8164.14.00 09/16  
06/02/2021                        14:59

Test	AC	Time
01 Diagnostic	OK	15:00
02 Room Air	0.000	15:01
03 Reference	IPS*	15:01
04 Room Air	0.000	15:01

\*Improper Sample

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                     Weight = NA  
Test = DUI                         Cit = IMPROPER SAMPLE  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                         Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *IMPROPER SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                    8164.14.00 09/16  
06/02/2021                            15:02

Test	AC	Time
01 Diagnostic	OK	15:03
02 Room Air	0.000	15:03
03 Subject Test 1	0.048*	15:04
04 Room Air	0.000	15:06
05 Std. Gas	0.081	15:07
06 Room Air	0.000	15:09
07 Subject Test 2	0.178*	15:09
08 Room Air	0.000	15:10

\*Difference Too Great

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = DUI                                Cit = DIFF TOO GREAT  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                                Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: **DIFFERENCE  
TOO GREAT**

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004191  
Location = TOXL                 8164.14.00 09/16  
06/02/2021                        15:11


Test	AC	Time
01 Diagnostic	OK	15:12
02 Room Air	0.000	15:13
03 Subject Test 1	0.188	15:13
04 Room Air	0.000	15:16
05 Std. Gas	0.079	15:17
06 Room Air	0.000	15:18
07 Subject Test 2	0.187	15:19
08 Room Air	0.000	15:20
09 Reported AC	0.187	15:19

Difference OK  
No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                     Weight = NA  
Test = DUI                         Cit = CMS MODE  
Dr. Lic. = ND/TES989643  
Lot No = 05620080A1  
Cyl No = 6  
Expiration Date = 04/05/2022  
County = 08                        Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*CMS MODE  
TEST*

Form 106-I8000





Operator Signature  
CHARLES EDER

Remarks: STABILITY TEST  
0.080 AC

Form 106-I8000

LOT: 05620080A1

CYL: 006

EXP: 4.5.22