

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	005946	Instrument Location	FINE	<u>/</u>
Reason for Install/Repair		**************************************		
Install After Receiving Fr	om Crime Laboratory	Install After Location Change		
Other (Specify)	Install af	ter download	at Sta	ate Lab
Check When Done:				
1. Surge Protector Instal	led/Property Grounded.			
2. Telephone Line Conne	ected to Intoxilyzer® 8000.			
3. Breath Tube Heated.				
4. Enter Preliminary Data	a (i.e. Date, Time, DST (Y), a	and Location; Level 2, Function	ı E).	
5. Scan/Enter Gas Cylind	der Information (Level 1, Fun	nction S).		
6. Run Tests:				
A. Print Test (Level1,	Function P).			
B. ACA Test (Level 1	, Function C).			
C. Radio Frequency I	nterference (RFI) Test (CMS	Mode or Level 1, Function B	or C; Key Radio Du	uring Test).
7. Repair and/or Mainten	ance Performed (if any):			
8. Complete the Top Por	tion of the Intoxilyzer® Reco	rd (SFN50496, Form 120-G) a	nd Place it by the I	ntoxilyzer® for Use.
		orm 120-G) at the Intoxilyzer®		
10. Send the Following to 104-G), Print Test, AC	the Crime Laboratory: Comp CA Test, and RFI Test. Rec	oleted Intoxilyzer® 8000 Instal CEIVED FORM MOD TE	ation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature	well Sand	4		Date 6/4/202
Crime Laboratory Use Only		*		
This installation has been reviewed the date the Field Inspector perform County of Burleigh, North Dakota,	med the installation. This red	cord on file at the Office of Atto	rney General, Crin	
Reviewed/Certified By	Cearles	Less .		Certified Date 6.7.202/

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FINL
06/04/2021

Alcohol Analyzer
SN 80-005946
8164.14.00 09/16
20:36

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 24 Standard Expiration? 11/05/2021 Oper No? 130363

Flow Cal. Date: 06/02/2021 Slope 642 Intercept -600927

IR Calibration Date: 06/02/2021 3um 9um

0 + b	Coef(*100):		0011	
			-20165	-16483
1st	Coef(*100):		290371	
			4903/1	141330
ZIIQ	Coef(*100):		953	910
H20	adj (mg/1*10)	- \ -	106	910
1120	adj (mg/ 1 ~ 10)	<):	496	404

******** Printer Test End *******

Operator Signature RICK CZAPIEWSKI

Remarks:

P-Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FINL
06/04/2021

Alcohol Analyzer SN 80-005946 8164.14.00 09/16 20:44

DRY CAL CHECK

AC	Time
0.000 0.080 0.000 0.079 0.000 0.079 0.000	20:45 20:45 20:46 20:46 20:47 20:47
	0.000 0.080 0.000 0.079 0.000 0.079

Lot No = 24119080A1

Cyl No = 24

Exp Date = 11/05/2021

County = 46

Oper No. = 130363

Operator Signature RICK CZAPIEWSKI

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FINL
06/04/2021

Alcohol Analyzer
SN 80-005946
8164.14.00 09/16
20:52

Test	AC	Time
01 Room Air	RFI*	20:54
02 Room Air	0.000	20:54

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 11/11/2000
Sub Sex = Male
Test = DUI
Dr. Lic. = ND/NA
Lot No = 24119080A1
Cyl No = 24
Expiration Date = 11/05/2021

Oper No. = 130363

Operator Signature RICK CZAPIEWSKI

Remarks:

County = 46

RFI

Form 106-I8000