# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

# **INTOXILYZER® 8000 INSPECTION**

Into	oxilyzer® 8000 Serial Number: <u>80-∞9206</u>	Location: ToXL		
A.	General Setup and Checks:  1. Diagnostic Tests Pass and Instrument in 2. Breath tube heated 3. Date, time and location code (Level 2,E). 4. Print test (Level 1,P). Sign and attach test 5. Tank monitor (Level 3,D,G). Display and a. Display: 189 psi Regulator: 200 psi	Re-set if necessary. st record. Regulator within 50 psi.		
B.	Tests (Sign and attach test records):  Configure simulator for the following tests (Level 1. Wet Calibration Check  a. December 1. Wet Calibration Check  a. December 1. Consum Value 1. Consu	e ≤ 0.03 AC: <u>0.020</u> AC <u>0</u> Exp. Date: <u>2.13.2</u> 2 e ≥ 0.25 AC: <u>0.300</u> AC 2A Exp. Date: <u>12.8.2</u> 2		
	Configure dry gas standard for the following tes  2. Interferent Check (Level 1,B): Known Value Sim. Ser #: DR 3847 Lot #: FCS 6 Explosed Explosed From Explosed Expl	Iue: <u>0./0</u> AC + 0.05% Acetone kp. Date: <u>//A</u> reads "RFI Detect"  n Value <u>0.080</u> AC Exp. Date: <u>//.5.2/</u> Test 7 <u>0.080</u> Test 8 <u>0.079</u> Test 9 <u>0.080</u>		
C.	Remarks/Maintenance: ANNUAL INSPEC	770N		
	nstrument is acceptable to be used in the field.  ath Analyst Signature	6/11/21 Date		
<i>A</i> Revi	riewed by	<u>φ - 14-71</u> Date		

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = TOXL 8164.14.00 09/16
06/11/2021 16:21

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? Standard Expiration? 11/05/2021 Oper No? 666666

Flow Cal. Date: 06/10/2009 Slope 690

Intercept -645883

IR Calibration Date: 04/09/2015

Oth Coef(\*100): -24157 -27270
1st Coef(\*100): 262661 128882
2nd Coef(\*100): 2137 1256
H2O adj(mg/l\*10k): 491 640

3um

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/11/2021

Alcohol Analyzer SN 80-004206 8164.14.00 09/16 16:23

#### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:24
02 Std. Sol.	0.018	16:25
03 Room Air	0.000	16:25
04 Std. Sol.	0.018	16:26
05 Room Air	0.000	16:26
06 Std. Sol.	0.018	16:27
07 Room Air	0.000	16:28

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3064 Std Sol No = 20070

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low AC

Remarks:

Form 106-I8000

0.020 Ac

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/11/2021 SN 80-004206 8164.14.00 09/16 16:38

#### WET CAL CHECK

Test		AC	Time
01 Ro	om Air	0.000	16:38
02 St	d. Sol.	0.294	16:39
03 Ro	om Air	0.000	16:40
04 St	d. Sol.	0.297	16:40
05 Ro	om Air	0.000	16:41
06 St	d. Sol.	0.297	16:42
07 Ro	om Air	0.000	16:42

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP5317 Std Sol No = 202012A

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004206 Location = TOXL 8164.14.00 09/16 06/11/2021 16:43

Test	AC	Time
01 Room Air	0.000	16:44
02 *Subject Test	INT*	16:44
03 Room Air	0.000	16:45

<sup>\*</sup>Invalid Test

Interferent Detected

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = MaleWeight = NA Cit = INTERFERENT Test = DUI

Dr. Lic. = ND/TES989643

Lot No = 24119080A1

Cyl No = 5

Expiration Date = 11/05/2021

County = 08Oper No. = 666666

> Operator Signature CHARLES EDER

Remarks:

Therferent Check 0.10AC + 0.05% ACETONE

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-004206 Location = TOXL 8164.14.00 09/16 06/11/2021

Test	AC	Time
01 Diagnostic	OK	16:46
02 Room Air	0.000	16:46
03 *Subject Test	RFI*	16:46
04 Room Air	0.000	16:47

\*Invalid Test

Inhibited - RFI

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = MaleWeight = NA Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/TES989643

Lot No = 24119080A1

Cyl No = 5

Expiration Date = 11/05/2021

County = 08

Oper No. = 666666

16:45

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

> Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/11/2021

SN 80-004206 8164.14.00 09/16 16:47

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	16:47
02	Std.	Gas	0.079	16:48
03	Room	Air	0.000	16:48
04	Std.	Gas	0.080	16:49
05	Room	Air	0.000	16:49
06	Std.	Gas	0.079	16:49
07	Room	Air	0.000	16:50

Lot No = 24119080A1

Cyl No = 5

Exp Date = 11/05/2021

County = 08

Oper No. = 666666

Operator Signature

Remarks:

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/11/2021

SN 80-004206 8164.14.00 09/16 16:50

# DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:51
02 Std. Gas	0.080	16:51
03 Room Air	0.000	16:52
04 Std. Gas	0.079	16:52
05 Room Air	0.000	16:52
06 Std. Gas	0.080	16:53
07 Room Air	0.000	16:53

Lot No = 24119080A1

Cyl No = 5

Exp Date = 11/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks: (ANBRATION CHECK

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/11/2021

Alcohol Analyzer SN 80-004206 8164.14.00 09/16 16:53

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:54
02 Std. Gas	0.080	16:54
03 Room Air	0.000	16:55
04 Std. Gas	0.079	16:55
05 Room Air	0.000	16:56
06 Std. Gas	0.080	16:56
07 Room Air	0.000	16:57

Lot No = 24119080A1

Cyl No = 5

Exp Date = 11/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Pemarks: ALBOADA

Form 106-I8000

0.080 AC