




**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-006506</b>	Instrument Location <b>Bottineau Co. S.O.</b>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	


Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): None
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <b>06-11-21</b>
--	-------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date <b>6-15-21</b>
--	----------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006506  
Location = BOTT               8164.14.00 09/16  
06/11/2021                     10:30

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Dry  
Standard Value?            0.080  
Standard Lot #?            02621080A1  
Standard Cyl #?            11  
Standard Expiration?       03/05/2023  
Oper No?                    039117

Flow Cal. Date:            08/20/2015  
Slope                       657  
Intercept                  -627881

IR Calibration Date:       08/20/2015  
                            3um                    9um

-----  
0th Coef(\*100):       -22470           -18919  
1st Coef(\*100):       266505           136109  
2nd Coef(\*100):       3702             1644  
H2O adj(mg/l\*10k):   709               510

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
SEAN MURPHY

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006506  
Location = BOTT      8164.14.00 09/16  
06/11/2021      10:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:31
02 Std. Gas	0.081	10:32
03 Room Air	0.000	10:32
04 Std. Gas	0.081	10:33
05 Room Air	0.000	10:33
06 Std. Gas	0.080	10:33
07 Room Air	0.000	10:34

Lot No = 02621080A1  
Cyl No = 11  
Exp Date = 03/05/2023  
County = 05      Oper No. = 039117

  
Operator Signature  
SEAN MURPHY

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006506  
Location = BOTT                    8164.14.00 09/16  
06/11/2021                            10:38

Test	AC	Time
01 Diagnostic	OK	10:40
02 Room Air	0.000	10:41
03 *Subject Test	RFI*	10:41
04 Room Air	0.000	10:41

\*Invalid Test  
Inhibited - RFI


I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
Operator Signature  
SEAN MURPHY

Remarks:

Form 106-I8000

NOTE: NO SUBJECT OR GAS CYLINDER INFORMATION APPEAR ON THIS SUBJECT TEST RECORD. THIS IS DUE TO THE OPERATOR NOT ENTERING A FIRST OR LAST NAME FOR A SUBJECT. THE INSTRUMENT PROGRAMMING IS DESIGNED TO PRODUCE A REPORT SUCH AS THIS WHEN NO SUBJECT INFORMATION IS ENTERED. IT IS NOT AN ERROR FROM THE OPERATOR. THIS IS A VALID TEST RECORD.

  
6.15.21