



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number	80 - 00 5940	Instrument Location	New Town (TATLES)
Reason for Install/Repair			
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____			

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	<i>Al Mat 5661</i>	Date	6-11-21 <del>6-11-21</del>
---------------------------	--------------------	------	-------------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	<i>Charles E. Eh</i>	Certified Date	6-16-21
-----------------------	----------------------	----------------	---------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005940  
Location = NEWT      8164.14.00 09/16  
06/11/2021      20:27

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

Current Instrument Setup


Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      24119080A1  
Standard Cyl #?      26  
Standard Expiration?      11/05/2021  
Oper No?      131739

Flow Cal. Date:      06/21/2013  
Slope      688  
Intercept      -659053

IR Calibration Date:      04/02/2015

	3um	9um
0th Coef(*100):	-18917	-25088
1st Coef(*100):	267616	136158
2nd Coef(*100):	3374	1493
H2O adj(mg/l*10k):	427	417

\*\*\*\*\* Printer Test End \*\*\*\*\*

 5661

Operator Signature  
ANDRES MARTINEZ

Remarks:

Install

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005940  
Location = NEWT      8164.14.00 09/16  
06/11/2021      20:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:37
02 Std. Gas	0.079	20:37
03 Room Air	0.000	20:38
04 Std. Gas	0.079	20:38
05 Room Air	0.000	20:39
06 Std. Gas	0.079	20:39
07 Room Air	0.000	20:40

Lot No = 24119080A1  
Cyl No = 26  
Exp Date = 11/05/2021  
County = 60      Oper No. = 131739

 5661  
Operator Signature  
ANDRES MARTINEZ

Remarks: *Install*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = NEWT                8164.14.00 09/16  
06/11/2021                        20:28

Test	AC	Time
01 Room Air	0.000	20:30
02 *Subject Test	RFI*	20:30
03 Room Air	0.000	20:31

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, SUBJECT NA  
Sub DOB = 09/09/1980  
Sub Sex = Male                    Weight = 200  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 24119080A1  
Cyl No = 26  
Expiration Date = 11/05/2021  
County = 60                        Oper No. = 131739

  
Operator Signature  
ANDRES MARTINEZ

Remarks:

Install

Form 106-I8000