




**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-005357</b>	Instrument Location <b>MOUNTRAIL BOOKING</b>
Reason for Install/Repair	
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	


Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function **B** or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <b>06-11-2021</b>
--	---------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date <b>7/6/21</b>
--	---------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005357  
Location = MOSO                8164.14.00 09/16  
06/10/2021                        18:46

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 031123



Operator Signature  
JOHN TRONSET

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005357  
Location = MOSO                8164.14.00 09/16  
06/11/2021                      20:20

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            05620080A1  
Standard Cyl #?            5  
Standard Expiration?       04/05/2022  
Oper No?                    031123

Flow Cal. Date:            05/16/2018  
Slope                      573  
Intercept                  140460

IR Calibration Date:        06/05/2018  
                                    3um                      9um

-----  
0th Coef(\*100):        -15117                -26541  
1st Coef(\*100):        266971                144348  
2nd Coef(\*100):        2581                  959  
H2O adj(mg/l\*10k):    293                    383

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
JOHN TRONSET

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005357  
Location = MOSO      8164.14.00 09/16  
06/11/2021      20:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:21
02 Std. Gas	0.081	20:21
03 Room Air	0.000	20:22
04 Std. Gas	0.081	20:22
05 Room Air	0.000	20:23
06 Std. Gas	0.082	20:23
07 Room Air	0.000	20:24

Lot No = 05620080A1  
Cyl No = 5  
Exp Date = 04/05/2022  
County = 31      Oper No. = 031123

  
\_\_\_\_\_  
Operator Signature  
JOHN TRONSET

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005357  
Location = MOSO                      8164.14.00 09/16  
06/11/2021    20:24

Test	AC	Time
01 Room Air	0.000	20:26
02 Reference	RFI*	20:26
03 Room Air	0.000	20:26

\*Invalid Test  
Inhibited - RFI

Sub Name = DOE, JOHN F  
Sub DOB = 01/01/1900  
Sub Sex = Male                                      Weight = 500  
Test = DUI                                      Cit = 00000000000  
Dr. Lic. = ND/000000000  
Lot No = 05620080A1  
Cyl No = 5  
Expiration Date = 04/05/2022  
County = 31                                      Oper No. = 031123

  
\_\_\_\_\_  
Operator Signature  
JOHN TRONSET

Remarks:

Form 106-I8000