INTOXILYZER® 8000 INSPECTION

Intoxilyzer[®] 8000 Serial Number: 80.0030-57

Location: 70XL

BrW-004

- A. General Setup and Checks:
 - 1. Diagnostic Tests Pass and Instrument in Ready Mode
 - 2. **X** Breath tube heated
 - 3. A Date, time and location code (Level 2,E). Re-set if necessary.
 - 4. Print test (Level 1,P). Sign and attach test record.
 - 5. A Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: <u>78</u> psi Regulator: <u>75</u> psi
- B. Tests (Sign and attach test records):
 - \boxtimes Configure simulator for the following tests (Level 1,S).
 - 1. Wet Calibration Check
 - a. If Low AC (Level 1,C): Known Value ≤ 0.03 AC: $O_{\bullet} O \overline{SO} AC$ Sim. Ser #: <u>MP5289</u> Lot #: $\overline{SO19/IE}$ Exp. Date: $\overline{II \cdot I9 \cdot 21}$
 - b. A High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.400 AC Sim. Ser #: MP3062 Lot #: $2^{0}2/03E$ Exp. Date: $3 \cdot 24.23$

Configure dry gas standard for the following tests (Level 1,S).

- Interferent Check (Level 1,B): Known Value: <u>0.10</u>AC + 0.05% Acetone Sim. Ser #: <u>DR3897</u> Lot #: <u>ICS6</u> Exp. Date: <u>NA</u>
 Display reads "Interferent Detect"
- 3. XRFI Check (CMS Mode) X Display reads "RFI Detect"
- 4. X Dry Calibration Check (Level 1, C): Known Value <u>0.080</u> AC Gas Cyl Lot #: <u>02621080A1</u> Cyl #: <u>26</u> Exp. Date: <u>3.5.23</u> Test 1 <u>0.080</u> Test 4 <u>0.080</u> Test 7 <u>0.080</u> Test 2 <u>0.079</u> Test 5 <u>0.080</u> Test 8 <u>0.079</u> Test 3 <u>0.080</u> Test 6 <u>0.080</u> Test 9 <u>0.080</u> Average 0.080

Remarks/Maintenance: <u>ANNUAL INSPECTION</u>, REPLACED FORTERNAL ATTERY, REPLACED SIMULATOR RETURN O-RING DURING C. ALIBRATION ADJUSTMENT NOTE: CLASSROOM ONLY UNTIL FLOW SENSOR REPAIRED

Instrument is acceptable to be used in the field. NO

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

Page 1 of 1

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-003057 Location = TOXL 8164.14.00 09/16 06/14/2021 12:46 abcdefqhijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 26 Standard Expiration? 03/05/2023 Oper No? 666666 Flow Cal. Date: 06/14/2021 Slope 516 Intercept 495821 IR Calibration Date: 06/14/2021 3um 9um _____ Oth Coef(*100):-207061st Coef(*100):2641692nd Coef(*100):1074 -12354 140515 859 H2O adj(mg/l*10k): 351239

Operator Signature CHARLES EDER

Remarks:

RINT TEST

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-003057Location = TOXL8164.14.00 09/1606/14/202112:29

	WET CAL CHECK	
Test	AC	Time
01 Room Air	0.000	12:29
02 Std. Sol.	0.030	12:30
03 Room Air	0.000	12:30
04 Std. Sol.	0.031	12:31
05 Room Air	0.000	12:32
06 Std. Sol.	0.030	12:32
07 Room Air	0.000	12:33
	24.000	

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP5289 Std Sol No = 201911E County = 08

Oper No. = 666666

U

Operator Signature CHARLES EDER

Low AC 0.030AC

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-003057Location = TOXL8164.14.00 09/1606/14/202112:39

	WET CAL CHECK	
Time	AC	est
12:40	0.000	Room Air
12:41	0.396	Std. Sol.
12:41	0.000	Room Air
12:42	0.396	Std. Sol.
12:43	0.000	Room Air
12:43	0.397	Std. Sol.
12:44	0.000	Room Air

08 Sim Temp = 34.0° C

Simul Ser No = MP3062 Std Sol No = 202103E County = 08

Oper No. = 666666

Cart 1 1

Operator Signature CHARLES EDER

Remarks:

Form 106-18000

0.4001

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-003057 North Dakota Model 8000 Location = TOXL 8164.14.00 09/16 06/14/2021 12:47 Time Test AC 12:48 01 Room Air 0.000 02 *Subject Test 12:48 INT* 03 Room Air 0.000 12:49 *Invalid Test Interferent Detected Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NATest = DUI Cit = INTERFERENT Dr. Lic. = ND/TES989643Lot No = 02621080A1Cyl No = 26Expiration Date = 03/05/2023County = 08Oper No. = 666666

Operator Signature CHARLES EDER

- Interferent Clock 0.10AC + 0.05% ACETONE C82

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-003057 Location = TOXL 8164.14.00 09/16 06/14/2021 12:49 Time Test AC 12:50 01 Diagnostic OK 02 Room Air 0.000 12:51 03 *Subject Test RFI* 12:51 04 Room Air 0.000 12:51 *Invalid Test Inhibited - RFI Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Cit = RFI CHECK Test = DUI Dr. Lic. = ND/TES989643Lot No = 0.02621080A1Cyl No = 26Expiration Date = 03/05/2023County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

10

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-003057Location = TOXL8164.14.00 09/1606/14/202112:55

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.079 0.000 0.080 0.000	12:55 12:56 12:56 12:57 12:57 12:57 12:58

Lot No = 02621080A1 Cyl No = 26 Exp Date = 03/05/2023 County = 08 Oper No. = 666666

Operator Signature CHARLES EDER A CALIBRATTON CHECK 0.080 AC

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-003057Location = TOXL8164.14.00 09/1606/14/202112:58

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	12:59
02 Std. Gas	0.080	12:59
03 Room Air	0.000	13:00
04 Std. Gas	0.080	13:00
05 Room Air	0.000	13:00
06 Std. Gas	0.080	13:01
07 Room Air	0.000	13:01

Lot No = 02621080A1 Cyl No = 26 Exp Date = 03/05/2023 County = 08

Oper No. = 666666

C

Operator Signature CHARLES EDER

ALIBRATION CHECK 0.080AC

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-003057 Location = TOXL 8164.14.00 09/16 06/14/2021 13:02

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.079 0.000 0.080 0.000	13:02 13:02 13:03 13:03 13:04 13:04 13:05

Lot No = 02621080A1Cyl No = 26Exp Date = 03/05/2023County = 08

Oper No. = 666666

Operator Signature CHARLES EDER ALIBRATION COLECK

Remarks:

Form 106-18000

0.080AC