

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005360

Location: TOXL

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 959 psi Regulator: 1000psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
 - a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC
Sim. Ser #: MP5069 Lot #: 20070 Exp. Date: 2.13.22
 - b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: MP5317 Lot #: 202012A Exp. Date: 12.8.22
 - Configure dry gas standard for the following tests (Level 1,S).
 2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR3847 Lot #: ICS6 Exp. Date: NA
 Display reads "Interferent Detect"
 3. RFI Check (CMS Mode) Display reads "RFI Detect"
 4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 05620080A1 Cyl #: 6 Exp. Date: 4.5.22

Test 1	<u>0.080</u>	Test 4	<u>0.080</u>	Test 7	<u>0.080</u>
Test 2	<u>0.079</u>	Test 5	<u>0.080</u>	Test 8	<u>0.081</u>
Test 3	<u>0.079</u>	Test 6	<u>0.080</u>	Test 9	<u>0.080</u>
Average	<u>0.080</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED
SIMULATOR RETURN O-RING

Instrument is acceptable to be used in the field.

Charles E. Ehr
Breath Analyst Signature

6/16/21
Date

Kali Y. Hieb
Reviewed by

6-21-21
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:53

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=-|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=-|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 24119080A1
Standard Cyl #? 5
Standard Expiration? 11/05/2021
Oper No? 666666

Flow Cal. Date: 06/16/2021
Slope 517
Intercept 187172

IR Calibration Date: 06/16/2021

	3um	9um
0th Coef(*100):	-7908	-18879
1st Coef(*100):	270905	148636
2nd Coef(*100):	959	696
H2O adj(mg/l*10k):	168	275

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:54
02 Std. Sol.	0.020	16:55
03 Room Air	0.000	16:56
04 Std. Sol.	0.021	16:56
05 Room Air	0.000	16:57
06 Std. Sol.	0.021	16:58
07 Room Air	0.000	16:58

08 Sim Temp = 34.0°C

Simul Ser No = MP3064
Std Sol No = 20070
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*
 0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:59

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:59
02 Std. Sol.	0.295	17:00
03 Room Air	0.000	17:00
04 Std. Sol.	0.296	17:01
05 Room Air	0.000	17:02
06 Std. Sol.	0.296	17:02
07 Room Air	0.000	17:03

08 Sim Temp = 34.0°C

Simul Ser No = MP5317
Std Sol No = 202012A
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 17:03

Test	AC	Time
01 Room Air	0.000	17:04
02 *Subject Test	INT*	17:04
03 Room Air	0.000	17:05

*Invalid Test
Interferent Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/TES989643
Lot No = 05620080A1
Cyl No = 6
Expiration Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10AC + 0.05% ACETONE
Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 17:05

Test	AC	Time
01 Diagnostic	OK	17:06
02 Room Air	0.000	17:07
03 *Subject Test	RFI*	17:07
04 Room Air	0.000	17:07

*Invalid Test
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/TES989643
Lot No = 05620080A1
Cyl No = 6
Expiration Date = 04/05/2022
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 17:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:08
02 Std. Gas	0.080	17:08
03 Room Air	0.000	17:09
04 Std. Gas	0.079	17:09
05 Room Air	0.000	17:10
06 Std. Gas	0.079	17:10
07 Room Air	0.000	17:11

Lot No = 05620080A1
Cyl No = 6
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 17:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:11
02 Std. Gas	0.080	17:12
03 Room Air	0.000	17:12
04 Std. Gas	0.080	17:13
05 Room Air	0.000	17:13
06 Std. Gas	0.080	17:13
07 Room Air	0.000	17:14

Lot No = 05620080A1
Cyl No = 6
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 17:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:15
02 Std. Gas	0.080	17:15
03 Room Air	0.000	17:16
04 Std. Gas	0.081	17:16
05 Room Air	0.000	17:16
06 Std. Gas	0.080	17:17
07 Room Air	0.000	17:17

Lot No = 05620080A1
Cyl No = 6
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000