

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-005360 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged ^{A DJUST}
 2. Flow Meter Serial Number: 40655 & 55260 ^{VERIFY}
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation ≥ 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. 175 L/S X 60 Sec/min = 10.50 L/min
 - b. 20 L/min: 0. 332 L/S X 60 Sec/min = 19.92 L/min
 - c. Flow Rates within ± 1 L/min of Expected Value

- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 920 psi Regulator: 900 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)

- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value ≥ 10
 - b. Auto Range Res Value ≥ 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

| Soln. | g/210 L | Lot No. | Exp. Date | Simulator SN |
|-------|-------------------|---------------------------------|------------------------------|--------------|
| 1 | 0.000 (ACTUAL) | NA - MilliQ H ₂ O | NA - MilliQ H ₂ O | MP3057 |
| 2 | 0.040 (0.041) | 20203AA ^{03A} | 3.10.22 | MP3059 |
| 3 | 0.080 (0.080) | 21050 | 2.15.23 | MP5318 |
| 4 | 0.100 (0.102) | 202010E | 10.20.22 | MP3003 |
| 5 | 0.300 (0.298) | 20030 | 1.21.22 | MP3069 |

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 07220100A1 Cyl No. 4 Exp. Date: 5.5.22
4. Atmospheric Pressure
 - a. 916 mbar Displayed by Intoxilyzer® 8000
 - b. 952 mbar Adjusted to using barometer
 - c. 952 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached
- a. Solution 1 Avg % Abs ≤ 0.2500
 - b. Solution 2-5 REL STD DEV ≤ 3.000
 - c. Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μm and 9 μm channels
 - d. Dry Gas H₂O Adjustment Sum for 3 μm and 9 μm channels within ± 10

| | Average | | H ₂ O Adjust | |
|-----------------|-------------|---|-------------------------|---------------|
| 3 μm | <u>4593</u> | + | <u>168</u> | = <u>4761</u> |
| 9 μm | <u>4486</u> | + | <u>275</u> | = <u>4761</u> |

7. Optical Bench Calibration Verification (Level 1, S and C)
- a. Wet Calibration Check
 - i. Low AC Known Value ≤ 0.03 AC: 0.030 AC
Sim. SN: MP5289 Lot No.: 201911E Exp. Date: 11.19.21
 - ii. High AC Known Value ≥ 0.25 AC: 0.400 AC
Sim. SN: MP3062 Lot No.: 202103E Exp. Date: 3.24.23
 - b. Dry Calibration Check: Known Value 0.08 AC
Lot No. 24119080 A1 Cyl No. 5 Exp. Date: 11.5.21
Test 1 0.080 AC Test 4 0.081 AC Test 7 0.080 AC
Test 2 0.080 AC Test 5 0.080 AC Test 8 0.081 AC
Test 3 0.080 AC Test 6 0.080 AC Test 9 0.080 AC
Average 0.080 AC
 - c. Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: CAL. ADJ. DUE TO ATMOSPHERIC
SENSOR (BAROMETER) READING 916 mbar WHEN ACTUAL
ATMOSPHERIC PRESSURE WAS 952 mbar.

Instrument is acceptable to be used in the field.

Charles E. Ehr
Breath Analyst Signature

6/16/21
Date

Kali L. Hines
Reviewed by

6-21-21
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 14:24

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
 SQRT(Diff) = 2.000
2: Rate (Liters/min) = 15
 SQRT(Diff) = 10.582
3: Rate (Liters/min) = 30
 SQRT(Diff) = 20.949

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 517

Rounded Intercept = 187172

Correlation = 0.99816



TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 06/16/2021 15:53:31

Auto Calibration

pg 1 of 2

```

    <<<<<      3um      >>>>>
    -----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs   (% Abs Ref)   % Abs   (% Abs Ref)
Sample #1  0.0360   (0.0190)     0.1490   (-0.0150)
Sample #2  0.0230   (0.0480)     0.1050   (0.0270)
Sample #3  0.0330   (0.0560)     0.1220   (0.0020)
Sample #4  0.0040   (0.1090)     0.1270   (0.0050)
Avg % Abs  0.0200   (0.0710)     0.1180   (0.0113)
STD DEV    0.0147   (0.0332)     0.0115   (0.0137)
REL STD DEV 73.655   (46.692)     9.773    (120.445)
  
```

```

    <<<<<      9um      >>>>>
    -----
Solution = 0.041 g/210L or 0.1952 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs   (% Abs Ref)   % Abs   (% Abs Ref)
Sample #1  0.7230   (0.0180)     1.4100   (0.0270)
Sample #2  0.7840   (0.0260)     1.4500   (0.0120)
Sample #3  0.7540   (0.0560)     1.4690   (0.0300)
Sample #4  0.7620   (0.0610)     1.4540   (0.0290)
Avg % Abs  0.7667   (0.0477)     1.4577   (0.0237)
STD DEV    0.0155   (0.0189)     0.0100   (0.0101)
REL STD DEV 2.026   (39.713)     0.687    (42.744)
  
```

```

    <<<<<      3um      >>>>>
    -----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs   (% Abs Ref)   % Abs   (% Abs Ref)
Sample #1  1.4190   (0.0120)     2.6860   (0.0390)
Sample #2  1.4620   (0.0090)     2.6790   (0.0460)
Sample #3  1.4010   (0.0540)     2.5980   (0.0750)
Sample #4  1.4180   (0.0450)     2.6110   (0.0800)
Avg % Abs  1.4270   (0.0360)     2.6293   (0.0670)
STD DEV    0.0315   (0.0238)     0.0435   (0.0184)
REL STD DEV 2.206   (66.144)     1.654    (27.399)
  
```

```

    <<<<<      9um      >>>>>
    -----
Solution = 0.102 g/210L or 0.4857 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs   (% Abs Ref)   % Abs   (% Abs Ref)
Sample #1  1.7630   (0.0130)     3.3680   (0.0140)
Sample #2  1.7850   (0.0210)     3.3530   (0.0390)
Sample #3  1.8420   (0.0130)     3.3580   (0.0370)
Sample #4  1.7770   (0.0450)     3.3450   (0.0600)
Avg % Abs  1.8013   (0.0263)     3.3520   (0.0453)
STD DEV    0.0354   (0.0167)     0.0066   (0.0127)
REL STD DEV 1.968   (63.240)     0.196    (28.105)
  
```

```

    <<<<<      3um      >>>>>
    -----
Solution = 0.298 g/210L or 1.4190 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs   (% Abs Ref)   % Abs   (% Abs Ref)
Sample #1  5.1700   (0.0180)     9.2570   (0.0000)
Sample #2  5.1980   (0.0530)     9.3100   (0.0330)
Sample #3  5.1720   (0.0260)     9.2650   (0.0110)
Sample #4  5.1500   (0.0480)     9.2400   (0.0200)
Avg % Abs  5.1733   (0.0423)     9.2717   (0.0213)
STD DEV    0.0240   (0.0144)     0.0355   (0.0111)
REL STD DEV 0.464   (33.931)     0.383    (51.846)
  
```

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 06/16/2021 15:53:31

Auto Calibration

pg 2 of 2

<<<< 3um >>>>

 Zero Order Coef -79.08
 First Order Coef 2709.05
 Second Order Coef 9.59

<<<< 9um >>>>

 -188.79
 1486.36
 6.96

| Act (g/210L) | Fit (g/210L) | Residual (g/210L) |
|-----------------|-----------------|----------------------|
| 0.000 | -0.001 | 0.0005 |
| 0.041 | 0.042 | -0.0011 |
| 0.080 | 0.080 | 0.0001 |
| 0.102 | 0.101 | 0.0005 |
| 0.298 | 0.298 | -0.0000 |

| Act (g/210L) | Fit (g/210L) | Residual (g/210L) |
|-----------------|-----------------|----------------------|
| 0.000 | -0.000 | 0.0003 |
| 0.041 | 0.042 | -0.0008 |
| 0.080 | 0.079 | 0.0009 |
| 0.102 | 0.102 | -0.0003 |
| 0.298 | 0.298 | -0.0000 |


<<<< 3um >>>> <<<< 9um >>>>

 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

| Sample | 3um | 9um |
|-----------------------|-----------|-----------|
| Sample #1 | 4640.00 | 4531.00 |
| Sample #2 | 4572.00 | 4505.00 |
| Sample #3 | 4632.00 | 4479.00 |
| Sample #4 | 4576.00 | 4475.00 |
| Avg | 4593.3335 | 4486.3335 |
| STD DEV | 33.5460 | 16.2891 |
| REL STD DEV | 0.730 | 0.363 |
| H2O adjust (mg/l*10k) | 168 | 275 |

Atmospheric Pressure = 952

*****CALIBRATION SUCCESSFUL*****



TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 06/16/2021 15:53:31

Auto Calibration
 Max Power Res Value = 38
 Auto Range Res Value = 20

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:32

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 16:33 |
| 02 Std. Sol. | 0.030 | 16:34 |
| 03 Room Air | 0.000 | 16:34 |
| 04 Std. Sol. | 0.029 | 16:35 |
| 05 Room Air | 0.000 | 16:36 |
| 06 Std. Sol. | 0.030 | 16:36 |
| 07 Room Air | 0.000 | 16:37 |

08 Sim Temp = 34.0°C

Simul Ser No = MP5289
Std Sol No = 201911E
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low Ac*
0.030 Ac

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:37

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 16:38 |
| 02 Std. Sol. | 0.397 | 16:39 |
| 03 Room Air | 0.000 | 16:39 |
| 04 Std. Sol. | 0.397 | 16:40 |
| 05 Room Air | 0.000 | 16:41 |
| 06 Std. Sol. | 0.397 | 16:41 |
| 07 Room Air | 0.000 | 16:42 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3062
Std Sol No = 202103E
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.400 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:42

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:43 |
| 02 Std. Gas | 0.080 | 16:43 |
| 03 Room Air | 0.000 | 16:44 |
| 04 Std. Gas | 0.080 | 16:44 |
| 05 Room Air | 0.000 | 16:45 |
| 06 Std. Gas | 0.080 | 16:45 |
| 07 Room Air | 0.000 | 16:46 |

Lot No = 24119080A1
Cyl No = 5
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:46

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:46 |
| 02 Std. Gas | 0.081 | 16:47 |
| 03 Room Air | 0.000 | 16:47 |
| 04 Std. Gas | 0.080 | 16:48 |
| 05 Room Air | 0.000 | 16:48 |
| 06 Std. Gas | 0.080 | 16:48 |
| 07 Room Air | 0.000 | 16:49 |

Lot No = 24119080A1
Cyl No = 5
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK*
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
06/16/2021 16:49

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:50 |
| 02 Std. Gas | 0.080 | 16:50 |
| 03 Room Air | 0.000 | 16:51 |
| 04 Std. Gas | 0.081 | 16:51 |
| 05 Room Air | 0.000 | 16:52 |
| 06 Std. Gas | 0.080 | 16:52 |
| 07 Room Air | 0.000 | 16:52 |

Lot No = 24119080A1
Cyl No = 5
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK*
0.080 AC

Form 106-I8000