

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION Toxicology Section/Breath Alcohol Program SFN 59281 (10-2016)

Serial Number	Instrument Location			
30.0064985	NTPP			
Reason for Install/Repair	_			
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Properly Ground	led.			
2. Telephone Line Connected to Intoxilyzer®	8000.			
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, D	ST (Y), and Location; Level 2,	Function E).		
5. Scan/Enter Gas Cylinder Information (Lev	el 1, Function S).			
6. Run Tests:				
A. Print Test (Level 1, Function P).				
B. ACA Test (Level 1, Function C).				
C. RFI Test (CMS Mode or Level 1, Func	tion B or C; Key Radio During	Test).		
7. Repair and/or Maintenance Performed (if	Any): NA			
8. Complete the Top Portion of the Intoxilyze Intoxilyzer® for Use.	r® Record (SFN50496, Form 1	120-G) and Place it by the		
9. File Previous Intoxilyzer® Record (SFN50	496, Form 120-G) at the Intoxil	yzer® Location at the Agency.		
10. Send the Following to the Crime Laborator (SFN59281, Form 104-G), Print Test, ACA	y: Completed Intoxilyzer®800 Test, and RFI Test.	0 Installation and Repair Checkou		
Field Inspector Signature		Date		
Jak H		June 18, 202		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is a alcohol concentration from the date the Field Inspector per Attorney General, Crime Laboratory Division, in the County copy of the documents received.	formed the installation. This re	cord on file at the Office of		
Reviewed/Certified by	660	Certified Date		
1 XInver		6.18-2021		

Form 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006498
Location = NTPD 8164.14.00 09/16
06/18/2021 12:26

********* Printer Test *********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&*()_+?

Current Instrument Setup Enabled Data Entry Mode: DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Yes Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 24119080A1 10 Standard Expiration? 11/05/2021 130826 Oper No?

Flow Cal. Date: 08/18/2015 Slope 669 Intercept -637469

IR Calibration Date: 08/18/2015
3um 9um

Oth Coef(*100): -28161 -22508
1st Coef(*100): 286657 143843
2nd Coef(*100): 2341 1281
H2O adj(mg/l*10k): 776 451

******* Printer Test End ********

Óperator Signature JACQUELYN HALONEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006498
Location = NTPD 8164.14.00 09/16 06/18/2021

12:26

DRY CAL CHECK

Test		AC	Time
01	Room Air	0.000	12:27
	Std. Gas	0.081	12:27
	Room Air	0.000	12:28
04	Std. Gas	0.080	12:28
05	Room Air	0.000	12:28
06	Std. Gas	0.081	12:29
	Room Air	0.000	12:29

Lot No = 24119080A1

Cyl No = 10

Exp Date = 11/05/2021

County = 31Oper No. = 130826

> Operator Signature JACQUELYN HALONEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = NTPD 06/18/2021

SN 80-006498 8164.14.00 09/16 12:30

Test	AC	Time
01 Room Air	RFI*	12:32
02 Room Air	0.000	12:32

*Invalid Test Inhibited - RFI

Sub Name = DOE, JOHN DAVID

Sub DOB = 12/12/1900

Weight = 210 Sub Sex = Male Cit = 000000Test = DUI

Dr. Lic. = ND/DOE190012

Lot No = 24119080A1

Cyl No = 10

Expiration Date = 11/05/2021

Oper No. = 130826County 31 Operator Signature JACQUELYN HALONEN

Remarks:

Form 106-I8000