

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005360	rument Location Grant For FS Police Department			
Reason for Install/Repair	Oran 2 TO: FS TO TEC Department			
Linstall After Receiving From Crime Laboratory	cation Change			
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
2. A. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Leve	l 1, Function B or C; Key Radio During Test).			
💢 7. Repair and/or Maintenance Performed (if any): non 🕰				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyze 104-G), Print Test, ACA Test, and RFI Test.	er® 8000 Installation and Repair Checkout (SFN59281, Form			
Field Inspector Signature	Date			
Hon Lanps	06-18-21			
Crime Laboratory Use Only				

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Chana El	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005360 8164.14.00 09/16 Location = GFPD 15:41 06/18/2021 ************* Printer Test ************** abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Enabled Data Entry Mode. Start Test Sequence: DABA Display Prelim Rslt? Yes Display Third Digit? Yes Data Entry Mode: DABACABA No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 35 Standard Expiration? 03/05/2023 039947 Oper No? 06/16/2021 Flow Cal. Date: 517 Slope 187172 Intercept IR Calibration Date: 06/16/2021 3um 9um _____ Oth Coef(*100):-79081st Coef(*100):2709052nd Coef(*100):959 -18879 148636 696 275 H2O adj(mg/l*10k): 168

Operator Signature JONATHAN LAMPI

Remarks:

Install

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005360Location = GFPD8164.14.00 09/1606/18/202115:42

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.080 0.000 0.079 0.000	15:43 15:43 15:43 15:44 15:44 15:45 15:45

Operator Signature JONATHAN LAMPI

Remarks:

Form 106-I8000

Install

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-005360 North Dakota Model 8000 8164.14.00 09/16 Location = GFPD15:49 06/18/2021

DRY CAL	CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Room Air	0.000 0.079 RFI* 0.000	15:49 15:49 15:50 15:50
*Invalid Test Inhibited - RFI		
Lot No = 02621080A1 Cyl No = 35 Exp Date = 03/05/2023 County = 18	Oper No.	= 039947

Operator Signature

JONATHAN LAMPI

Remarks:

Form 106-18000

Install