

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 65 05	Instrument Location She Malls	Office Mohl
Reason for Install/Repair		
Install After Receiving From Crime Laboratory Install	After Location Change	
Other (Specify)		
Check When Done;		
💢 1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
🐧 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loc	ation, Level 2, Function E).	
5 Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:		
A. Print Test (Level1, Function P)		
B. ACA Test (Level 1, Function C).		
Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1, Function B or C; Key Radio D	ouring Test).
7. Repair and/or Maintenance Performed (if any):	esued from hah	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120	-G) at the Intoxilyzer® Location at the Ag	ency.
10. Send the Following to the Crime Laboratory: Completed I 104-G), Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair	Checkout (SFN59281, Form
Field Inspector Signature		Date 8 - 2 - 71)
Crime Laboratory Use Only		·
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct of the county of Burleigh.	file at the Office of Attorney General, Cri	ermine alcohol concentration from me Laboratory Division, in the
Réviewed/Certified By	El	Certified Date 8 • 2 • 2

FORM 104-G

SN 80-006505 Page 1 of 4 Uploaded 8/2/2021 CEE

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006505
Location = MOHL 8164.14.00 09/16
08/02/2021 16:02

********** Printer Test ********

abcdefghijklmnopgrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ1@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ1@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Display Third Digit? Yes Yes Inhib Printer(Y/N)? No Display Volume? Νo Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 24119080A1 3.4 11/05/2021 Standard Expiration? 131291 Oper No?

Flow Cal. Date: 08/20/2015 Slope 671 Intercept -642269

IR Calibration Date: 08/20/2015
3um 9um

Oth	Coef(*100):	-29025	-17051
1st	Coef(*100):	252705	133345
	Coef(*100):	3840	1546
H20	adj(mg/1*10k):	817	461

****** Printer Test End ********

Operator Signature ROGER HUTCHINSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006505 Location = MOHL 08/02/2021

8164.14.00 09/16 16:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:05
02 Std. Gas	0.080	16:05
03 Room Air	0.000	16:06
04 Std. Gas	0.079	16:06
05 Room Air	0.000	16:06
06 Std. Gas	0.080	16:07
07 Room Air	0.000	16:07

Lot No = 24119080A1

Cyl No = 34

Exp Date = 11/05/2021

County = 38

Oper No. = 131291

perator signature ROCER HUTCHINSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = MOHL 08/02/2021 Alcohol Analyzer SN 80-006505 8164.14.00 09/16 16:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:25
02 Std. Gas	RFI*	16:25
03 Room Air	0.000	16:26

*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 34

Exp Date = 11/05/2021

County = 38

Oper No. = 131291

Operator Signature ROGER HUTCHINSON

Remarks:

Form 106-I8000