

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-66 5943	Instrument Location N CSO		
Reason for Install/Repair			
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change		
Other (Specify) re-install offer storage ((Covid)		
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.	
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature		Date 10 / Cl /a a /	
		10/14/2021	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be			
the date the Field Inspector performed the installation. This record on fil County of Burleigh, North Dakota, is certified to be a true and correct copies.		ne Laboratory Division, in the	

FORM 104-G

Reviewed/Certified By

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005943
Location = NCSO 8164.14.00 09/16
10/09/2021 19:13

****** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Prelim Rsit:
Display Third Digit?
Yes
Third Digit?
No Display Volume? No Disable On Memfull? # of Print Copies? Yes Dry Select Std (D/W/I)? Standard Value? 0.080 Standard Lot #? Standard Cyl #? 02621080A1 15 03/05/2023 Standard Expiration? Oper No? 100754

Flow Cal. Date: 06/20/2013 Slope 686 Intercept -735246

IR Calibration Date: 08/26/2016
3um 9um
0th Coef(*100): -13589 -19547

Oth Coef(*100): -13589 -19547
1st Coef(*100): 284566 139614
2nd Coef(*100): 2717 1336
H2O adj(mg/l*10k): 381 404

****** Printer Test End *******

Operator Signature
DANIELLE BJORLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005943 Location = NCSO 10/14/2021

8164.14.00 09/16 09:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:09
02 Std. Gas	0.079	09:09
03 Room Air	0.000	09:10
04 Std. Gas	0.078	09:10
05 Room Air	0.000	09:11
06 Std. Gas	0.079	09:11
07 Room Air	0.000	09:11

Lot No = 02621080A1

Cyl No = 15

Exp Date = 03/05/2023

County = 32

Oper No. = 100754

Operator Signature DANIELLE BJORLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
10/14/2021

Alcohol Analyzer SN 80-005943 8164.14.00 09/16 09:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:22
02 Std. Gas	0.079	09:22
03 Room Air	RFI*	09:22
04 Room Air	0.000	09:23

*Invalid Test Inhibited - RFI

Lot No = 02621080A1

Cyl No = 15

Exp Date = 03/05/2023

County = 32

Oper No. = 100754

Operator Signature DANJELLE BJORLIE

Remarks:

Form 106-I8000