




**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number 80-005953	Instrument Location DICKEY COUNTY COURTHOUSE, ELLENDALE ND
Reason for Install/Repair	
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	


Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date 10/28/2021
--	--------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date 11/2/21
--	---------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005953  
Location = DCSO                    8164.14.00 09/16  
10/28/2021                            15:37

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            24119080A1  
Standard Cyl #?            15  
Standard Expiration?        11/05/2021  
Oper No?                    020504

Flow Cal. Date:            06/16/2015  
Slope                      678  
Intercept                  -668987

IR Calibration Date:        04/30/2015  
                                    3um                      9um

-----  
0th Coef(\*100):            -13435                      -19950  
1st Coef(\*100):            257165                      129445  
2nd Coef(\*100):            3995                        1631  
H2O adj(mg/l\*10k):        505                         548

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHRIS ESTES

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005953  
Location = DCSO      8164.14.00 09/16  
10/28/2021      15:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:38
02 Std. Gas	0.081	15:39
03 Room Air	0.000	15:39
04 Std. Gas	0.082	15:40
05 Room Air	0.000	15:40
06 Std. Gas	0.081	15:41
07 Room Air	0.000	15:41

Lot No = 24119080A1  
Cyl No = 15  
Exp Date = 11/05/2021  
County = 11      Oper No. = 020504

  
\_\_\_\_\_  
Operator Signature  
CHRIS ESTES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005953  
Location = DCSO      8164.14.00 09/16  
10/28/2021      15:42

Test	AC	Time
01 Room Air	RFI*	15:47
02 Room Air	0.000	15:47

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST RADIO  
Sub DOB = 10/28/2020  
Sub Sex = Male      Weight = 100  
Test = OTH      Cit = 0000  
Dr. Lic. = ND/00000000  
Lot No = 24119080A1  
Cyl No = 15  
Expiration Date = 11/05/2021  
County = 11      Oper No. = 020504



Operator Signature  
CHRIS ESTES

Remarks:

Form 106-I8000