

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004945	Instrument Location BMDC		
Reason for Install/Repair			
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change			
Other (Specify) Calibration(s) after power went out at BMDC in early mornings of 10/28/21			
Check When Done:			
▼ 1. Surge Protector Installed/Property Grounded.			
ヌ 3. Breath Tube Heated.			
X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
▼ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
★ 6. Run Tests:			
■ A. Print Test (Level1, Function P).			
☑ B. ACA Test (Level 1, Function C).			
▼ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
X 7. Repair and/or Maintenance Performed (if any):			
▼ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.	
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature 4834		Date 10/28/2021	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.			
Reviewed/Certified By		Certified Date	

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004945
Location = BMDC 8164.14.00 09/16
10/28/2021 10:03

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 02621080A1 Standard Expiration? 03/05/2023 Oper No? 130315

Flow Cal. Date: 09/26/2012 Slope 622 Intercept -485836

IR Calibration Date: 04/22/2016
3um 9um

Oth Coef(\*100): -23427 -26406 1st Coef(\*100): 266989 144341 2nd Coef(\*100): 2085 1078 H2O adj(mg/l\*10k): 539 400

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
JARED LEMIEUX

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = BMDC 8164.14.00 09/16 10/28/2021

SN 80-004945 09:59

## DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas	0.000 0.081 0.000 0.081 0.000 0.080	10:00 10:00 10:01 10:01 10:01
07 Room Air	0.000	10:02

Lot No = 02621080A1

Cyl No = 7

Exp Date = 03/05/2023

County = 08

Remarks:

Oper No. = 130315

Form 106-I8000

Operator Signature JARED LEMIEUX

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
10/28/2021

SN 80-004945 8164.14.00 09/16 10:04

## DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Room Air	RFI*	10:04 10:04

\*Invalid Test Inhibited - RFI

Lot No = 02621080A1

Cyl No = 7

Exp Date = 03/05/2023

County = 08

Oper No. = 130315

Remarks:

Form 106-I8000

Operator Signaty
JARED LEMIEUX