

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004945	Instrument Location BMDC		
Reason for Install/Repair			
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	r Location Change		
☑ Other (Specify) December Routine Check			
Check When Done:			
▼ 1. Surge Protector Installed/Property Grounded.			
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.			
ヌ 3. Breath Tube Heated.			
▼ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
▼ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
▼ 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
▼ C. Radio Frequency Interference (RFI) Test (CMS Mode or)	Level 1, Function B or C; Key Radio Du	ring Test).	
X 7. Repair and/or Maintenance Performed (if any):			
▼ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50)	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.	
🗵 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.	
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature 4834		Date 12/5/2021	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct copies.	e at the Office of Attorney General, Crin	mine alcohol concentration from ne Laboratory Division, in the	
Reviewed/Certified By		Certified Date 2/25 2021	

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004945
Location = BMDC 8164.14.00 09/16
12/05/2021 01:03

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes No Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? 02621080A1 Standard Lot #? Standard Cyl #? Standard Expiration? 03/05/2023 130315 Oper No?

Flow Cal. Date: 09/26/2012 Slope 622 Intercept -485836

IR Calibration Date: 04/22/2016
3um 9um

 0th Coef(*100):
 -23427
 -26406

 1st Coef(*100):
 266989
 144341

 2nd Coef(*100):
 2085
 1078

 H2O adj(mg/l*10k):
 539
 400

******* Printer Test End ********

Operator Signature
JARED LEMIEUX

Remarks:

December checks (2001)

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
12/05/2021

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 01:03

DRY CAL CHECK

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std. 07 Room	Gas Air Gas Air Gas	0.000 0.080 0.000 0.081 0.000 0.081 0.000	01:05 01:05 01:05 01:06 01:06 01:07

Lot No = 02621080A1

Cyl No = 7

Exp Date = 03/05/2023

County = 08

Oper No. = 130315

Remarks:

Dec. charks (2021)

JARED LEMIEUX

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
12/05/2021

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 01:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	01:11
02 Std. Gas	0.082	01:12
03 Room Air	RFI*	01:12
04 Room Air	RFI*	01:13

*Invalid Test Inhibited - RFI

Lot No = 02621080A1

Cyl No = 7

Exp Date = 03/05/2023

County = 08

Oper No. = 130315

Operator Signature JARED LEMIEUX

Remarks:

Dre. Routine Chark (2031)

Form 106-I8000