



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-004950</i>	Instrument Location <i>Cakes PD</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): *None*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Matt C. Brown</i>	Date <i>3/9/22</i>
---	-----------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>Charles E. Eide</i>	Certified Date <i>3/10/2022</i>
---	------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004950  
Location = OKPD                8164.14.00 09/16  
03/09/2022                      21:40

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            26021080A1  
Standard Cyl #?            021  
Standard Expiration?       10/05/2023  
Oper No?                    130690

Flow Cal. Date:            01/18/2010  
Slope                      774  
Intercept                  -908904

IR Calibration Date:        05/05/2016  
                                  3um                            9um

-----  
0th Coef(\*100):            -13717                        -19495  
1st Coef(\*100):            271916                        140513  
2nd Coef(\*100):            3504                            1503  
H2O adj(mg/l\*10k):        370                            426

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
MATTHEW O'BRIEN

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004950  
Location = OKPD      8164.14.00 09/16  
03/09/2022      21:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:43
02 Std. Gas	0.080	21:43
03 Room Air	0.000	21:44
04 Std. Gas	0.081	21:44
05 Room Air	0.000	21:45
06 Std. Gas	0.081	21:45
07 Room Air	0.000	21:45

Lot No = 26021080A1  
Cyl No = 021  
Exp Date = 10/05/2023  
County = 11      Oper No. = 130690

  
\_\_\_\_\_  
Operator Signature  
MATTHEW O'BRIEN

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004950  
Location = OKPD                    8164.14.00 09/16  
03/09/2022                            21:53

Test	AC	Time
01 Room Air	RFI*	21:54
02 Room Air	0.000	21:55

\*Invalid Test  
Inhibited - RFI

Sub Name = RFI TEST, NA NA  
Sub DOB = 01/01/2001  
Sub Sex = Unknown                    Weight = NA  
Test = OTH                              Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 26021080A1  
Cyl No = 021  
Expiration Date = 10/05/2023  
County = 11                              Oper No. = 130690

  
\_\_\_\_\_  
Operator Signature  
MATTHEW O'BRIEN

Remarks:

Form 106-I8000