

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004943	Instrument Location BCCC		
Reason for Install/Repair			
Install After Receiving From Crime Laboratory	er Location Change		
Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
(X) 2. Telephone Line Connected to Intoxilyzer® 8000.			
∑ 3. Breath Tube Heated.			
(i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).		
[X] 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
(A) 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or I	Level 1, Function B or C; Key Radio During Test).		
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	496, Form 120-G) and Place it by the Intoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	xilyzer® 8000 Installation and Repair Checkout (SFN59281, Form		
Field Inspector Signature	Date _		
De XIII	2/10/2022		
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be he date the Field Inspector performed the installation. This record on file County of Burleigh, North Dakota, is certified to be a true and correct copy	at the Office of Attorney General, Crime Laboratory Division, in the		
Reviewed/Certified By	Certified Date		
	7 3/1/2420		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = BCCC 8164.14.00 09/16
03/10/2022 13:45

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 26021080A1 Standard Cyl #? 12 Standard Expiration? 10/05/2023 100012 Oper No?

Flow Cal. Date: 04/06/2020 Slope 679 Intercept -694345

IR Calibration Date: 04/06/2020 3um 9um

Oth Coef(*100): -23601 -21377
1st Coef(*100): 255715 132279
2nd Coef(*100): 1319 1274
H2O adj(mg/l*10k): 424 419

****** Printer Test End ********

Operator Signature DANA RUSTEBAKKE

Remarks:

Print Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BCCC
03/10/2022

Alcohol Analyzer SN 80-004943 8164.14.00 09/16 13:45

DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas	0.000 0.082 0.000 0.081 0.000 0.081	13:46 13:46 13:47 13:47 13:47
07 Room Air	0.000	13:48

Lot No = 26021080A1

Cyl No = 12

Exp Date = 10/05/2023

County = 02

Oper No. = 100012

Operator Signature
DANA RUSTEBAKKE

Remarks:

ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BCCC
03/10/2022

Alcohol Analyzer SN 80-004943 8164.14.00 09/16 13:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:49
02 Std. Gas	RFI*	13:50
03 Room Air	0.000	13:50

*Invalid Test Inhibited - RFI

Lot No = 26021080A1

Cyl No = 12

Exp Date = 10/05/2023

County = 02

Oper No. = 100012

Operator Signature DANA RUSTEBAKKE

Remarks:

RFI Tot

Form 106-I8000