




INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59281 (06/2018)

| | | |
|--|--|------|
| Serial Number 80-006502 | Instrument Location Stutsman Co SO - downstairs | STVT |
| Reason for Install/Repair | | |
| <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____ | | |

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): NA
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

| | | |
|---------------------------|--|-----------------|
| Field Inspector Signature |  358 (010319) | Date 3/11/22 |
|---------------------------|--|-----------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | | |
|-----------------------|---|-----------------------------|
| Reviewed/Certified By |  | Certified Date 3/11/2022 |
|-----------------------|---|-----------------------------|

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006502
Location = STUT 8164.14.00 09/16
03/11/2022 01:48

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 26021080A1
Standard Cyl #? 5
Standard Expiration? 10/05/2023
Oper No? 010319

Flow Cal. Date: 08/19/2015
Slope 674
Intercept -654491

IR Calibration Date: 08/19/2015
 3um 9um

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -24446 | -18103 |
| 1st Coef(*100): | 267787 | 135417 |
| 2nd Coef(*100): | 3387 | 1610 |
| H2O adj(mg/l*10k): | 740 | 513 |

***** Printer Test End *****


Operator Signature
RICK WIDICKER

Remarks:

Print test good!

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006502
Location = STUT 8164.14.00 09/16
03/11/2022 01:51

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 01:52 |
| 02 Std. Gas | 0.082 | 01:52 |
| 03 Room Air | 0.000 | 01:53 |
| 04 Std. Gas | 0.082 | 01:53 |
| 05 Room Air | 0.000 | 01:54 |
| 06 Std. Gas | 0.082 | 01:54 |
| 07 Room Air | 0.000 | 01:55 |

Lot No = 26021080A1
Cyl No = 5
Exp Date = 10/05/2023
County = 47 Oper No. = 010319


Operator Signature
RICK WIDICKER

Remarks: *ACA TEST - Good*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006502
Location = STUT 8164.14.00 09/16
03/11/2022 01:57

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | RFI* | 02:01 |
| 02 Room Air | 0.000 | 02:02 |

*Invalid Test
Inhibited - RFI

Sub Name = RFI TEST INSTALL, INSTALL RFI
TEST
Sub DOB = 01/01/1980
Sub Sex = Unknown Weight = NA
Test = OTH Cit = NA
Dr. Lic. = ND/NA
Lot No = 26021080A1
Cyl No = 5
Expiration Date = 10/05/2023
County = 47 Oper No. = 010319


Operator Signature
RICK WIDICKER

Remarks: *RFI TEST - Good*

Form 106-I8000