

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 004 944	Instrument Location
Reason for Install/Repair	WILL
Install After Receiving From Crime Laboratory	☐ Install After Location Change
Other (Specify)	
heck When Done:	
1. Surge Protector Installed/Property Grounded	ı.
2. Telephone Line Connected to Intoxilyzer® 80	000.
☑ 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST	(Y), and Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1	1, Function S).
☑ 6. Run Tests:	
A. Print Test (Level1, Function P).	
☑ B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test	(CMS Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any	/):
8. Complete the Top Portion of the Intoxilyzer®	Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN50409	96, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: 104-G), Print Test, ACA Test, and RFI Test.	Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
ield Inspector Signature	Date 3-/7- 3032-
frank sin	
rime Laboratory Use Only	
	s approved to be used for the analysis of breath to determine alcohol concentration his record on file at the Office of Attorney General, Crime Laboratory Division, in the

County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	lealer Et	Certified Date 3/22	/22
		/	

Page 1 of 4

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = WILL 8164.14.00 09/16
03/17/2022 19:05

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 26021080A1 Standard Cyl #? 17 10/05/2023 Standard Expiration? Oper No? 130755

Flow Cal. Date: 07/01/2011 Slope 685 Intercept -634740

IR Calibration Date: 07/01/2011
3um 9ur

Oth Coef(*100): -22591 -17234
1st Coef(*100): 265915 132691
2nd Coef(*100): 2841 1562
H2O adj(mg/l*10k): 559 425

****** Printer Test End *******

Operator Signature
JUSTIN PELZL

Remarks: Install

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL

03/17/2022

Alcohol Analyzer SN 80-004944 8164.14.00 09/16 19:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:28
02 Std. Gas	0.082	19:29
03 Room Air	0.000	19:29
04 Std. Gas	0.082	19:29
05 Room Air	0.000	19:30
06 Std. Gas	0.082	19:30
07 Room Air	0.000	19:31

Lot No = 26021080A1

Cyl No = 17

Exp Date = 10/05/2023

County = 53

Oper No. = 130755

WARY #190

Operator Signature
JUSTIN PELZL

Remarks:

Install

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = WILL 8164.14.00 09/16
03/17/2022 19:31

Т	est	AC	Time
01	Room Air	0.000	19:33
02	*Subject T	'est RFI*	19:33
03	Room Air	0.000	19:34

*Invalid Test Inhibited - RFI

Sub Name = PELZL, JUSTIN WILLIAM

Sub DOB = 06/06/1986

Sub Sex = Male Weight = 245 Test = DUI Cit = NA

Dr. Lic. = ND/PEL865188

Lot No = 26021080A1

Cyl No = 17

Expiration Date = 10/05/2023

County = 53 Oper No. = 130755

Operator Signature
JUSTIN PELZL

Remarks: Inctall

Form 106-I8000