



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-004184</b>	Instrument Location <b>Divide Co. So Intake/Processing</b>
Reason for Install/Repair	
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <b>6/2/2022</b>
-------------------------------	-------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By 	Certified Date <b>07/11/2022</b>
---------------------------	-------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000           SN 80-004184  
Location = CROS                   8164.14.00 09/16  
06/02/2022                               15:08

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Dry  
Standard Value?            0.080  
Standard Lot #?            26021080A1  
Standard Cyl #?            24  
Standard Expiration?       10/05/2023  
Oper No?                    131836

Flow Cal. Date:            07/30/2010  
Slope                       804  
Intercept                  -925008

IR Calibration Date:       12/07/2015

	3um	9um
0th Coef(*100):	-15679	-18680
1st Coef(*100):	263558	136948
2nd Coef(*100):	2594	1146
H2O adj(mg/l*10k):	405	327

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
SHAWN HOLM

Remarks:

The undersigned, having custody of the original record, certifies that the information hereon is a true and correct copy of the original document maintained as part of the activities of this office.

Form 106-I8000

  
Signature

6/2/22  
Date


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004184  
Location = CROS      8164.14.00 09/16  
06/02/2022      15:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Gas	0.081	15:11
03 Room Air	0.000	15:11
04 Std. Gas	0.082	15:11
05 Room Air	0.000	15:12
06 Std. Gas	0.082	15:12
07 Room Air	0.000	15:13


Lot No = 26021080A1  
Cyl No = 24  
Exp Date = 10/05/2023  
County = 12      Oper No. = 131836

  
\_\_\_\_\_  
Operator Signature  
SHAWN HOLM

Remarks:

Form 106-I8000

The undersigned, having custody of the original record, certifies that the information hereon is a true and correct copy of the original document maintained as part of the activities of this office.

  
\_\_\_\_\_  
Signature

6/2/22  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004184  
Location = CROS      8164.14.00 09/16  
06/02/2022      15:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	15:15
02 Room Air	0.000	15:15

\*Invalid Test  
Inhibited - RFI

Lot No = 26021080A1  
Cyl No = 24  
Exp Date = 10/05/2023  
County = 12      Oper No. = 131836




Operator Signature  
SHAWN HOLM

Remarks:

Form 106-I8000

The undersigned, having custody of the original record, certifies that the information hereon is a true and correct copy of the original document maintained as part of the activities of this office.

      6/2/22  
Signature      Date