



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80 - 006490</b>	Instrument Location <b>Missotash</b>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

*Note: All forms were received by the Crime Laboratory on 20 June 22. 20 July 22 JPL*

Field Inspector Signature 	Date <b>6/17/22</b>
-------------------------------	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <b>Janelle Purtscheller</b>	Certified Date <b>6/22/22</b>
--	----------------------------------

*Note: AJP I noted the location code was wrong on the attached printouts. I contacted Steve Delzer on 6/22/22 to change the instrument location code from TOLX to ASHL. Steve verified no subject tests were run on the #59281 instrument since the instrument install by Mark Pearson on 6/17/22 until the location code was changed on 6/22/22.*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006490  
Location = TOXL                    8164.14.00 09/16  
06/17/2022                         18:53

*Actual Location  
ASHL 6/22/22 gup*

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                 No  
Display Volume?                     No  
Disable On Memfull?                Yes  
# of Print Copies?                  1  
Select Std (D/W/I)?                 Dry  
Standard Value?                     0.080  
Standard Lot #?                      02621080A1  
Standard Cyl #?                      36  
Standard Expiration?                03/05/2023  
Oper No?                              130117

Flow Cal. Date:                      09/04/2020  
Slope                                 649  
Intercept                            -598817

IR Calibration Date:                09/04/2020

	3um	9um
0th Coef(*100):	-25190	-16392
1st Coef(*100):	271645	136906
2nd Coef(*100):	2533	1222
H2O adj(mg/l*10k):	460	309

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MARK PEARSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006490  
Location = TOXL      8164.14.00 09/16  
06/17/2022      18:57

*Actual location  
ASHL 6/22/22 gen*

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:58
02 Std. Gas	0.077	18:58
03 Room Air	0.000	18:59
04 Std. Gas	0.078	18:59
05 Room Air	0.000	19:00
06 Std. Gas	0.078	19:00
07 Room Air	0.000	19:01

Lot No = 02621080A1  
Cyl No = 36  
Exp Date = 03/05/2023  
County = 26      Oper No. = 130117



Operator Signature  
MARK PEARSON

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006490  
Location = TOXL 8164.14.00 09/16  
06/17/2022 18:54  
*Actual Location 6/22/22 JWP*

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	18:56
02 Room Air	0.000	18:57

\*Invalid Test  
Inhibited - RFI

Lot No = 02621080A1  
Cyl No = 36  
Exp Date = 03/05/2023  
County = 26                                  Oper No. = 130117

  
\_\_\_\_\_  
Operator Signature  
MARK PEARSON

Remarks: