

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004206

Location: TOXL

A. General Setup and Checks:

- 1.  Diagnostic Tests Pass and Instrument in Ready Mode
- 2.  Breath tube heated
- 3.  Date, time and location code (Level 2,E). Re-set if necessary.
- 4.  Print test (Level 1,P). Sign and attach test record.
- 5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 883 psi Regulator: 850 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

a.  Low AC (Level 1,C): Known Value  $\leq$  0.03 AC: 0.020 AC

Sim. Ser #: MP3070 Lot #: 220P41184 Exp. Date: 3/1/24

b.  High AC (Level 1,C): Known Value  $\geq$  0.25 AC: 0.400 AC

Sim. Ser #: MP3065 Lot #: 202103E Exp. Date: 3/24/23

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone

Sim. Ser #: DR38478 Lot #: 1CS68 Exp. Date: NA

Display reads "Interferent Detect" 1CS7

3.  RFI Check (CMS Mode)  Display reads "RFI Detect"

4.  Dry Calibration Check (Level 1,C): Known Value 0.08 AC

Gas Cyl Lot #: 34921080A1 Cyl #: 44 Exp. Date: 2/5/24

Test 1 0.080 Test 4 0.079 Test 7 0.080

Test 2 0.080 Test 5 0.079 Test 8 0.078

Test 3 0.079 Test 6 0.078 Test 9 0.078

Average 0.079

C. Remarks/Maintenance: Annual inspection.

Replaced breath tube + simulator return o-rings.  
Note: initial high + low Ac checks were low due to a bad simulator solution.

Instrument is acceptable to be used in the field.

Jamille Perichiller  
Breath Analyst Signature

06/30/2022  
Date

[Signature]  
Reviewed by

07/01/22  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/30/2022    10:51

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                  1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    24119080A1  
Standard Cyl #?                    5  
Standard Expiration?                11/05/2021  
Oper No?                            132128

Flow Cal. Date:                    06/10/2009  
Slope                                690  
Intercept                            -645883

IR Calibration Date:                04/09/2015  
    3um                                9um

-----  
0th Coef(\*100):            -24157                                -27270  
1st Coef(\*100):            262661                                128882  
2nd Coef(\*100):            2137                                    1256  
H2O adj(mg/l\*10k):        491                                    640

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Print test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
06/30/2022      10:56

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:57
02 Std. Sol.	0.016	10:58
03 Room Air	0.000	10:58
04 Std. Sol.	0.017	10:59
05 Room Air	0.000	10:59
06 Std. Sol.	0.017	11:00
07 Room Air	0.000	11:01

08 Sim Temp = 34.0°C

Simul Ser No = MP3070  
Std Sol No = 22CP41184  
County = 08      Oper No. = 132128

*Janelle Portscheller*  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *#1 Low AC check 0.020*  
*Bad simulator solution*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/30/2022                            12:10

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:11
02 Std. Sol.	0.019	12:12
03 Room Air	0.000	12:12
04 Std. Sol.	0.019	12:13
05 Room Air	0.000	12:13
06 Std. Sol.	0.019	12:14
07 Room Air	0.000	12:15

08 Sim Temp = 34.0°C

Simul Ser No = MP3060  
Std Sol No = 202201A  
County = 08                            Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

#2  
Remarks: Low Act check 0.020

New Simulator ~~Sott~~ Solution.

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/30/2022                            11:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:03
02 Std. Sol.	0.387	11:04
03 Room Air	0.000	11:05
04 Std. Sol.	0.389	11:05
05 Room Air	0.000	11:06
06 Std. Sol.	0.390	11:07
07 Room Air	0.000	11:07

08 Sim Temp = 34.0°C

Simul Ser No = MP3065  
Std Sol No = 202103E  
County = 08                            Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: #1  
High AC check 0.400  
Bad simulator solution

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/30/2022                            12:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:16
02 Std. Sol.	0.295	12:17
03 Room Air	0.000	12:18
04 Std. Sol.	0.297	12:19
05 Room Air	0.000	12:19
06 Std. Sol.	0.298	12:20
07 Room Air	0.000	12:20

08 Sim Temp = 34.0°C

Simul Ser No = MP3064  
Std Sol No = 202012A  
County = 08                            Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *High AC check 0.300*  
*New simulator solution.*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-004206  
Location = TOXL                              8164.14.00 09/16  
06/30/2022    11:18

Test	AC	Time
01 Room Air	0.000	11:19
02 *Subject Test	INT*	11:19
03 Room Air	0.000	11:20

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                              Weight = 150  
Test = DUI                                      Cit = INTERFERENT CK  
Dr. Lic. = ND/DIS821456  
Lot No = 34921080A1  
Cyl No = 44  
Expiration Date = 02/05/2024  
County = 08                                      Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Interferent check*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/30/2022                            11:21

Test	AC	Time
01 Diagnostic	OK	11:21
02 Room Air	0.000	11:22
03 *Subject Test	RFI*	11:23
04 Room Air	0.000	11:24

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHK  
Dr. Lic. = ND/DIS821456  
Lot No = 34921080A1  
Cyl No = 44  
Expiration Date = 02/05/2024  
County = 08                            Oper No. = 132128

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *RFI check*




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                      8164.14.00 09/16  
06/30/2022                                      11:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:24
02 Std. Gas	0.080	11:25
03 Room Air	0.000	11:25
04 Std. Gas	0.080	11:26
05 Room Air	0.000	11:26
06 Std. Gas	0.079	11:26
07 Room Air	0.000	11:27

Lot No = 34921080A1  
Cyl No = 44  
Exp Date = 02/05/2024  
County = 08                                      Oper No. = 132128



Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check #1-3*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
06/30/2022      11:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:28
02 Std. Gas	0.079	11:28
03 Room Air	0.000	11:29
04 Std. Gas	0.079	11:29
05 Room Air	0.000	11:29
06 Std. Gas	0.078	11:30
07 Room Air	0.000	11:30

Lot No = 34921080A1  
Cyl No = 44  
Exp Date = 02/05/2024  
County = 08      Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check # 4-6*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                      8164.14.00 09/16  
06/30/2022                              11:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:32
02 Std. Gas	0.080	11:32
03 Room Air	0.000	11:33
04 Std. Gas	0.078	11:33
05 Room Air	0.000	11:34
06 Std. Gas	0.078	11:34
07 Room Air	0.000	11:34

Lot No = 34921080A1  
Cyl No = 44  
Exp Date = 02/05/2024  
County = 08                              Oper No. = 132128

  
\_\_\_\_\_  
Operator Signature  
JANELLE PORTSCHELLER

Remarks: *Dry Calibration Check # 7-9*