

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location			
80-004944 Will			
Reason for Install/Repair			
☐ Install After Receiving From Crime Laboratory 💢 Install After Location Change			
Other (Specify)			
Check When Done:	Malane ang an Maddalagay ang panggan naga panggangan ang ang ang ang ang ang ang an		
1. Surge Protector Installed/Property Grounded.			
Z 2. Telephone Line Connected to Intoxilyzer® 8000.			
🔀 3. Breath Tube Heated.			
X4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
5 € 6. Run Tests:			
🗖 A. Print Test (Level1, Function P).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test	).		
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer	® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (104-G), Print Test, ACA Test, and RFI Test.	(SFN59281, Form		
Field Inspector Signature  HP 21   Date	23		
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alc the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Labora County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.			
Reviewed/Certified By Certified	i Date		
	11122		

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = WILL 8164.14.00 09/16
01/03/2023 13:19

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 26021080A1 Standard Cyl #? 13 Standard Expiration? 10/05/2023 Oper No? 130739

Flow Cal. Date: 07/01/2011 Slope 685 Intercept -634740

IR Calibration Date: 07/01/2011

3um
9um

Oth Coef(\*100): -22591 -17234
1st Coef(\*100): 265915 132691
2nd Coef(\*100): 2841 1562
H2O adj(mg/l\*10k): 559 425

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator/Signature
MICHELLE ROMANS

Remarks:

Print test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
01/03/2023

Alcohol Analyzer SN 80-004944 8164.14.00 09/16 13:20

## DRY CAL CHECK

Test		AC	Time
01 Room	n Air	0.000	13:21
02 Std.	Gas	0.083	13:21
03 Room	ı Air	0.000	13:22
04 Std.	Gas	0.083	13:22
05 Room	n Air	0.000	13:23
06 Std.	Gas	0.082	13:23
07 Room	a Air	0.000	13:24

Lot No = 26021080A1

Cyl No = 13

Exp Date = 10/05/2023

County = 53

Oper No. = 130739

Operator Signature MICHELLE ROMANS

Remarks:

ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
01/03/2023

Alcohol Analyzer

SN 80-004944
8164.14.00 09/16
13:24

Te	est	AC	Time
	Room Air	RFI*	13:26
	Room Air	0.000	13:27

\*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI

Sub DOB = 01/01/2023

Sub Sex = Male

Weight = NA

Test = DUI

Cit = NA

Dr. Lic. = ND/NA

Lot No = 26021080A1

Cyl No = 13

Expiration Date = 10/05/2023

County = 53

Oper No. = 130739

Operator gignature MICHELLE ROMANS

Remarks:

RFI test

Form 106-I8000