

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 30 - 60 4939	Instrument Location Werd Compty Ja: 1			
Reason for Install/Repair				
☑ Install After Receiving From Crime Laboratory	er Location Change			
Other (Specify)				
Check When Done:				
1 Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
☑ 3 Breath Tube Heated.				
4 Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5, Scan/Enter Gas Cylinder Information (Level 1, Function S).				
🗹 6. Bun Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
S Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9, File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date		
ale		61/13/2023		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By		$\left(\right)$	Certified Date
	1 che	X	IJAn 23
	~ into		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004939 Location = WARD 8164.14.00 09/16 01/13/2023 17:16 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 # OI Princ Copies.Select Std (D/W/I)?DryStandard Value?0.080 Standard Lot #? Standard Cyl #? 34921080A1 33 Standard Expiration? 02/05/2024 Oper No? 131334 01/18/2010 Flow Cal. Date: Slope 724 Intercept -779308 IR Calibration Date: 04/05/2017 3um 9um _____ 0th Coef(*100): -24406 -20805 lst Coef(*100): 261875 2nd Coef(*100): 2828 135290 1427 H2O adj(mg/l*10k): 563393

282

Operator Signature JOSHUA GUDVANGEN


Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004939
Location = WARD	8164.14.00 09/16
01/13/2023	17:25

DRY CAI	L CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.085 0.000 0.085 0.000 0.085 0.000	17:25 17:25 17:26 17:26 17:27 17:27 17:27
Lot No = 34921080A1 Cyl No = 33 Exp Date = 02/05/2024 County = 51	4 Oper No. :	= 131334



Operator Signature JOSHUA GUDVANGEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004939 Location = WARD 8164.14.00 09/16 01/13/2023 17:35 Test AC Time 01 Room Air RFI* 17:37 02 Room Air 0.000 17:37 *Invalid Test Inhibited - RFI Sub Name = RFI TEST, RFI TEST NA Sub DOB = 01/01/1900Sub Sex = Unknown Weight = NA Test = OTHCit = NA Dr. Lic. = ND/NALot No = 34921080A1Cyl No = 33Expiration Date = 02/05/2024County = 51Oper No. = 131334

282

Operator Signature JOSHUA GUDVANGEN

Remarks:

Form 106-18000