

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 1188	Instrument Location Ward County Ja: 1		
Reason for Install/Repair			
✓ Install After Receiving From Crime Laboratory ☐ Install After Location Change			
Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
Z. Telephone Line Connected to Intoxilyzer® 8000.			
3/ Breath Tube Heated.			
4/ Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.			
Field Inspector Signature		Date	
lesson		01/2412023	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from			
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.			
Reviewed/Certified By		Certified Date	
		asimondo	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = WARD 8164.14.00 09/16
01/24/2023 16:13

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? NoDisplay Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34921080A1 Standard Cyl #? 37 Standard Expiration? 02/05/2024 Oper No? 131334

Flow Cal. Date: 04/30/2020 Slope 645 Intercept -637424

IR Calibration Date: 05/06/2020

3um 9um

Oth Coef(*100): -30092 -33449

1st Coef(*100): 293453 143538

2nd Coef(*100): 151 602

H2O adj(mg/l*10k): 425 462

****** Printer Test End ********

Operator Signature
JOSHUA GUDVANGEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WARD
01/24/2023

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 16:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:18
02 Std. Gas	0.080	16:18
03 Room Air	0.000	16:18
04 Std. Gas	0.079	16:19
05 Room Air	0.000	16:19
06 Std. Gas	0.079	16:20
07 Room Air	0.000	16:20

Lot No = 34921080A1

Cyl No = 37

Exp Date = 02/05/2024

County = 51

Oper No. = 131334

Operator Signature
JOSHUA GUDVANGEN

282

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = WARD 8164.14.00 09/16
01/24/2023 16:21

Test	AC	Time
01 Room Air	0.000	16:22
02 *Subject Test	RFI*	16:24
03 Room Air	0.000	16:24

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST NA

Sub DOB = 01/01/2023

Sub Sex = Unknown Weight = NA Test = OTH Cit = NA

Dr. Lic. = ND/NA
Lot No = 34921080A1

Cyl No = 37

Expiration Date = 02/05/2024

> Operator Signature JOSHUA GUDVANGEN

282

Remarks:

Form 106-I8000