



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number: <u>80-003059</u>	Instrument Location: <u>WAPD</u>
Reason for Install/Repair: <input checked="" type="checkbox"/> Install After Receiving from Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify)	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location). (Level 2, Function E)
- 5. Scan/Enter Gas Cylinder Information. (Level 1, Function S)
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode – Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN 50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN 50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN 59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Troy Heseck
 Field Inspector Signature

5-30-11
 Date

Debi Kashin
 Reviewed By

01 June 2011
 Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
05/30/2011 11:05

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 05911080A1
Standard Cyl #? 34
Standard Expiration? 02/01/2013
Oper No? 049524

Flow Cal. Date: 04/07/2008
Slope 713
Intercept -692349

IR Calibration Date: 06/04/2008

	3um	9um
0th Coef(*100):	-15808	-12583
1st Coef(*100):	279401	139580
2nd Coef(*100):	1046	760
H2O adj(mg/l*10k):	478	208

***** Printer Test End *****


Operator Signature
TROY A. HISCHER

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
05/30/2011 11:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:06
02 Std. Gas	0.078	11:07
03 Room Air	0.000	11:07
04 Std. Gas	0.078	11:08
05 Room Air	0.000	11:08
06 Std. Gas	0.078	11:09
07 Room Air	0.000	11:09

Lot No = 05911080A1
Cyl No = 34
Exp Date = 02/01/2013
County = 09 Oper No. = 049524



Operator Signature
TROY A. HISCHER

Remarks:

Form 106-I8000

80-003059
P. 131 of 132
AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
05/30/2011 11:13

Test	AC	Time
01 Diagnostic	OK	11:14
02 Room Air	RFI*	11:14
03 Room Air	RFI*	11:14

*Invalid Test
Inhibited - RFI

Sub Name = RFI, RFI RFI
Sub DOB = 01/01/1990
Sub Sex = Male Weight = 200
Test = DUI Cit = 000000001
Dr. Lic. = ND/RFI000000
Lot No = 05911080A1
Cyl No = 34
Expiration Date = 02/01/2013
County = 09 Oper No. = 049524

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
TROY A. HISCHER

Remarks:

Form 106-I8000

80-003059
p. 132 of 132
VM