



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 2667

A. General Setup and Checks:

- 1. Initial start up
  - Power On
  - 5-minute countdown during initial warm-up
  - 1-minute countdown to begin test after standby mode
- 2.  Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
  - Level 1 (Operator)
  - Level 2 (Field Inspector)
  - Level 3 (Lab Personnel)
- 4.  Printer Setup (Level 3,S,P)
- 5.  Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

- 1.  Internal Printer Test (Level 1,P). Sign and attach record.
  - Paper Feed
- 2.  External Printer Test (Level 1,P). Sign and attach record.
- 3. Barcode Reader (Level 1,B) (Tone sounds after each)
  - Operator ID
  - Invalid ID
  - ND License
  - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

- 1.  Continuous Air Blank (Level 2,A)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 1 of 4

2.  Diagnostic (Level 2,D). Sign and attach test record.

3.  Prelim Data Entry (Level 2,E)

4. Ser Number/Version Info (Level 2,V):

Inst Ser#: 80-002667 DSP Ver: 1362.14

Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12508 s: 0 C: 47.00

9: 12488 s: 1 B: 45

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0 → 1

2. Function (F) Keys:

F5 – Tone

F9 - Fan *← does not turn on.*

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 954 CF: 1.06

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 74 F: 0.000L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 45.0 I: 0.0

B: 44.9 V: 0.0

Cell Temperature (C) between  $47 \pm 0.02$  °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 349 psi Regulator: 350 psi

Display and Regulator  $\pm 50$  psi of each other Yes or No

Gas tank tare necessary? Yes or No

If Yes, display readings after tare (Level 3,M,C,G):

Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

7.  General Setup (Level 3,S,G)

8.  Configure Start Test (Level 3,S,T)

9.  Network (Level 3,S,N)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 2 of 4

E. Tests (Sign and attach all test records):

1.  Configure wet bath simulator for the following test (Level 1,S)  
 Zero AC-Wet Calibration Check (Level 1,C):  
 Sim ser #: MP3066 Lot #: 32833023 Exp date: N/A  
 Known value: 0.000 AC  
 Results are Zero AC
2.  Configure wet bath simulator for the following test (Level 1,S)  
 Low AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\leq 0.03$  AC: 0.020 AC  
 Sim ser #: MP5318 Lot #: 202201A Exp date: 11Jan24  
 Results  $\pm 0.005$  of known AC
3.  Configure wet bath simulator for the following test (Level 1,S)  
 High AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\geq 0.25$  AC: 0.300 AC  
 Sim ser #: MP3067 Lot #: 202201F Exp date: 18Jan24  
 Results  $\pm 5\%$  of known AC
4.  Configure dry gas standard for the following tests (Level 1,S)  
 Known Value: 0.080 AC  
 Gas Cylinder Lot #: 26021080A1 Cylinder #: 24 Exp. Date: 10/5/23
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7345 Lot #: 1CS7 Exp. Date: N/A  
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7349 Lot #: ACE3 Exp. Date: N/A  
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):  
 Sim. Ser #: DR7347 Lot #: MET3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):  
 Sim. Ser #: DR7350 Lot #: IPA3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
9.  Ambient Fail (Level 1,B): Display read and test record correct.
10.  Invalid Sample (Level 1,B): Display read and test record correct.
11.  Sequence Aborted (Level 1,B): Display read test record correct.
12.  ABA (Level 1,B): Test record correct.

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 3 of 4

- 13.  Calibration check out of tolerance (CMS mode): Display read and test record correct.
- 14.  RFI Detected (B,C, or CMS mode): Display read and test record correct.
- 15.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 16.  Deficient Samples (CMS mode): Display read and test record correct.
- 17.  Improper Sample (CMS mode): Display read and test record correct.
- 18.  Difference Too Great (CMS mode): Display read and test record correct.
- 19.  CMS (CMS mode): Test record correct.
- 20. Stability Test (Level 3,T,S):  
 Ave: 0.0800 Std Dev: 0.0005 Rel Std Dev: 0.58%

Maintenance/Repairs/Notes: N/A

---



---



---

Instrument is acceptable to be used in the field. Yes or No

If No, state reason(s) why:  
Classroom only instrument b/c fan motor  
does not turn on.

---



---

*[Signature]*  
 Inspector Signature

24 May 2023  
 Date

*[Signature]*  
 Reviewer

25 May 2023  
 Date

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 4 of 4

Intoxilyzer Test Record and Checklist  
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-002667  
Location = TOXL 8164.14.00 09/16  
05/24/2023 12:01

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMN

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMN

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMN

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef

Cu  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/1)? Dry  
Standard Value? 0.080  
Standard Lot #? 34921080A1  
Standard Cyl #? 33  
Standard Expiration? 02/05/2024  
Oper No? 133237

Flow Cal. Date: 04/19/2021  
Slope 615  
Intercept 113911

IR Calibration Date: 05/08/2017  
3um 9um

0th Coef(\*100): -16077 -8130

1st Coef(\*100): 268518 134907

2nd Coef(\*100): 2115 1175

H2O adj(mg/l\*10k): 380 220

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks:  
Internal Print  
Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-002667  
Location = TOXL                    8164.14.00 09/16  
05/24/2023                            12:03

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34921080A1  
Standard Cyl #?            33  
Standard Expiration?       02/05/2024  
Oper No?                    133237

Flow Cal. Date:            04/19/2021  
Slope                        615  
Intercept                    113911

IR Calibration Date:        05/08/2017  
                                  3um                    9um

	3um	9um
0th Coef(*100):	-16077	-8130
1st Coef(*100):	268518	134907
2nd Coef(*100):	2115	1175
H2O adj(mg/l*10k):	380	220

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks: External Print Test


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                    8164.14.00 09/16  
05/24/2023                            12:03

Test	AC	Time
01 Room Air	0.000	12:04
02 *Subject Test	REF*	12:04
03 Room Air	0.000	12:05

\*Subject Test Refused

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                            Cit = REFUSAL  
Dr. Lic. = ND/TES989643  
Lot No = 34921080A1  
Cyl No = 33  
Expiration Date = 02/05/2024  
County = 08                          Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Refusal*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:02

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks:

Diagnostic

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:09

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:09
02 Std. Sol.	0.000	12:10
03 Room Air	0.000	12:11
04 Std. Sol.	0.000	12:11
05 Room Air	0.000	12:12
06 Std. Sol.	0.000	12:12
07 Room Air	0.000	12:13

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 33023  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: *Zero AC Check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:15

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:16
02 Std. Sol.	0.023	12:17
03 Room Air	0.000	12:17
04 Std. Sol.	0.023	12:18
05 Room Air	0.000	12:18
06 Std. Sol.	0.022	12:19
07 Room Air	0.000	12:20

08 Sim Temp = 34.0°C

Simul Ser No = MP5318  
Std Sol No = 202201A  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: LOW AC - 0.020Ac

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:21
02 Std. Sol.	0.300	12:22
03 Room Air	0.000	12:22
04 Std. Sol.	0.303	12:23
05 Room Air	0.000	12:24
06 Std. Sol.	0.304	12:24
07 Room Air	0.000	12:25

08 Sim Temp = 34.0°C

Simul Ser No = MP3067  
Std Sol No = 202201F  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: High AC - 0.300Ac

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:27

Test	AC	Time
01 Room Air	0.000	12:28
02 *Subject Test	INT*	12:28
03 Room Air	0.000	12:29

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
ETOH + Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                    8164.14.00 09/16  
05/24/2023                            12:29

Test	AC	Time
01 Room Air	0.000	12:30
02 *Subject Test	INT*	12:30
03 Room Air	0.000	12:30

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:31

Test	AC	Time
01 Room Air	0.000	12:32
02 *Subject Test	INT*	12:32
03 Room Air	0.000	12:33

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
Methanol

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:33

Test	AC	Time
01 Room Air	0.000	12:34
02 *Subject Test	INT*	12:34
03 Room Air	0.000	12:35

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
IPA

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                    8164.14.00 09/16  
05/24/2023                            12:35

Test	AC	Time
01 Room Air	AMB*	12:36
02 Room Air	0.000	12:36

\*Invalid Test  
Check Ambient Conditions

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                            Cit = AMBIENT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks:

Ambient Fail

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:39

Test	AC	Time
01 Room Air	0.000	12:40
02 Invalid Sample	X.XXX	12:40
03 Room Air	0.000	12:41

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = INVALID  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Invalid Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      12:41

Test	AC	Time
01 Room Air	0.000	12:42
02 Reference	ABT*	12:42
03 Room Air	0.000	12:43

\*Sequence Aborted

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = SEQUENCEABORTED  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Sequence Aborted*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      13:03

Test	AC	Time
01 Room Air	0.000	13:03
02 Subject Test 1	0.000	13:04
03 Room Air	0.000	13:04
04 Reported AC	0.000	13:04

No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = ABA TEST  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: ABA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                  Alcohol Analyzer  
North Dakota Model 8000                  SN 80-002667  
Location = TOXL                  8164.14.00 09/16  
05/24/2023                                  13:05

Test	AC	Time
01 Diagnostic	OK	13:05
02 Room Air	0.000	13:06
03 Subject Test 1	0.000	13:06
04 Room Air	0.000	13:09
05 Std. Gas	0.050*	13:10
06 Room Air	0.000	13:10

\*Cal Check Out of Tolerance

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                  Weight = NA  
Test = OTH                                  Cit = CAL CHECK OOT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                                  Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Calibration Check  
Out of Tolerance

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                 8164.14.00 09/16  
05/24/2023                         13:12

Test	AC	Time
01 Room Air	0.000	13:13
02 *Subject Test	RFI*	13:13
03 Room Air	0.000	13:13

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = RFI CHECK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: *RFI check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-002667  
Location = TOXL                                      8164.14.00 09/16  
05/24/2023    13:14

Test	AC	Time
01 Diagnostic	OK	13:15
02 Room Air	0.000	13:15
03 *Subject Test	0.000*	13:18
04 Room Air	0.000	13:20
05 Std. Gas	0.080	13:22
06 Room Air	0.000	13:23
07 *Subject Test	0.000*	13:26
08 Room Air	0.000	13:28

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                      Weight = NA  
Test = OTH    Cit = DEFICIENT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                                      Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002667  
Location = TOXL                    8164.14.00 09/16  
05/24/2023                            13:29

Test	AC	Time
01 Diagnostic	OK	13:30
02 Room Air	IPS*	13:30
03 Room Air	0.000	13:31

\*Improper Sample

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = IMPROPER SAMPLE  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Improper Sample*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-002667  
Location = TOXL                                      8164.14.00 09/16  
05/24/2023    13:31

Test	AC	Time
01 Diagnostic	OK	13:32
02 Room Air	0.000	13:32
03 Subject Test 1	0.042*	13:33
04 Room Air	0.000	13:35
05 Std. Gas	0.080	13:36
06 Room Air	0.000	13:37
07 Subject Test 2	0.000*	13:38
08 Room Air	0.000	13:39

\*Difference Too Great

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                      Weight = NA  
Test = OTH                                      Cit = DIFF TOO GREAT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 24  
Expiration Date = 10/05/2023  
County = 08                                      Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000





Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = TOXL      8164.14.00 09/16  
05/24/2023      13:49

Test	AC	Time
01 Room Air	0.000	13:49
02 Std. Gas	0.080	13:50
03 Room Air	0.000	13:50
04 Std. Gas	0.080	13:50
05 Room Air	0.000	13:51
06 Std. Gas	0.080	13:51
07 Room Air	0.000	13:52
08 Std. Gas	0.079	13:52
09 Room Air	0.000	13:53
10 Std. Gas	0.080	13:53
11 Room Air	0.000	13:53
12 Std. Gas	0.081	13:54
13 Room Air	0.000	13:54
14 Std. Gas	0.080	13:55
15 Room Air	0.000	13:55
16 Std. Gas	0.080	13:56
17 Room Air	0.000	13:56
18 Std. Gas	0.079	13:56
19 Room Air	0.000	13:57
20 Std. Gas	0.080	13:57
21 Room Air	0.000	13:58
22 Std. Gas	0.080	13:58
23 Room Air	0.000	13:59
24 Std. Gas	0.080	13:59
25 Room Air	0.000	13:59
26 Std. Gas	0.080	14:00
27 Room Air	0.000	14:00
28 Std. Gas	0.081	14:01
29 Room Air	0.000	14:01
30 Std. Gas	0.080	14:01
31 Room Air	0.000	14:02
32 Std. Gas	0.080	14:02
33 Room Air	0.000	14:03
34 Std. Gas	0.080	14:03
35 Room Air	0.000	14:04
36 Std. Gas	0.080	14:04
37 Room Air	0.000	14:04
38 Std. Gas	0.080	14:05
39 Room Air	0.000	14:05
40 Std. Gas	0.080	14:06
41 Room Air	0.000	14:06
42 Cal Check Stats:		
43 Average	0.0800	
44 Std Dev	0.0005	
45 Rel Std Dev	0.58%	



Operator Signature  
ANNA NAREHOOD

Remarks:

Stability Test  


Form 106-I8000