

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 70 95	Instrument Location	
Reason for Install/Repair	1/12	
	er Location Change	
φ	er Education Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loca	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B., ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode of	or Level 1, Function B or C; Key Radio Dur	ring Test).
7. Repair and/or Maintenance Performed (if any):	L Field Insp	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN	50496, Form 120-G) and Place it by the In	toxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120		
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature Field Kolik		Date 2/16/23
		/ /
Crime Laboratory Use Only	he would far the analysis of brooth to date	rming alcohol concentration from
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record or	be used for the analysis of breath to deter In file at the Office of Attorney General, Crin	me Laboratory Division, in the
County of Burleigh, North Dakota, is certified to be a true and correct	copy of the documents received.	

Reviewed/Certified By

Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = ADAM 8164.16.00 09/18
02/16/2023 13:56

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34921080A1 29 Standard Expiration? 02/05/2024 Oper No? 132666

Flow Cal. Date: 09/27/2018 Slope 683

Intercept -637452

IR Calibration Date: 09/27/2018
3um 9um

Oth Coef(*100): -59629 -23097
1st Coef(*100): 252379 137307
2nd Coef(*100): 2489 1549
H2O adj(mg/l*10k): 1414 594

******* Printer Test End *******

Operator Signature JEDEDIAH KOHLER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ADAM
02/16/2023

Alcohol Analyzer SN 80-007095 8154.16.00 09/18 14:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:28
02 Std. Gas	0.082	14:29
03 Room Air	0.000	14:29
04 Std. Gas	0.082	14:29
05 Room Air	0.000	14:30
06 Std. Gas	0.082	14:30
07 Room Air	0.000	14:31

Lot No = 34921080A1

Cy1 No = 29

Exp Date = 02/05/2024

County = 01

Oper No. = 132666

Operator Signature JEDEDIAH KOHLER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = ADAM 02/16/2023

SN 80-007095 8164.16.00 09/18 14:33

Test	AC	Time
01 Room Air	RFI*	14:35
02 Room Air	0.000	14:35

*Invalid Test Inhibited - RFI

Sub Name = KOHLER, JEDEDIAH SCOTT

Sub DOB = 08/20/1987

Sub Sex = Male

Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/KOH872036

Lot No = 34921080A1

Cyl No = 29

Expiration Date = 02/05/2024

County = 01

Oper No. = 132666

Operator Signature JEDEDIAH KOHLER

Remarks:

Form 106-I8000