



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-000678</b>	Instrument Location <b>GFPD</b>
Reason for Install/Repair	
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <b>Moved Machine to UNPPD for Data File/Install back at GFPD</b>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i> #611	Date <b>03/07/23</b>
--	-------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <b>21 Mar 2023</b>
---	--------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006678  
Location = GFPD                8164.14.00 09/16  
03/17/2023                      11:50

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            01923080A3  
Standard Cyl #?            57  
Standard Expiration?       02/05/2025  
Oper No?                    132718

Flow Cal. Date:            08/22/2016  
Slope                      686  
Intercept                  -611596

IR Calibration Date:        08/22/2016  
                                    3um                      9um

-----  
0th Coef(\*100):            -15445                      -21650  
1st Coef(\*100):            251854                      131766  
2nd Coef(\*100):            3094                        1378  
H2O adj(mg/l\*10k):        549                         513

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
-----  
Operator Signature  
DYLAN BARBER

Remarks: *Print Test Install 03/17/23*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006678  
Location = GFPD      8164.14.00 09/16  
03/17/2023      11:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:51
02 Std. Gas	0.083	11:51
03 Room Air	0.000	11:52
04 Std. Gas	0.083	11:52
05 Room Air	0.000	11:53
06 Std. Gas	0.084	11:53
07 Room Air	0.000	11:54

Lot No = 01923080A3  
Cyl No = 57  
Exp Date = 02/05/2025  
County = 18      Oper No. = 132718

  
Operator Signature  
DYLAN BARBER

Remarks:

ACA Test

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006678  
Location = GFPD      8164.14.00 09/16  
03/17/2023      11:54

Test	AC	Time
01 Room Air	0.000	11:56
02 *Subject Test	RFI*	11:57
03 Room Air	0.000	11:57

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST NA  
Sub DOB = 01/01/2001  
Sub Sex = Male      Weight = 150  
Test = DUI      Cit = 00  
Dr. Lic. = ND/00  
Lot No = 01923080A3  
Cyl No = 57  
Expiration Date = 02/05/2025  
County = 18      Oper No. = 132718

  
Operator Signature  
DYLAN BARBER

Remarks: RFI TEST

Form 106-I8000