

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004201 Instrument Location 5
Reason for Install/Repair
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change
Other (Specify) Annual Inspection at Chine Laboratory
Check When Done:
1. Surge Protector Installed/Property Grounded.
2. Telephone Line Connected to Intoxilyzer® 8000.
☐ 3. Breath Tube Heated.
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
∑ 6. Run Tests:
A. Print Test (Level1, Function P).
B. ACA Test (Level 1, Function C).
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any): NONE
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.
Field Inspector Signature July Dill Date 04/17/2023
Crime Laboratory Use Only
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration fro the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.
Reviewed/Certified By Certified Date 17 Apr 2023

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004201 Location = CARS 8164.14.00 09/16 04/17/2023 08:47

******** Printer Test ********

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes

Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? Standard Lot #? Standard Cyl #? 0.080 26021080A1

10/05/2023

Standard Expiration? Oper No? 132318

04/25/2018 Flow Cal. Date: Slope 580 Intercept -164747

IR Calibration Date: 04/26/2018 3um 9um

Oth Coef(*100): -26833 1st Coef(*100): 272526 2nd Coef(*100): 3966 -28328 132392 1397 H20 adj(mq/l*10k): 561514

****** Printer Test End *******

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CARS
04/17/2023

Alcohol Analyzer SN 80-004201 8164.14.00 09/16 08:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:54
02 Std. Gas	0.081	08:54
03 Room Air	0.000	08:55
04 Std. Gas	0.081	08:55
05 Room Air	0.000	08:55
06 Std. Gas	0.080	08:56
07 Room Air	0.000	08:56

Lot No = 26021080A1

Cyl No = 6

Exp Date = 10/05/2023

County = 19

Oper No. = 132318

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Al North Dakota Model 8000 Location = CARS 81 04/17/2023

Alcohol Analyzer SN 80-004201 8164.14.00 09/16 08:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:49
02 Std. Gas	0.080	08:50
03 Room Air	0.000	08:50
04 Std. Gas	0.080	08:50
05 Room Air	0.000	08:51
06 Std. Gas	RFI*	08:51
07 Room Air	0.000	08:52

*Invalid Test Inhibited - RFI

Lot No = 26021080A1

Cyl No = 6

Exp Date = 10/05/2023

County = 19

Oper No. = 132318

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000