



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4198

A. General Setup and Checks:

- 1. Initial start up
  - Power On
  - 5-minute countdown during initial warm-up
  - 1-minute countdown to begin test after standby mode
- 2.  Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
  - Level 1 (Operator)
  - Level 2 (Field Inspector)
  - Level 3 (Lab Personnel)
- 4.  Printer Setup (Level 3,S,P)
- 5.  Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

- 1.  Internal Printer Test (Level 1,P). Sign and attach record.
  - Paper Feed
- 2.  External Printer Test (Level 1,P). Sign and attach record.
- 3. Barcode Reader (Level 1,B) (Tone sounds after each)
  - Operator ID
  - Invalid ID
  - ND License
  - R key: Refusal. Sign and attach record.

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 1
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/20/2023
UNCONTROLLED WHEN PRINTED	Page 1 of 4

C. Level 2 Setup and Checks:

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Sign and attach test record.
- Prelim Data Entry (Level 2,E)
- Ser Number/Version Info (Level 2,V):

Inst Ser#: 80-004198 DSP Ver: 1362.14  
Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

- DVM Monitor (Level 3,D,D):

3: 12762 s: 1 C: 47.00  
9: 12733 s: 2 B: 45

RFI indicator with keyed radio: 071

- Function (F) Keys:

F5 – Tone                       F9 - Fan  
 F6 – Green LED                 F10 – Pump  
 F7 – Red LED                     F11 – Solenoid

- Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 950 CF: 1.06

- Flow Monitor (Level 3,D,F): Flow Sensor, R: 214 F: 0.000 L/s

- Temperature Monitor (Level 3,D,T): C: 47.0 A: 48.5 I: 0.0  
B: 45.0 V: 0.0

- Tank Monitor (Level 3,D,G):

Display: 688 psi Regulator: 700 psi  
Display and Regulator ± 50 psi of each other  Yes or No  
Gas tank tare necessary? Yes or  No  
If Yes, display readings after tare (Level 3,M,C,G):  
Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

- General Setup (Level 3,S,G)
- Configure Start Test (Level 3,S,T)
- Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

- Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP3066 Lot #: 33023 Exp date: N/A  
Known value: 0.000 AC

Results are Zero AC

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 1
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/20/2023
UNCONTROLLED WHEN PRINTED	Page 2 of 4

2.  Configure wet bath simulator for the following test (Level 1,S)  
 Low AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\leq 0.03$  AC: 0.020 AC  
 Sim ser #: MP5318 Lot #: 202201A Exp date: 11Jan24  
 Results  $\pm 0.005$  of known AC
3.  Configure wet bath simulator for the following test (Level 1,S)  
 High AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\geq 0.25$  AC: 0.300 AC  
 Sim ser #: MP6039 Lot #: 202201F Exp date: 18Jan24  
 Results  $\pm 5\%$  of known AC
4.  Configure dry gas standard for the following tests (Level 1,S)  
 Known Value: 0.080 AC  
 Gas Cylinder Lot #: 26021080A1 Cylinder #: 25 Exp. Date: 10/5/23
5. Ethanol + 0.05% Acetone (Level 1,B):  
 Known Value: 0.10% EtC + .05% ACE  
 Sim. Ser #: DR3841 Lot #: 1CS7 Exp. Date: NIA  
 Display reads "Interferent Detect"
6. Acetone (Level 1,B):  
 Known Value: 0.05% ACE  
 Sim. Ser #: DR7349 Lot #: ACE3 Exp. Date: NIA  
 Display reads "Interferent Detect"
7. Methanol (Level 1,B):  
 Known Value: 0.10% Methanol  
 Sim. Ser #: DR7347 Lot #: MET3 Exp. Date: NIA  
 Display Reads "Interferent Detect"
8. Isopropanol (Level 1,B):  
 Known Value: 0.10% IPA  
 Sim. Ser #: DR7350 Lot #: IPA3 Exp. Date: NIA  
 Display Reads "Interferent Detect"
9.  Ambient Fail (Level 1,B): Display read and test record correct.
10.  Invalid Sample (Level 1,B): Display read and test record correct.
11.  Sequence Aborted (Level 1,B): Display read test record correct.
12.  ABA (Level 1,B): Test record correct.
13.  Calibration check out of tolerance (CMS mode): Display read and test record correct.
14.  RFI Detected (B,C, or CMS mode): Display read and test record correct.
15.  Tank Pressure Below Minimum (CMS mode): Display reads correct.

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 1
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/20/2023
UNCONTROLLED WHEN PRINTED	Page 3 of 4



- 16.  Deficient Samples (CMS mode): Display read and test record correct.
- 17.  Improper Sample (CMS mode): Display read and test record correct.
- 18.  Difference Too Great (CMS mode): Display read and test record correct.
- 19.  CMS (CMS mode): Test record correct.
- 20. Stability Test (Level 3,T,S):  
 Ave: 0.0793 Std Dev: 0.006 Rel Std Dev: 0.73%

Maintenance/Repairs/Notes: N/A The internal print test, external print test, refusal and diagnostic all state 04/08/2023 because they were ran before C.3. The preliminary data entry.

Instrument is acceptable to be used in the field.  Yes or  No

If No, state reason(s) why:

---



---



---

  
 Inspector Signature

21 Apr 2023  
 Date

  
 Reviewer

21 Apr 2023  
 Date

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 1
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/20/2023
UNCONTROLLED WHEN PRINTED	Page 4 of 4

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004198  
Location = TOXL 8164.14.00 09/16  
04/08/2023 13:30

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMN

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMN

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMN

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMN


Current Test Mode: 0000  
Data Entry Mode: 0000  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Wet  
Standard Value? 0.000  
Standard Solution #? 33023  
Sim Serial #? MP3066  
Oper No? 133237

Flow Cal. Date: 01/09/2012  
Slope: 637  
Intercept: -467208

IR Calibration Date: 05/18/2016  
3um 9um

0th Coef(*100):	-19233	-23050
1st Coef(*100):	267031	132119
2nd Coef(*100):	2449	1142
H2O adj(ng/l*10k):	329	398

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks:  
Internal Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-004198  
Location = TOXL               8164.14.00 09/16  
04/08/2023                       13:38

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&\*()\_+?

Current Instrument Setup

Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?           No  
Disable On Memfull?       Yes  
# of Print Copies?         1  
Select Std (D/W/I)?       Wet  
Standard Value?           0.000  
Standard Solution #?       33023  
Sim Serial #?             MP3066  
Oper No?                  133237

Flow Cal. Date:            01/09/2012  
Slope                      637  
Intercept                 -467208

IR Calibration Date:       05/18/2016  
                            3um                   9um

0th Coef(*100):	-19233	-23050
1st Coef(*100):	267031	132119
2nd Coef(*100):	2449	1142
H2O adj(mg/l*10k):	329	398

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks: External Printer Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/08/2023                            13:39

Test	AC	Time
01 Room Air	0.000	13:40
02 *Subject Test	REF*	13:40
03 Room Air	0.000	13:41

\*Subject Test Refused

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = REFUSAL  
Dr. Lic. = ND/TES989643  
Simul Ser No = MP3066  
Std Sol No = 33023  
County = 08                                Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Refusal

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = TOXL                8164.14.00 09/16  
04/08/2023                        13:41

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      15:59

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:59
02 Std. Sol.	0.000	16:00
03 Room Air	0.000	16:00
04 Std. Sol.	0.000	16:01
05 Room Air	0.000	16:01
06 Std. Sol.	0.000	16:02
07 Room Air	0.000	16:03

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 33023  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: Zero AC Check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      16:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:04
02 Std. Sol.	0.018	16:05
03 Room Air	0.000	16:05
04 Std. Sol.	0.020	16:06
05 Room Air	0.000	16:06
06 Std. Sol.	0.020	16:07
07 Room Air	0.000	16:08

08 Sim Temp = 34.0°C

Simul Ser No = MP5318  
Std Sol No = 202201A  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Low Ae Check - 0.020Ac*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      16:08

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:09
02 Std. Sol.	0.297	16:10
03 Room Air	0.000	16:10
04 Std. Sol.	0.302	16:11
05 Room Air	0.000	16:11
06 Std. Sol.	0.301	16:12
07 Room Air	0.000	16:13

08 Sim Temp = 34.0°C

Simul Ser No = MP6039  
Std Sol No = 202201F  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: High AC Check - 0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004198
Location = TOXL	8164.14.00 09/16
04/20/2023	16:14

Test	AC	Time
01 Room Air	0.000	16:15
02 *Subject Test	INT*	16:15
03 Room Air	0.000	16:16

\*Invalid Test  
 Interferent Detected

Sub Name = TEST, DONOR2 NONE  
 Sub DOB = 07/25/1998  
 Sub Sex = Male                      Weight = NA  
 Test = OTH                      Cit = ETOH INTFT CK  
 Dr. Lic. = ND/TES989643  
 Lot No = 26021080A1  
 Cyl No = 25  
 Expiration Date = 10/05/2023  
 County = 08                      Oper No. = 133237

Operator Signature  
 ANNA NAREHOOD

Remarks: ETOH Interferent CK

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      16:16

Test	AC	Time
01 Room Air	0.000	16:17
02 *Subject Test	INT*	16:17
03 Room Air	0.000	16:18

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = ACETONE CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Acetone Interferent check

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      16:18

Test	AC	Time
01 Room Air	0.000	16:19
02 *Subject Test	INT*	16:19
03 Room Air	0.000	16:19

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = METOH CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: MeOH Interferent Check

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = TOXL      8164.14.00 09/16  
04/20/2023      16:20

Test	AC	Time
01 Room Air	0.000	16:20
02 *Subject Test	INT*	16:21
03 Room Air	0.000	16:21

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = IPA INTFT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: IPA Interferent Check

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/20/2023                            16:22

Test	AC	Time
01 Room Air	AMB*	16:22
02 Room Air	0.000	16:23

\*Invalid Test  
Check Ambient Conditions

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = AMBIENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: *Ambient Fail*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/20/2023                            16:25

Test	AC	Time
01 Room Air	0.000	16:26
02 Invalid Sample	X.XXX	16:26
03 Room Air	0.000	16:27

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = INVALID SAMPLE  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: *Invalid Sample*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = TOXL                8164.14.00 09/16  
04/20/2023                        16:27

Test	AC	Time
01 Room Air	0.000	16:28
02 *Subject Test	ABT*	16:28
03 Room Air	0.000	16:28

\*Sequence Aborted

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = SEQUENCE ABT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Sequence Aborted*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = TOXL                8164.14.00 09/16  
04/20/2023                        16:29

Test	AC	Time
01 Room Air	0.000	16:29
02 Subject Test 1	0.011	16:30
03 Room Air	0.000	16:30
04 Reported AC	0.011	16:30

No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = ABA  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: ABA TEST

Form 106-I8000





Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/21/2023                            08:53

Test	AC	Time
01 Diagnostic	OK	08:54
02 Room Air	0.000	08:54
03 *Subject Test	0.000*	08:57
04 Room Air	0.000	08:59
05 Std. Gas	0.079	09:01
06 Room Air	0.000	09:02
07 *Subject Test	0.000*	09:05
08 Room Air	0.000	09:06

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = DEFICIENT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/21/2023                            09:07

Test	AC	Time
01 Diagnostic	OK	09:08
02 Room Air	0.000	09:08
03 Subject Test 1	0.000	09:09
04 Room Air	IPS*	09:10
05 Room Air	0.000	09:10

\*Improper Sample

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = IMPROPER  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Improper Sample*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004198  
Location = TOXL                    8164.14.00 09/16  
04/21/2023                            09:11

Test	AC	Time
01 Diagnostic	OK	09:12
02 Room Air	0.000	09:12
03 Subject Test 1	0.000*	09:12
04 Room Air	0.000	09:15
05 Std. Gas	0.079	09:16
06 Room Air	0.000	09:17
07 Subject Test 2	0.042*	09:17
08 Room Air	0.000	09:18

\*Difference Too Great

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = DIFF TOO GREAT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-004198  
Location = TOXL                              8164.14.00 09/16  
04/21/2023    09:19


Test	AC	Time
01 Diagnostic	OK	09:19
02 Room Air	0.000	09:20
03 Subject Test 1	0.000	09:20
04 Room Air	0.000	09:22
05 Std. Gas	0.079	09:24
06 Room Air	0.000	09:25
07 Subject Test 2	0.000	09:25
08 Room Air	0.000	09:26
09 Reported AC	0.000	09:20

Difference OK  
No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                              Weight = NA  
Test = OTH                                      Cit = CMS  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                              Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: CMS test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = TOXL                8164.14.00 09/16  
04/21/2023                      09:28

Test	AC	Time
01 Room Air	0.000	09:28
02 Std. Gas	0.079	09:28
03 Room Air	0.000	09:29
04 Std. Gas	0.079	09:29
05 Room Air	0.000	09:30
06 Std. Gas	0.079	09:30
07 Room Air	0.000	09:31
08 Std. Gas	0.078	09:31
09 Room Air	0.000	09:32
10 Std. Gas	0.080	09:32
11 Room Air	0.000	09:32
12 Std. Gas	0.079	09:33
13 Room Air	0.000	09:33
14 Std. Gas	0.079	09:34
15 Room Air	0.000	09:34
16 Std. Gas	0.079	09:34
17 Room Air	0.000	09:35
18 Std. Gas	0.079	09:35
19 Room Air	0.000	09:36
20 Std. Gas	0.079	09:36
21 Room Air	0.000	09:37
22 Std. Gas	0.080	09:37
23 Room Air	0.000	09:37
24 Std. Gas	0.079	09:38
25 Room Air	0.000	09:38
26 Std. Gas	0.079	09:39
27 Room Air	0.000	09:39
28 Std. Gas	0.079	09:39
29 Room Air	0.000	09:40
30 Std. Gas	0.080	09:40
31 Room Air	0.000	09:41
32 Std. Gas	0.080	09:41
33 Room Air	0.000	09:42
34 Std. Gas	0.080	09:42
35 Room Air	0.000	09:43
36 Std. Gas	0.080	09:43
37 Room Air	0.000	09:43
38 Std. Gas	0.079	09:44
39 Room Air	0.000	09:44
40 Std. Gas	0.080	09:45
41 Room Air	0.000	09:45
42 Cal Check Stats:		
43 Average	0.0793	
44 Std Dev	0.0006	
45 Rel Std Dev	0.73%	



Operator Signature  
ANNA NAREHOOD

Remarks:

Stability check

Form 106-I8000