

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006685	Instrument Location GARR	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory	After Location Change	
Other (Specify) ANNUAL INSPECTION AT CRIME LABORATORY		

Check When Done:

- X 1. Surge Protector Installed/Property Grounded.
- X 2. Telephone Line Connected to Intoxilyzer® 8000.
- X 3. Breath Tube Heated.
- X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- X 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).

X 6. Run Tests:

- X A. Print Test (Level1, Function P).
- X B. ACA Test (Level 1, Function C).
- X C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- X 7. Repair and/or Maintenance Performed (if any): NOT APPLICABLE
- X 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature

Date			
25	ADD	2023	

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	nAll	Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006685 Location = GARR 8164.14.00 09/16 04/25/2023 12:28 *********** Printer Test *********** abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 26021080A1 32 Standard Expiration? 10/05/2023 Oper No? 131749 Flow Cal. Date: 09/06/2016 Slope 705 Intercept -687370 IR Calibration Date: 09/06/2016 3um 9um _____ 0th Coef(*100):-141551st Coef(*100):2714362nd Coef(*100):2016 -22915135296 1251 H2O adj(mg/l*10k): 372 505 ********** Printer Test End **********

> Operator Signature JONATHAN GALLAGHER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-006685Location = GARR8164.14.00 09/16
12:23

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	12:24
02 Std. Gas	0.082	12:24
03 Room Air	0.000	12:24
04 Std. Gas	0.082	12:25
05 Room Air	0.000	12:25
06 Std. Gas	0.081	12:26
07 Room Air	0.000	12:26

Lot No = 26021080A1 Cyl No = 32 Exp Date = 10/05/2023 County = 28

Oper No. = 131749

Operator Signature JONATHAN GALLAGHER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006685Location = GARR8164.14.00 09/1604/25/202312:32

DRY CAL	CHECK	
Test	AC Time	;
01 Room Air 02 Room Air	RFI* 12:32 RFI* 12:32	
*Invalid Test Inhibited - RFI		
Lot No = 26021080A1 Cyl No = 32 Exp Date = 10/05/2023 County = 28	Oper No. = 131749	

Operator Signature JONATHAN GALLAGHER

Remarks:

Form 106-18000