



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4189

A. General Setup and Checks:

- 1. Initial start up
 - Power On
 - 5-minute countdown during initial warm-up
 - 1-minute countdown to begin test after standby mode
- 2. Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
 - Level 1 (Operator)
 - Level 2 (Field Inspector)
 - Level 3 (Lab Personnel)
- 4. Printer Setup (Level 3,S,P)
- 5. Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

- 1. Internal Printer Test (Level 1,P). Sign and attach record.
 - Paper Feed
- 2. External Printer Test (Level 1,P). Sign and attach record.
- 3. Barcode Reader (Level 1,B) (Tone sounds after each)
 - Operator ID
 - Invalid ID
 - ND License
 - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

- 1. Continuous Air Blank (Level 2,A)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 1 of 4

2. Diagnostic (Level 2,D). Sign and attach test record.

3. Prelim Data Entry (Level 2,E)

4. Ser Number/Version Info (Level 2,V):

Inst Ser#: 80-004189 DSP Ver: 1302.14

Core Ver: 8104.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12830 s: 1 C: 47.00

9: 13259 s: 2 B: 45

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0→1

2. Function (F) Keys:

F5 – Tone

F9 - Fan

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 959 CF: 1.05

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 206 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 48.5 I: 0.0

B: 45.0 V: 0.0

Cell Temperature (C) between 47 ± 0.02 °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 387 psi Regulator: 400 psi

Display and Regulator ± 50 psi of each other Yes or No

Gas tank tare necessary? Yes or No

If Yes, display readings after tare (Level 3,M,C,G):

Display: _____ psi Regulator: _____ psi

7. General Setup (Level 3,S,G)

8. Configure Start Test (Level 3,S,T)

9. Network (Level 3,S,N)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 2 of 4

E. Tests (Sign and attach all test records):

1. Configure wet bath simulator for the following test (Level 1,S)
 Zero AC-Wet Calibration Check (Level 1,C):
 Sim ser #: MP3066 Lot #: 33023 Exp date: NIA
 Known value: 0.000 AC
 Results are Zero AC
2. Configure wet bath simulator for the following test (Level 1,S)
 Low AC-Wet Calibration Check (Level 1,C):
 Known Value ≤ 0.03 AC: 0.020 AC
 Sim ser #: MP5318 Lot #: 202201A Exp date: 11Jan24
 Results ± 0.005 of known AC
3. Configure wet bath simulator for the following test (Level 1,S)
 High AC-Wet Calibration Check (Level 1,C):
 Known Value ≥ 0.25 AC: 0.300 AC
 Sim ser #: MP6039 Lot #: 202201F Exp date: 18Jan24
 Results $\pm 5\%$ of known AC
4. Configure dry gas standard for the following tests (Level 1,S)
 Known Value: 0.080 AC
 Gas Cylinder Lot #: 34921080A1 Cylinder #: 49 Exp. Date: 2/5/24
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR3841 Lot #: ICS 7 Exp. Date: NIA
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7349 Lot #: ACE 3 Exp. Date: NIA
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):
 Sim. Ser #: DR7347 Lot #: MET 3 Exp. Date: NIA
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):
 Sim. Ser #: DR7350 Lot #: IPA 3 Exp. Date: NIA
 Display Reads "Interferent Detect"
9. Ambient Fail (Level 1,B): Display read and test record correct.
10. Invalid Sample (Level 1,B): Display read and test record correct.
11. Sequence Aborted (Level 1,B): Display read test record correct.
12. ABA (Level 1,B): Test record correct.

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 3 of 4

- 13. Calibration check out of tolerance (CMS mode): Display read and test record correct.
- 14. RFI Detected (B,C, or CMS mode): Display read and test record correct.
- 15. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 16. Deficient Samples (CMS mode): Display read and test record correct.
- 17. Improper Sample (CMS mode): Display read and test record correct.
- 18. Difference Too Great (CMS mode): Display read and test record correct.
- 19. CMS (CMS mode): Test record correct.
- 20. Stability Test (Level 3,T,S):
 Ave: 0.0815 Std Dev: 0.0005 Rel Std Dev: 0.63%

Maintenance/Repairs/Notes: NIA

Instrument is acceptable to be used in the field. Yes or No

If No, state reason(s) why:

Ann Mill
 Inspector Signature

26 Apr 2023
 Date

Janelli Pertschiller
 Reviewer

27 Apr 2023
 Date

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 4 of 4

Intoxilyzer Test Record and Checklist
NDODAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8154.14.00 09/16
04/25/2023 09:41

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890


Data Entry Mode: enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34921080A1
Standard Cyl #? 47
Standard Expiration? 02/05/2024
Oper No? 133237

Flow Cal. Date: 02/02/2011
Slope 628
Intercept -452608

IR Calibration Date: 05/19/2016
3um 9um

0th Coef(*100): -17767 -18761
1st Coef(*100): 263439 133606
2nd Coef(*100): 1960 1140
H2O adj(ng/l)*10K: 341 362

***** Test End *****


Operator Signature
ANNA NAREHOOD

Remarks:
Internal Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/25/2023 09:46

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34921080A1
Standard Cyl #? 47
Standard Expiration? 02/05/2024
Oper No? 133237

Flow Cal. Date: 02/02/2011
Slope 628
Intercept -452608

IR Calibration Date: 05/19/2016
 3um 9um

0th Coef(*100): -17767 -18761
1st Coef(*100): 263439 133606
2nd Coef(*100): 1960 1140
H2O adj(mg/l*10k): 341 362

***** Printer Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: External Print Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/25/2023 09:48

Test	AC	Time
01 Room Air	0.000	09:49
02 *Subject Test	REF*	09:49
03 Room Air	0.000	09:49

*Subject Test Refused

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = REFUSAL
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 47
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Refusal*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/25/2023 10:06

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/25/2023 13:19

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:19
02 Std. Sol.	0.000	13:20
03 Room Air	0.000	13:21
04 Std. Sol.	0.000	13:21
05 Room Air	0.000	13:22
06 Std. Sol.	0.000	13:23
07 Room Air	0.000	13:23

08 Sim Temp = 34.0°C

Simul Ser No = MP3066
Std Sol No = 33023
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Zero AC check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:04

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:04
02 Std. Sol.	0.018	10:05
03 Room Air	0.000	10:05
04 Std. Sol.	0.020	10:06
05 Room Air	0.000	10:07
06 Std. Sol.	0.021	10:07
07 Room Air	0.000	10:08

08 Sim Temp = 34.0°C

Simul Ser No = MP5318
Std Sol No = 202201A
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Low Ac Check - 0.020Ac*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:11
02 Std. Sol.	0.301	10:12
03 Room Air	0.000	10:13
04 Std. Sol.	0.303	10:13
05 Room Air	0.000	10:14
06 Std. Sol.	0.302	10:15
07 Room Air	0.000	10:15

08 Sim Temp = 34.0°C

Simul Ser No = MP6039
Std Sol No = 202201F
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: High Ac Check - 0.300Ac

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:17

Test	AC	Time
01 Room Air	0.000	10:18
02 *Subject Test	INT*	10:18
03 Room Air	0.000	10:19

*Invalid Test
Interferent Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = INTERFERENT CK
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
 ETOH + Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:19

Test	AC	Time
01 Room Air	0.000	10:20
02 *Subject Test	INT*	10:20
03 Room Air	0.000	10:21

*Invalid Test
Interferent Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = INTERFERENT CK
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:21

Test	AC	Time
01 Room Air	0.000	10:22
02 *Subject Test	INT*	10:22
03 Room Air	0.000	10:23

*Invalid Test
Interferent Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = INTERFERENT CK
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *interferent check*
MeOH

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:23

Test	AC	Time
01 Room Air	0.000	10:24
02 *Subject Test	INT*	10:25
03 Room Air	0.000	10:25

*Invalid Test
Interferent Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = INTERFERENT CK
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent Check
IPA

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:27

Test	AC	Time
01 Room Air	AMB*	10:27
02 Room Air	0.000	10:28

*Invalid Test
Check Ambient Conditions

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = AMBIENT FAIL
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Ambient Fail

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:29

Test	AC	Time
01 Room Air	0.000	10:29
02 Invalid Sample	X.XXX	10:30
03 Room Air	0.000	10:30

*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = INVALID
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Invalid Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:31

Test	AC	Time
01 Room Air	0.000	10:31
02 *Subject Test	ABT*	10:32
03 Room Air	0.000	10:32

*Sequence Aborted

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = SEQUENCE ABORT
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Sequence Aborted*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:33

Test	AC	Time
01 Room Air	0.000	10:33
02 Subject Test 1	0.000	10:36
03 Room Air	0.000	10:37
04 Reported AC	0.000	10:36

No RFI Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = ABA
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *aba test*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:37

Test	AC	Time
01 Diagnostic	OK	10:38
02 Room Air	0.000	10:38
03 Subject Test 1	0.000	10:39
04 Room Air	0.000	10:41
05 Std. Gas	0.058*	10:42
06 Room Air	0.000	10:43

*Cal Check Out of Tolerance

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = CAL CHECK OOT
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Calibration check out
of tolerance

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:43

Test	AC	Time
01 Room Air	0.000	10:44
02 *Subject Test	RFI*	10:44
03 Room Air	0.000	10:45

*Invalid Test
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = RFI CHECK
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 10:45

Test	AC	Time
01 Diagnostic	OK	10:46
02 Room Air	0.000	10:46
03 *Subject Test	0.000*	10:49
04 Room Air	0.000	10:54
05 Std. Gas	0.080	10:55
06 Room Air	0.000	10:56
07 *Subject Test	0.000*	10:59
08 Room Air	0.000	11:00

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = DEFICIENT SAMPL
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.

Operator Signature
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 11:01

Test	AC	Time
01 Diagnostic	OK	11:02
02 Room Air	IPS*	11:02
03 Room Air	0.000	11:02

*Improper Sample

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = IMPROPER
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: *Improper Sample*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 11:03

Test	AC	Time
01 Diagnostic	OK	11:03
02 Room Air	0.000	11:04
03 Subject Test 1	0.000*	11:04
04 Room Air	0.000	11:07
05 Std. Gas	0.080	11:08
06 Room Air	0.000	11:09
07 Subject Test 2	0.042*	11:09
08 Room Air	0.000	11:10

*Difference Too Great

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = DIFF TOO GREAT
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Difference Too Great

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 11:11

Test	AC	Time
01 Diagnostic	OK	11:12
02 Room Air	0.000	11:12
03 Subject Test 1	0.000	11:12
04 Room Air	0.000	11:14
05 Std. Gas	0.077	11:16
06 Room Air	0.000	11:17
07 Subject Test 2	0.000	11:17
08 Room Air	0.000	11:19
09 Reported AC	0.000	11:12

Difference OK
No RFI Detected

Sub Name = TEST, DONOR2 NONE
Sub DOB = 07/25/1998
Sub Sex = Male Weight = NA
Test = OTH Cit = CMS TEST
Dr. Lic. = ND/TES989643
Lot No = 34921080A1
Cyl No = 49
Expiration Date = 02/05/2024
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: CMS Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/26/2023 11:20

Test	AC	Time
01 Room Air	0.000	11:20
02 Std. Gas	0.081	11:21
03 Room Air	0.000	11:21
04 Std. Gas	0.082	11:22
05 Room Air	0.000	11:22
06 Std. Gas	0.081	11:22
07 Room Air	0.000	11:23
08 Std. Gas	0.081	11:23
09 Room Air	0.000	11:24
10 Std. Gas	0.082	11:24
11 Room Air	0.000	11:25
12 Std. Gas	0.082	11:25
13 Room Air	0.000	11:26
14 Std. Gas	0.082	11:26
15 Room Air	0.000	11:26
16 Std. Gas	0.082	11:27
17 Room Air	0.000	11:27
18 Std. Gas	0.081	11:28
19 Room Air	0.000	11:28
20 Std. Gas	0.081	11:29
21 Room Air	0.000	11:29
22 Std. Gas	0.081	11:29
23 Room Air	0.000	11:30
24 Std. Gas	0.082	11:30
25 Room Air	0.000	11:31
26 Std. Gas	0.082	11:31
27 Room Air	0.000	11:32
28 Std. Gas	0.082	11:32
29 Room Air	0.000	11:33
30 Std. Gas	0.081	11:33
31 Room Air	0.000	11:33
32 Std. Gas	0.081	11:34
33 Room Air	0.000	11:34
34 Std. Gas	0.082	11:35
35 Room Air	0.000	11:35
36 Std. Gas	0.082	11:35
37 Room Air	0.000	11:36
38 Std. Gas	0.081	11:36
39 Room Air	0.000	11:37
40 Std. Gas	0.081	11:37
41 Room Air	0.000	11:38
42 Cal Check Stats:		
43 Average	0.0815	
44 Std Dev	0.0005	
45 Rel Std Dev	0.63%	

Operator Signature
ANNA NAREHOOD

Remarks:

Stability Test

Form 106-I8000