



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59281 (06/2018)

Serial Number <i>80-004942</i>	Instrument Location <i>Stanton</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u><i>Annual Inspection at Crime Laboratory</i></u>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): *N/A*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>4/26/2023</i>
---	--------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date 28Apr2023
---	------------------------------------


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = STAT 8164.14.00 09/16
04/26/2023 18:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:12
02 Std. Gas	0.084	18:13
03 Room Air	0.000	18:13
04 Std. Gas	0.083	18:14
05 Room Air	0.000	18:14
06 Std. Gas	0.083	18:14
07 Room Air	0.000	18:15

Lot No = 26021080A1
Cyl No = 1
Exp Date = 10/05/2023
County = 29 Oper No. = 120906



Operator Signature
PAMELA PEHL

Remarks:

ACA Test

Form 106-18000

- High ACA values. Contacted Pam,
she ran a continuous air
blank for 1/2 hour and reran.
Results were lower. See attached
document. AEN 28Apr2023

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = STAT 8164.14.00 09/16
04/27/2023 19:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:39
02 Std. Gas	0.082	19:39
03 Room Air	0.000	19:40
04 Std. Gas	0.082	19:40
05 Room Air	0.000	19:41
06 Std. Gas	0.082	19:41
07 Room Air	0.000	19:42

Lot No = 26021080A1
Cyl No = 1
Exp Date = 10/05/2023
County = 29 Oper No. = 120906



Operator Signature
PAMELA PEHL

Remarks:

ACA Test

Form 106-I8000