

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	Λ σ
80-006490	mstrument Location	ASHL
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install A	After Location Change	
Dother (Specify) Annual Inspection	at Crime Labor	atory.
Check When Done:		O
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	ocation; Level 2, Function E).	
✓ 5. Scan/Enter Gas Cylinder Information (Level 1, Function)	S).	
☑ 6. Run Tests:		
A Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mod	le or Level 1, Function B or C;	Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (Sf	FN50496, Form 120-G) and Pla	ace it by the Intoxilyzer® for Use.
	20-G) at the Intoxilyzer® Locati	on at the Agency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	d Intoxilyzer® 8000 Installation	and Repair Checkout (SFN59281, Form
Field Inspector Signature		Date 4-26-23
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved	to be used for the analysis of h	preath to determine alcohol concentration from
the date the Field Inspector performed the installation. This record County of Burleigh, North Dakota, is certified to be a true and corre	on file at the Office of Attorney	General, Crime Laboratory Division, in the
Reviewed/Certified By		Certified Date
(iu/hl		27Apr2023

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006490
Location = ASHL 8164.14.00 09/16
04/26/2023 18:34

\*\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes 1 # of Print Copies? Select Std (D/W/I)? Dry 0.080 Standard Value? 01923080A3 Standard Lot #? Standard Cyl #? 62 Standard Expiration? 02/05/2025 133262 Oper No?

Flow Cal. Date: 09/04/2020 Slope 649 Intercept -598817

IR Calibration Date: 09/04/2020 3um 9um

Coef(*100):	-25190	-16392
	271645	136906
Coef(*100):	2533	1222
adj (mg/l*10k):	460	309
	Coef(*100): Coef(*100): Coef(*100): adj(mg/l*10k):	Coef(*100): 271645

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature VICTOR FISCHER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ASHL
04/26/2023

Alcohol Analyzer SN 80-006490 8164.14.00 09/16 18:52

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:53
02 Std. Gas	0.080	18:53
03 Room Air	0.000	18:54
04 Std. Gas	0.080	18:54
05 Room Air	0.000	18:55
06 Std. Gas	0.079	18:55
07 Room Air	0.000	18:56

Lot No = 01923080A3

Cyl No = 62

Exp Date = 02/05/2025

County = 26

Oper No. = 133262

Operator Signature VICTOR FISCHER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = ASHL
04/26/2023

Alcohol Analyzer SN 80-006490 8164.14.00 09/16 18:58

Test	AC	Time
01 Diagnostic	OK	19:00
02 Room Air	RFI*	19:00
03 Room Air	RFI*	19:00

\*Invalid Test Inhibited - RFI

Sub Name = BOSS, BOSS BOSS

Sub DOB = 02/02/2022

Sub Sex = Female

Weight = 111 Cit = 11

Test = DUI

Dr. Lic. = ND/11 Lot No = 01923080A3

Cyl No = 62

Expiration Date = 02/05/2025

County = 26

Oper No. = 133262

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature VICTOR FISCHER

Remarks:

Form 106-I8000