



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-004935</i>	Instrument Location <i>EMMONS COUNTY</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <i>Annual Inspection at Crime Laboratory</i>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>4/26/23</i>
---	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.


Reviewed/Certified By <i>Aime Nahl</i>	Certified Date 28Apr2023
---	-----------------------------

Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-004935  
 Location = LINT 8164.14.00 09/16  
 04/26/2023 17:52

Test	AC	Time
01 Room Air	0.000	17:54
02 *Subject Test	RFI*	17:54
03 Room Air	0.000	17:55

\*Invalid Test  
 Inhibited - RFI

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-004935  
 Location = LINT 8164.14.00 09/16  
 04/26/2023 17:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:49
02 Std. Gas	0.080	17:49
03 Room Air	0.000	17:50
04 Std. Gas	0.079	17:50
05 Room Air	0.000	17:51
06 Std. Gas	0.078	17:51
07 Room Air	0.000	17:52

Lot No = 34921080A1  
 Cyl No = 56  
 Exp Date = 02/05/2024  
 County = 15  
 Oper No. = 130117

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-I8000

Current Instrument Setup  
 Data Entry Mode: Enabled  
 Start Test Sequence: DABACABA  
 Display Prelim Rslt? Yes  
 Display Third Digit? Yes  
 Inhib Printer(Y/N)? No  
 Display Volume? No  
 Disable On MenFull? Yes  
 # of Print Copies? 1  
 Select Std (D/W/I)? Dry  
 Standard Value? 0.080  
 Standard Lot #? 34921080A1  
 Standard Cyl #? 56  
 Standard Expiration? 02/05/2024  
 Oper No? 130117

Flow Cal. Date: 07/12/2012  
 Slope: 655  
 Intercept: -51.0289

IR Calibration Date: 05/20/2015  
 3um 9um

0th Coef(\*100): -7871 -19684  
 1st Coef(\*100): 264505 133108  
 2nd Coef(\*100): 2865 1320  
 H2O adj(mg/l\*10K): 200 298

\*\*\*\*\* Test End \*\*\*\*\*

  
 Operator Signature  
 MARK PEARSON

Remarks:

Form 106-I8000

42 characters per line  
 abcdefghijklmnopqrstuvwxyz1234567890abcde  
 fghijklmnopqrstuvwxyz1234567890abcdef

32 characters per line  
 abcdefghijklmnopqrstuvwxyz123456  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ123456

24 characters per line  
 abcdefghijklmnopqrstuvw  
 ABCDEFGHIJKLMNOPQRSTUWX

21 chars per line  
 abcdefghijklmnopqrstu  
 ABCDEFGHIJKLMNOPQRSTU

16 char per line  
 abcdefghijklmnop  
 ABCDEFGHIJKLMNOP

12 char/line  
 abcdefghijkl  
 ABCDEFGHIJKL

\*\*\*\*\* Printer Test \*\*\*\*\*  
 Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501  
 CMI, Inc. Intoxilyzer Alcohol Analyzer  
 North Dakota Model 8000 SN 80-004935  
 Location = LINT 8164.14.00 09/16  
 04/26/2023 17:44