

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	1	
Serial Number 80 - 905 36 l	Instrument Location FTOT \$16.4	1,14.60
Reason for Install/Repair	770.	,71,00
☐ Install After Receiving From Crime Laboratory ☐ Install Aft	er Location Change	
	. 1	
**XOther (Specify) Annual Inspection at	Crime Laboratory	
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Local	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode of	Level 1, Function B or C; Key Radio Du	uring Test).
☑ 7. Repair and/or Maintenance Performed (if any): ✓ IA		
	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date
(1) "		4-28-2023
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on f County of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Crir	
Reviewed/Certified By		Certified Date
		200 200

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005361
Location = FTOT	8164.14.00 09/16
04/28/2023	10:25

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current	Instrument	Setup
Markov Co. Markov Co. Residence Co.	and the second s	200

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	34921080A1
Standard Cyl #?	41
Standard Expiration?	02/05/2024
Oper No?	132596

Flow Cal. Date:	05/15/2020
Slope	634
Intercept	-488053

IR	Calibration	Date:	05/15/2	020
		3um		9um
0 1 1	C (-1-100)	2520	2.2	07000

Oth	Coef(*100):	-25383	-27323
1st	Coef(*100):	288661	141732
2nd	Coef(*100):	458	822
H20	adj(mg/l*10k):	335	405

****** Printer Test End *******

Operator Signature DAVID LAWSON

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FTOT
04/28/2023

Alcohol Analyzer SN 80-005361 8164.14.00 09/16 10:26

DRY CAL CHECK

$T\epsilon$	est		AC	Time
01	Room	Air	0.000	10:26
02	Std.	Gas	0.081	10:27
03	Room	Air	0.000	10:27
04	Std.	Gas	0.081	10:27
05	Room	Air	0.000	10:28
06	Std.	Gas	0.080	10:28
07	Room	Air	0.000	10:29

Lot No = 34921080A1

Cyl No = 41

Exp Date = 02/05/2024

County = 62

Oper No. = 132596

Operator Signature DAVID LAWSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005361
Location = FTOT	8164.14.00 09/16
04/28/2023	10:29

Test	AC	Time
01 Room Air 02 *Subject Test	0.000 RFT*	10:31
03 Room Air	0.000	10:32

^{*}Invalid Test Inhibited - RFI

Oper No. = 132596

Operator Signature DAVID LAWSON

Remarks:

County = 62

Form 106-I8000