

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location Minot Police Department
Reason for Install/Repair	
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change
Decity) Annual Inspection at	- Crime laboratory
	0
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
☑ 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).
☑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
⊠ 6. Run Tests:	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio During Test).
X 7. Repair and/or Maintenance Performed (if any):	
▼ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.
▼ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0)	3) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Inf 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date 4/28/23
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct or	
Reviewed/Certified By	Certified Date
Une I'm	03May2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006503
Location = MIPD	8164.14.00 09/16
04/28/2023	19:24

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current	Instrument	Setup

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	34921080A1
Standard Cyl #?	31
Standard Expiration?	02/05/2024
Oper No?	132527
<del></del>	

Flow Cal. D	ate:	08/20/2015
Slope		663
Intercept		-633300

IR	Calibration	Date:	06/12/2017
		3um	9um

	-		
0th	Coef(*100):	-6741	-15632
1st	Coef(*100):	256583	132194
2nd	Coef(*100):	3028	1131
H20	adj(mg/l*10}	(): 354	358

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature SETH LABODDA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MIPD
04/28/2023

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 19:25

## DRY CAL CHECK

T	est		AC	Time
01	Room	Air	0.000	19:25
02	Std.	Gas	0.082	19:26
03	Room	Air	0.000	19:26
04	Std.	Gas	0.081	19:26
05	Room	Air	0.000	19:27
06	Std.	Gas	0.082	19:27
07	Room	Air	0.000	19:28

Lot No = 34921080A1

Cyl No = 31

Exp Date = 02/05/2024

County = 51

Oper No. = 132527

Operator Signature SETH LABODDA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MIPD
04/28/2023

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 19:29

Test	AC	Time
01 Room	 RFI*	19:30
02 Room	0.000	19:30

\*Invalid Test Inhibited - RFI

Sub Name = LABODDA, SETH AARON

Sub DOB = 09/30/1993

Sub Sex = Male

Weight = 160

Test = OTH

Cit = NA

Dr. Lic. = ND/LAB931101

Lot No = 34921080A1

Cyl No = 31

Expiration Date = 02/05/2024

County = 51

Oper No. = 132527

Operator Signature SETH LABODDA

Remarks:

Form 106-I8000