

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location			
80-004944 Reason for local Ulliston			
Reason for Install/Repair			
Install After Receiving From Crime Laboratory Install After Location Change			
Dether (Specify) Annual Inspection at Chine Laboratory.			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio D	During Test).		
7. Repair and/or Maintenance Performed (if any): N/A			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Ag			
 Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair 104-G), Print Test, ACA Test, and RFI Test. 	Checkout (SFN59281, Form		
Field Inspector Signature	Date 5-8-23		
Crime Laboratory Use Only			

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	<u></u>	- <u>A</u> A	
interviewed/Certilied by	N	N 0	Certified Date
	-Um	That	09May2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004944Location = WILL8164.14.00 09/1605/08/202320:07

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 26021080A1 Standard Cyl #? 13 Standard Expiration? 10/05/2023 Oper No? 130739 Flow Cal. Date: 07/01/2011 Slope 685 Intercept -634740

IR Calibration Date: 07/01/2011 3um 9um Oth Coef(*100): -22591 -17234 1st Coef(*100): 265915 132691 2nd Coef(*100): 2841 1562 H2O adj(mg/l*10k): 559 425

********** Printer Test End **********

Operator Signature MICHELLE ROMANS

Remarks:

Print test

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004944Location = WILL8164.14.00 09/1605/08/202320:08

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	20:09 20:10 20:10 20:10 20:11 20:11 20:11

Operator Signature MICHELLE ROMANS

Remarks:

ACA - valiz

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004944 Location = WILL 8164.14.00 09/16 05/08/2023 20:13 Test AC Time 01 Room Air RFI* 20:14 02 Room Air 0.000 20:15 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 01/01/2000Sub Sex = Unknown Weight = NA Test = DUI Cit = NADr. Lic. = ND/NA Lot No = 26021080A1Cyl No = 13Expiration Date = 10/05/2023County = 53Oper No. = 130739

Operator Signature

MICHELLE ROMANS

Remarks:

Form 106-18000