	NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION				
	INTOXILYZER <sup>®</sup> 8000 ANNUAL	INSPECTION			
Intox	kilyzer <sup>®</sup> 8000 Serial Number: 80-00 0507	Inspection Location: TOXL			
Α.	<ul> <li>Pre-Inspection</li> <li>1. Items with Instrument: Gas Cylinder Yes or No (If Yes, Lot # Keys Yes or No Power Cord Yes or No</li> <li>2. Download Data</li> <li>3. Upload Operator File</li> <li>4. Current Location Code: LAMO</li> <li>5. Battery Check Was the external battery pack replace</li> <li>6. O-Rings Replaced Simulator O-Ring</li> </ul>				
	Replaced Breath Tube O-Ring	Yes or No			
Β.	<ul> <li>General Setup and Checks:</li> <li>1. A Diagnostics passed and instrument in "Re</li> <li>2. Breath tube heated</li> <li>3. Date, time and location code (Level 2,E). Time Zone: CST or MDT (Time on te</li> <li>4. Print test (Level 1,P). Sign and attach tes</li> <li>5. Tank monitor (Level 3,D,G). Display: <u>192</u> psi Regulator: Display and Regulator ± 50 psi of eac Gas tank tare necessary? If Yes, display readings after tare (Lev Display: psi Regulator:</li> </ul>	Re-set if necessary. est records will be in time zone circled). st record. <u>200</u> psi th other Yes or No Yes or No Yes or No			
C.	1. Configure simulator for the following test (Level 1,S). Wet Calibration Check - Low AC (Level 1,C) Known Value $\leq 0.03$ AC: $0.020$ AC Sim. Ser #: <u>MP5318</u> Lot #: <u>2022014</u> Exp. Date: <u>ILJan X4</u> Results ± 0.005 of known AC				
Toxicology Unit - Breath Alcohol Section Qualtrax ID: 11698 Revision: 1					
Intoxilyzer 8000 Annual Inspection Status: Published: 03/24/					
UNCO	NTROLLED WHEN PRINTED	Page 1 of 2			

Uploaded 10May2023

2.	Configure simulator for the follow Wet Calibration Check - High Known Value ≥ 0.25 AC: Sim. Ser #:NP300 Lot #:202201 Exp. Date:18Jan 20 Results ± 5% AC of known	AC (Level 1,C) <u>300</u> AC <del>7</del> <del>F</del>		
3.	Configure dry gas standard for the Known Value : <u>0.080</u> Ac Gas Cylinder Lot #: <u>260</u> Cylinder #: <u>16</u> Exp. Date: <u>10523</u>	e remaining tests (Level 1,S). C <u>A IO 80 IA I</u>		
4. 5.	Interferent Check (Level 1,B) Known Value: <u>0.10</u> AC + 0.05 Sim. Ser #: <u>DR7345</u> Lot #: <u>ICS 7</u> Exp. Date: <u>NIA</u> Display reads "Interferent D RFI Check (CMS Mode)	% Acetone		
6.	Display reads "RFI Detect" Display reads "RFI Detect" Dry Calibration Check (Level 1,C) Test 1 <u>0.081</u> Test 4 <u>0</u> Test 2 <u>0.080</u> Test 5 <u>0</u> Test 3 <u>0.081</u> Test 6 <u>0</u> Average	<u>080</u> Test 8 <u>0.081</u> <u>.080</u> Test 9 <u>0.080</u>		
D. Rem	arks/Maintenance: NA			
Instrument i	s acceptable to be used in the field. Ye	as or No		
	reason(s) why:			
in No, State				
If Yes, chan	ge location code back to A.4.	000000000		
Inspector Signature Date				
Reviewer	Portschiller	lo May 2023 Date		
Toxicology Ur	nit - Breath Alcohol Section	Qualtrax ID: 11698 Revision: 1		
Intoxilyzer 800	0 Annual Inspection	Status: Published Published: 03/24/2023		
UNCONTROLI	LED WHEN PRINTED	Page 2 of 2		

Uploaded 10May2023

AEN

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202313:28

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= | ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	33321080A4
Standard Cyl #?	10
Standard Expiration?	01/05/2024
Oper No?	133237
Flow Cal. Date:	08/20/2015
Slope	659
Intercept	-658315
IR Calibration Date:	08/20/2015
IR CALIDIALION Date:	00/20/2015

3um9um0th Coef(\*100):-9452-117221st Coef(\*100):2588081395882nd Coef(\*100):44761674H2O adj(mg/l\*10k):499376

Operator Signature ANNA NAREHOOD

Remarks:

Print Test

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202313:42

		WET CAL CHECK	
Te	est	AC	Time
01	Room Air	0.000	13:43
02	Std. Sol.	0.022	13:44
03	Room Air	0.000	13:45
04	Std. Sol.	0.022	13:45
05	Room Air	0.000	13:46
06	Std. Sol.	0.022	13:46
07	Room Air	0.000	13:47
08	Sim Temp =	34.0°C	

Simul Ser No = MP5318 Std Sol No = 202201A County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

LOW an cherk-0.020AC

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202313:48

Time
13:48
13:49
13:50
13:51
13:51
13:52
13:52

08 Sim Temp = 34.0°C

Operator Signature ANNA NAREHOOD

High AC Check- 0.300AC

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006507 Location = TOXL 8164.14.00 09/16 05/09/2023 13:56

Test	AC	Time
01 Room Air	0.000	13:56
02 *Subject Test	INT*	13:57
03 Room Air	0.000	13:57

\*Invalid Test Interferent Detected

Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998 Sub Sex = Male Weight = NA Test = OTH Cit = INTERFERENT CK Dr. Lic. = ND/TES989643 Lot No = 26021080A1 Cyl No = 16 Expiration Date = 10/05/2023 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Interferent Check

NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006507 Location = TOXL 8164.14.00 09/16 05/09/2023 13:58 Test AC Time 01 Diagnostic OK 13:58 02 Room Air 0.000 13:59 03 \*Subject Test RFI\* 13:59 04 Room Air 0.000 13:59 \*Invalid Test Inhibited - RFI Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Test = OTHCit = RFI CHECK Dr. Lic. = ND/TES989643 Lot No = 26021080A1Cyl No = 16Expiration Date = 10/05/2023 County = 08Oper No. = 133237

Intoxilyzer Test Record and Checklist

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks:

RFI Chelk

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202314:00

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	14:00
02 Std. Gas	0.081	14:01
03 Room Air	0.000	14:01
04 Std. Gas	0.080	14:02
05 Room Air	0.000	14:02
06 Std. Gas	0.081	14:03
07 Room Air	0.000	14:03

Operator Signature ANNA NAREHOOD

Remarks:

Dry cal Cher #1-3

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202314:04

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	14:05
02 Std. Gas	0.080	14:05
03 Room Air	0.000	14:06
04 Std. Gas	0.080	14:06
05 Room Air	0.000	14:06
06 Std. Gas	0.080	14:07
07 Room Air	0.000	14:07

Operator Signature ANNA NAREHOOD

Remarks:

Cal check #4-6

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006507Location = TOXL8164.14.00 09/1605/09/202314:07

		ĺ	DRY CAL CHECK	]
Т	est		AC	Time
01	Room	Air	0.000	14:08
02	Std.	Gas	0.081	14:08
03	Room	Air	0.000	14:09
04	Std.	Gas	0.081	14:09
05	Room	Air	0.000	14:10
06	Std.	Gas	0.080	14:10
07	Room	Air	0.000	14:11

Lot No = 26021080A1 Cyl No = 16 Exp Date = 10/05/2023 County = 08 Oper

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Dry cal check #7-9